USER MANUAL HEMODIALYSIS CHAIR

Attention: Carefully read this manual
before operating the Olsen's equipment!

Olsen
Equipment made to last

Index

1 - Introduction	5
2 - Package Contents	5
3 - Equipment Presentation	5
3.1 - Standard Items	5
3.2 - Optional Items	5
3.3 - Hemodialysis Chair	6
3.4 - Synchronized Hemodialysis Chair	6
4 - Parts Identification	6
4.1 - Applied Parts	6
4.2 - Accessories and Detachable Parts	6
5 - Description and Operating of the Equipment	7
5.1 - Turning On the Equipment	7
5.2 - Auxiliary Pillow	7
5.3 - Hand Control	7
5.4 - Chair Movements	8
5.4.1 - Synchronized Hemodialysis Chair	8
5.4.2 - Hemodialysis Chair	8
5.4.3 - Positioning of the Patient, the Operator and Others	8
5.4.4 - Movements Interruption	9
5.5 - Independent Backrest	9
5.6 - Upholstered Armrests	9
5.7 - Armrest for Blood Donation (Optional)	9
5.8 - Left and Right Bracket Locks	10
5.9 - Stainless-Steel Tray (Optional)	10
5.10 - Rear Castor with Lock	11
6 - General Features	11
7 - Installation Requirements	11
7.1 - Pre-installation	11
7.2 - Electrical Installation	11
8 - Installation	12
8.1 - Position of the Equipment	12
8.2 - Olsen Certified Technical Assistance Network	12
9 - Cleaning and Disinfection	12
9.1 - Upholsteries and Plastic Covers	12
9.2 - Painted Parts	12

9.3 - Sterilization by Autoclave.....	13
10 - Technical Features	14
10.1 - Electromagnetic Compatibility (EMC)	15
11 - Dimensional	17
12 - Symbology.....	18
13 - Important Notes.....	18
13.1 - General Cares - Compulsory Reading	19
13.2 - Disposal.....	20
13.3 - Transport and Storage.....	20
13.4 Contraindication.....	20
13.5 - Finishing.....	20
14 - Troubleshooting.....	20
15 - Preventive Review	21
16 - Warranty Terms	22
17 - Message from the President.....	23

1 - Introduction

Congratulations on the great choice!

You have acquired an equipment developed for the outpatient environment, aimed primarily at the specialty of nephrology, oncology and hematology with the aim of providing the maximum comfort to the patient during the procedures.

This manual is supplemented to the Hemodialysis Chair Quick Operation Guide provided with your equipment and provides all the information you need to get the most out of your equipment, so read them carefully before using it.

2 - Package Contents

Check out the equipment package contents:

Standard Items:

- 1 Hemodialysis Chair
- 1 Left Bracket Lock
- 1 Right Bracket Lock
- 2 Upholstered Armrests
- 1 Auxiliary Pillow
- 1 Hemodialysis Chair Quick Operation Guide

Optional Items:

- 1 Armrest for Blood Donation
- 1 Stainless-Steel Tray

3 - Equipment Presentation

3.1 - Standard Items

- Automatic Chair with 4 Commands
- Upholstered Armrests
- Hand Control
- Auxiliary Pillow
- 2 Oil-Free Bosch Motors
- Synchronized Movement of Seat and Backrest **OR** Independent Movements
- Swivel Castors with Individual Locking

3.2 - Optional Items

- Armrest for Blood Donation
- Stainless-Steel Tray

3.3 - Hemodialysis Chair

Independent Backrest and Seat Movements.

3.4 - Synchronized Hemodialysis Chair

Synchronous Seat and Backrest Movement.

4 - Parts Identification

- A - Swivel Castors with Individual Locking
- B - Hand Control
- C - Auxiliary Pillow
- D - Backrest
- E - Bracket Lock
- F - Upholstered Armrests
- G - Seat
- H - Electrical Panel

4.1 - Applied Parts

The following items of the Hemodialysis Chair are considered as applied parts to the patient:

Standard Items:

- Upholstery;
- Auxiliary Pillow;
- Upholstered Armrests;

Optional Items:

- Armrest for Blood Donation.

4.2 - Accessories and Detachable Parts

The following items of the Hemodialysis Chair are considered as accessories or detachable parts or:

Detachable parts:

- Auxiliary Pillow;
- Upholstered Armrests.

Accessories:

- Armrest for Blood Donation;
- Stainless-Steel Tray.

5 - Description and Operating of the Equipment

Technical Name: Hemodialysis Chair

Trade Name: Hemodialysis Chair

The Hemodialysis Chair is an automated equipment designed for patient accommodation on clinical and outpatient exams, mainly in the areas of nephrology, oncology and hematology.

It is equipped with 2 electrical motors to perform 4 individual movements, controlled by hand control.

It can easily be moved because it has high resistant castors manufactured with molded polyurethane. These castors protect the floor and have a very low noise when they are on displacement.

This equipment should not be used to patient transport.

Before starting use your equipment make sure that it is properly installed, connected to the mains, the circuit breaker is on and the mains is energized.

5.1 - Turning On the Equipment

Located on the chair base, under the footrest and easy access, the electrical panel provides access to the protection fuses (3), the On / Off switch (1) and indicator LED (2)


To turn On the equipment: press the On/Off switch (1) to the "I" position. The LED (2) next to the On/Off switch (1) will turn green, indicating that the equipment is ready for use

To turn Off the equipment: press the On/Off switch (1) to the "O" position.

5.2 - Auxiliary Pillow


The auxiliary pillow manufactured with the same upholstery of the Hemodialysis Chair, for the purpose of increasing patient comfort. has 2 strips with Velcro® (1) for easy positioning and removal for cleaning.

5.3 - Hand Control

It presents 4 commands with identified icons for assist the operation. The control has a holder (5) for allocation in the side rails of the chair.

- 1 - Seat Up
- 2 - Seat Down
- 3 - Backrest Down
- 4 - Backrest Up

5.4 - Chair Movements

5.4.1 - Synchronized Hemodialysis Chair

Commands 1 and 2: seat height adjustment movement.

Commands 3 and 4: simultaneous movement of backrest and seat.

5.4.2 - Hemodialysis Chair

Commands 1 and 2: seat height adjustment movement.

Commands 3 and 4: backrest adjustment movement.

5.4.3 - Positioning of the Patient, the Operator and Others

During the equipment movement, the operator and other persons close to him/her should position themselves on the equipment sides, respecting the minimum distance of 50 centimeters, remaining out of the equipment movement area and its components.

To ensure safe and proper patient positioning, the operator should advise him to keep his hands in sight during the entire equipment movement on the equipment upholstery or on the specific accessories for performing the procedure. If the operator sees any situation that may generate any kind of risk to the patient, the chair must be stopped immediately.

The presence of the operator, persons or objects in the equipment movement areas and its components may cause damage to the equipment and/or impair its correct operation.

All equipment's movements must be supervised by the operator. If the operator sees any situation that may generate any kind of risk to the patient, the chair must be stopped immediately.

5.4.4 - Movements Interruption

All Hemodialysis Chair movements are of continuous activation, i.e., while the operator presses the control button, the equipment realizes the movement. The movement is interrupted immediately when the operator releases the control button.

5.5 - Independent Backrest

The Hemodialysis Chair can present the independent backrest, that means, the movements of the backrest are performed without the synchronized movement of the seat (Synchronized Hemodialysis).

In this configuration the equipment offers the commands Seat Up, Seat Down, Backrest Up and Backrest Down. There is no independent seat tilt control.


5.6 - Upholstered Armrests

Manufactured with the same upholstery of the Hemodialysis Chair, they feature soft foam to provide greater comfort to the patient. With the aid of the bracket locks, the upholstered armrest (1) feature height and depth adjustment.

The armrest features upholsteries (1) fitted on the metal support (2), allowing a free spin on the fitting shaft (3).

The fitting shaft has joint (4), which gives the piece the inclination adjustment.

To adjust the inclination: position the upholstery on the requested position.



5.7 - Armrest for Blood Donation (Optional)

Are used to position the patient's arms to facilitate the administration of intravenous drugs or blood donation, the armrest for blood donation are concave shaped and are made of ABS for better asepsis.

The armrest for blood donation has a 36 cm rod.

Can be adjusted in the same way as the armrest (section 5.6), that means, through the metal supports. Allowing adjustment for height, distance from center of seat and distance from backrest.

5.8 - Left and Right Bracket Locks

Used to attach accessories to the side rails, the supports are identified according to the side of the chair to be installed. The support on the left side (1) receives the letter L and the right (2) receives the letter R.

To adjust the height (6) or turning (7) of the accessory: turn the plastic handle (3) counterclockwise. Adjust the accessory to the desired position and tighten the plastic handle (3) by turning it clockwise.

To adjust the position on the rails (8): pull the metal lever (4) counterclockwise by loosening the bracket on the side rail. Adjust the accessory to the desired position and turn the metal lever (4) clockwise.

Note: the metal lever (4) moves only in the area indicated in the figure to the side. To turn it the metal lever (4) clockwise, push it upwards causing the lever to pass through the metallic knob's axis (5), and then rotate it again.


5.9 - Stainless-Steel Tray (Optional)

Attached to the chair base (4), the stainless-steel tray support assembly remains stable while the chair executes the movements.

The set has metal support (2) and stainless-steel tray (1), which allows its cleaning by autoclave, as well as allows free rotation on its axis and can be adjusted for height and distance in relation to the backrest.

To adjust the height or distance from the seat: turn the plastic handle (3) counterclockwise. Set the stainless-steel tray to the desired position and tighten the plastic handle (3) by turning it clockwise.

5.10 - Rear Castor with Lock

The rear castors lock can be activated with the feet.

To lock the castor: press the lever (1) down.

To release the castor lock: push the lever (2) down.

The lock must remain active during the equipment use. Only disable the lock to move the equipment.

6 - General Features

- **Upholsteries:** are assembled on a very resistant frame, with soft foam and coated with flexible laminated PVC or leather, which besides being resistant, providing easy asepsis for the whole.
- **Mechanical structure:** manufactured in SAE 1020 rolled steel profiles and welded by MIG process, to guarantee the strength and durability.
- **Electrostatic painting:** it is applied to all structural metal parts of the equipment. The special polyurethane paint provides the equipment with high durability coating and has antibacterial properties according to *JIS Z 2801: 2000*, where in 24 hours, the bacterial reduction is greater than 99.9%.
- **Plastic covers:** made of high-strength ABS with acrylic cover, the plastic covers do not require painting and allowing polishing for recovery in case of wear or scratches.
- **Electrical system:** the equipment can work in frequencies of 50 or 60 Hz and can be configured to be connected to the voltages of 118/127/220/230 Volts by an authorized technician. The maximum supply voltage of electronic boards, motors, and other controls is 24 V. The electrical system has On/Off switch and protection fuses.
- **Spindle drive motor:** all motors used in the equipment are produced by Robert Bosch of Brazil and have advantages such as low noise, the absence of oil reservoir, uniformity in displacement, reduction of energy consumption and low cost for maintenance. The motors also have protection systems that act on the mechanical and electrical parts of the set.

7 - Installation Requirements

7.1 - Pre-installation

The pre-installation must be guided by the Olsen authorized service to ensure that the environment is adequate the specified requirements for operation of the equipment (*chapter 8 - Installation*) and the positioning in which it is to be used. In this step, the electrical connections must be prepared to supply the equipment.

This equipment is not designed to be installed or operated in a surgical center.

7.2 - Electrical Installation

The power grid must be single-phase, has specific grounding and 10 A/3 mA DR circuit breaker. The circuit breaker must exclusively supply for the Hemodialysis Chair and must be easy and quick to disconnect from the mains. If the power grid presents voltage variation, the installation of a surge protection device is required. The power switch must be turned off at the end of the workday.

There is a table below for sizing the electrical installation:

Tension (V)	Wire Diameter (mm ²)	Distance (m)	Current (A)
118/127/220/230	2,5	Up to 20	10

This equipment should only be connected to a power source with protective grounding. There is a risk of electric shock!

8 - Installation

The Hemodialysis Chair installation must be done by Olsen certified technician consists in assembling upholstered armrests and checking the mains voltage which will supply the equipment and, if necessary, make the tension adjustment. In the case of optional items have been purchased, these should also be installed at this time.

The technician will also verify if the equipment complies with the purchasing order, if it maintains its integrity, and will provide operation, cleaning and maintenance guidance.

When the installation is complete, observe the following items:

- Both right and left bracket locks are installed with the upholstered armrests;
- Armrest handles adjust support properly;
- When the On/Off switch is pressed, the green indicator LED lights up;
- All hand control commands work perfectly;
- The upholsteries of backrest, seat, auxiliary pillow and armrest are intact;
- The castors should be attached to the chair base and be free to rotate.

In the case of optional items have been purchased, check the items below according to the acquired optional item:

- The armrests for blood donation are installed with the accessory brackets;
- The auxiliary stainless-steel tray is installed correctly and makes all adjustments perfectly;

For instructions on the proper functioning of the equipment and its components, read *chapter 5 (Description and Operating of the Equipment)* of this manual.

Installation should only be made by an authorized technician. Installation performed by an unauthorized person will result in loss of warranty.

8.1 - Position of the Equipment

For proper positioning of the equipment, recline the backrest and place the Hemodialysis Chair where it will be used. Perform the up and down commands of the backrest and the seat, making sure there is enough room for the movements to be performed.

8.2 - Olsen Certified Technical Assistance Network

To access the Olsen certified technical assistance network for installation and maintenance contact us by e-mail export3@olsen.odo.br or if you prefer +55 48 2106 6000.

9 - Cleaning and Disinfection

The whole sanitize process must be done with the use of gloves suitable for cleaning and protection, in addition to a mask and protective glasses, according to biosafety standards.

9.1 - Upholsteries and Plastic Covers

The plastic covers and upholsteries must be cleaned with a damp cloth containing neutral soap or detergent only. Olsen advises against using any chemical product to clean these parts, but in the case of disinfection products, it is important to check if it has suitable compatibility and specifications before use on these materials.

Never use hypochlorite or alcohol-based products.

9.2 - Painted Parts

It must be cleaned with a slightly damp cloth containing only soap or mild detergent.

Never use hypochlorite or alcohol-based products.

9.3 - Sterilization by Autoclave

The stainless steel tray can be sterilized by steam autoclaving at the following values (sterilization according to SS-374 resolution of 12/15/1995 and ISO 17665):

- a) 130° C, 2 bar, 15 minutes;
- b) 120° C, 1 bar, 30 minutes;
- c) 134° C, 2,2 bar, 4 minutes.

Note: these items withstand up to 1.000 autoclaving cycles.

All equipment items mentioned in this chapter must be sanitized and sterilized (when appropriate) prior to use.

Olsen is not responsible for defects, deformities, stains or abnormalities caused by improper use of chemical products, contact with tissues, leather, disposable gloves, inks, pigmented detergents and other organic or synthetic products.

10 - Technical Features

Power supply: 118/127/220/230 VC~.

A certified technician when installing the equipment must select the voltage;

Note: all equipment is adjusted in the factory to 220 VC~.

Number of phases: single phase.

Frequency: 50/60 Hz.

Power 118/127 V: 260 VA.

Power 220/230 V: 290 VA.

Protection fuses:

- For 220/230 VC~: 2 fuses of 1,5 AH (5 x 20 mm).
- For 118/127 VC~: 2 fuses of 3 AH (5 x 20 mm).

Power cable specification (in compliance with requirements 6.1 and 6.2 of IEC 60601-1-2:2010):

- Flexible Cable PP Circular 500 V 3 x 1 mm 247-5 NM 53-C5.
- Tri-polar Male Plug 10 A - 250 V (NBR 14136).

Electrical shock protection type: Class I equipment, according to IEC 60601-1-1 and IEC 60601-1-2 standards.

Protection degree: Type B;

Operation mode: non-continuous.

Time On: 30 s.

Time Off: 5 min.

Operating environment conditions:

- **Temperature:** between 15° and 30° C.
- **Relative humidity:** 30% to 70% non-condensing.
- **Pressure:** 75 kPa to 106 kPa.

Harmful water penetration's protection:

- **Equipment:** IPX0.

Thermal protection of the transformer: aperture with 130° C ± 3%.

Lifting capacity (maximum patient weight + accessories): up to 200 kg.

Equipment weight:

- **Net:** 108 kg (with accessories: 122,1 kg).
- **Gross:** 131,5 kg (with accessories: 146,5 kg).

The weight of accessories and optional items:

- Stainless steel tray set: 5,50 kg;
- Armrest for blood donation set: 8,6 kg.

Seat height in relation to the floor:

- Minimum: 550 mm;
- Maximum: 850 mm.

10.1 - Electromagnetic Compatibility (EMC)

The Hemodialysis Chair needs special attention regarding electromagnetic compatibility and must be installed and put into use in accordance with the electromagnetic compatibility information presented in this chapter.

Radio frequency (RF) communication equipment, portable and mobile, can affect the Hemodialysis Chair.

Guidance and Manufacturer's Declaration - Electromagnetic Emissions		
The Hemodialysis Chair is intended for use in the electromagnetic environment specified below. The customer or the user should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic Environment - Guidance
RF emissions CISPR 11	Grupo 1	The Hemodialysis Chair is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings use for domestic purposes.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Non-applicable	
Emissions due to voltage fluctuations/scintillation emissions IEC 61000-3-3	Non-applicable	

Guidance and Manufacturer's Declaration - Electromagnetic Emissions - I			
The Hemodialysis Chair is intended for use in the electromagnetic environment specified below. The customer or the user should assure that it is used in such an environment.			
IMMUNITY Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact	± 6 kV contact	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
	± 8 kV air	± 8 kV air	
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines	± 2 kV for power supply lines	Mains power quality should be that a typical commercial or hospital environment.
	±1 kV for input/output lines	±1 kV for input/output lines	
Surge IEC 61000-4-5	± 1 kV differential mode	± 1 kV differential mode	Mains power quality should be that a typical commercial or hospital environment.
	± 2 kV common mode	± 2 kV common mode	
Voltage dips, short interruptions and voltage variations on the power supply input lines IEC 61000-4-11	<5% U_T (dip > 95% in U_T) for 0.5 cycle	<5% U_T (dip > 95% in U_T) for 0.5 cycle	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Hemodialysis Chair requires continued operation during power mains interruptions, it is recommended that the Hemodialysis Chair be powered from an uninterruptible power supply or a battery.
	40% U_T (60% dip in U_T) for 5 cycles	40% U_T (60% dip in U_T) for 5 cycles	
	70% U_T (30% dip in U_T) for 25 cycles	70% U_T (30% dip in U_T) for 25 cycles	
	<5% U_T (dip > 95% in U_T) for 5 s	<5% U_T (dip > 95% in U_T) for 5 s	
Power frequency (50/60 Hz) Magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE U_T is the a.c. mains voltage prior to applications of the test level.			

Guidance and Manufacture's Declaration - Electromagnetic Immunity

The Hemodialysis Chair is intended for use in the electromagnetic environment specified below. The customer or the user should assure that it is used in such an environment.

IMMUNITY Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Conducted RF IEC 61000-4-6	3 V _{rms} 150 kHz to 80 MHz	3V _{rms}	Portable and mobile RF communications equipment should not be used closer to any part of the Hemodialysis Chair, including cables, than the recommended separation distance calculated by the equation applicable to the transmitter frequency. Recommended separation distance $d = [1,2]^{\frac{1}{2}}\sqrt{P}$ 80 MHz to 800 MHz $d = [1,2]^{\frac{1}{2}}\sqrt{P}$ 800 MHz to 2,5 GHz $d = [2,3]^{\frac{1}{2}}\sqrt{P}$ Where <i>P</i> is the transmitter maximum output power rating in watts (W) according to the transmitter manufacturer and <i>d</i> is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range. ^b Interference may occur in the surroundings of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

a	Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measure filed strength in the location in which the Hemodialysis Chair is used exceeds the applicable RF COMPLIANCE LEVEL above, the Hemodialysis Chair should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Hemodialysis Chair.
b	Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended Separation Distances Between Portable and Mobile RF Communications Equipment and the Hemodialysis Chair

The Hemodialysis Chair is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Hemodialysis Chair as recommended below, according to the maximum output power of the communications equipment.

Rated Maximum Output the Power of the Transmitter (W)	Separation Distance According to Frequency of Transmitter (m)		
	150 kHz to 80 MHz	80 kHz to 800 MHz	800 kHz to 2,5 GHz
	$d = [1,2]^{\frac{1}{2}}\sqrt{P}$	$d = [1,2]^{\frac{1}{2}}\sqrt{P}$	$d = [2,3]^{\frac{1}{2}}\sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,37	0,37	0,74
1	1,2	1,2	2,3
10	3,7	3,7	7,4
100	12	12	23


For transmitters rated at a maximum output power not listed above, the recommended separation distance *d* in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where *P* is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

11 - Dimensional

Measured in Millimeters

12 - Symbology

The symbology complies with standards IEC 60601-1 and IEC 60878.

	Seat Up		Seat Down		Recline Seat
	Backrest Up		Backrest Down		Incline Seat
	Non-Sterile		«Off» (Power)		“On” (Power)
	Caution		General Warning Sign		Warning, Electricity
	Operating Instructions		Refer to the Instruction Manual		General Mandatory Action Sign
	General Prohibition Sign		Stepping Prohibited		Manufacturer
	Type B Applied Part		Serial Number		Alternating Current
	Earth (Ground)		Protective Earth (Ground)		Sterilizable Up to the Temperature Specified
	Humidity Limitation		Temperature Limitation		Fragile, Handle with Care
	Keep Away from Heat		Keep Dry		This Side Up
	Heaping Up				
		Authorized Representative in the European Community			

13 - Important Notes

The reproduction and distribution of these instructions can only be made with prior permission from Olsen Indústria e Comércio S.A.

The technical features of the products described in this manual correspond to the time of its publication. Future technical improvements do not result in any right to update existing products.

The images presented in this manual are illustrative.

This equipment is designed to be free from interference from magnetic fields, external electrical influences, electrostatic discharges, pressure or pressure variation, provided that the equipment is transported, installed, operated and sanitized in accordance with the instructions for use contained in this manual.

13.1 - General Cares - Compulsory Reading

 Follow the instructions in *chapter 7* of this manual (*Installation Requirements*) to suit the electrical network where the equipment will be installed.

 Follow instructions for proper use of the equipment and its accessories as described in *chapter 5* (*Description and Operating of the Equipment*). Improper use could be harmful to equipment, which would not be covered by warranty.

 Follow the instructions in *chapter 9* (*Cleaning and Disinfection*) of this manual for daily cleaning of your equipment.

 Protect your equipment from direct exposure to sunlight. The equipment direct exposure to sunlight may cause premature aging of plastic covers and upholstery.

 Turn off the circuit breaker or unplug the equipment at the end of the working day.

 In case of hand control damage, discontinue use of the equipment, turn it off, and contact Olsen accredited assistance.

 The electric mains cable was developed for use only on the Hemodialysis Chair. The use of this component in other equipment may compromise the emissions and electromagnetic immunity of these.

 Only the authorized technician can replace the equipment's mains cable and internal fuses.

 This equipment is not suitable for use in the presence of a flammable anesthetic mixture with air, O₂ or Nitrous Oxide.

 This equipment should only be operated by physicians and nurses for clinical and outpatient procedures of the areas nephrology, oncology and hematology.

 Use only the mains cable supplied with the equipment for mains connection. Use of cables other than those specified (*chapter 10 - Technical Features*) may result in increased emissions or reduced immunity of the Hemodialysis Chair.

 Do not remove the plastic covers from the equipment. The risk of electric shock. Only the accredited technician is authorized to perform this procedure.

 In the case of damage to the electrical panel, adjacent plastic covers and motors plastic covers disconnect the equipment from the mains and contact the Olsen accredited service. The use of the equipment must be interrupted until maintenance is complete. The use of the equipment in these conditions offers a risk of electric shock!

 Do not perform maintenance or cleaning on the Hemodialysis Chair while it is in use or turned On.

 Do not install or use any electrical equipment over or near the Hemodialysis Chair. If necessary, the Hemodialysis Chair must be checked to see if it is functioning normally in the configuration in which it will be used.

 Do not perform the following procedures if it is possible to touch the patient, even if unintentionally, during the procedure:

- Fuse replacement;
- Do not perform the patient load or other procedures on the Hemodialysis Chair if the castors are not properly locked.

13.2 - Disposal

Debris, residues and infectious materials resulting from the procedures performed on this equipment must be deposited in biological waste duly identified and in accordance with current legislation.

For proper disposal of this equipment and its components and accessories, we recommend that it be sent to specialized recycling companies to ensure the best destination of each component without harm to the environment.

13.3 - Transport and Storage

It is recommended that the equipment's transportation and storage to be made in its original package.

Transport carefully protecting the equipment from falls and impacts.

Protect from moisture, rain exposure and direct contact with liquids.

Keep it sheltered from the sun.

Do not heap more than 5 volumes.

Do not move or store the equipment on uneven surfaces.

Temperature range for transportation and storage: -10° to +45° C.

Moisture limits for transport and storage: 20% to 70%.

13.4 Contraindication

This equipment is contraindicated for any use other than that for which it is intended, or to be operated by unauthorized personnel.

13.5 - Finishing

At the end of the day, switch Off the equipment, unplug it from the mains and clean all upholstered parts, armrests or arms for blood donation and sterilize the stainless-steel tray as described in *chapter 9* of this manual (*Cleaning and Disinfection*).

14 - Troubleshooting

For solving possible problems in a simple and practical way, just follow the instructions in the following table:

Item	Problem	Causes	Solutions
1	The equipment doesn't perform any command.	1° - Equipment is not connected to the electric mains	1° - Connect the equipment to the main electricity
		2° - Electrical circuit breaker is switched Off	2° - Switch on the electrical circuit breaker
		3° - There is no power on the main electricity	3° - Call the Power Supply Company
		4° - Protection fuse is blown	5° - Contact the Olsen authorized service

If you are in doubt or find a problem with the equipment that is not mentioned in this chapter, stop using the equipment immediately and contact your authorized service center or contact us by e-mail export3@olsen.odo.br or by calling +55 48 2106 6000.

15 - Preventive Review

In order to extend your equipment's lifespan, Olsen has prepared a list of the main equipment items for which it recommends performing semi-annual preventive maintenance.

Performing a preventive maintenance by an accredited technician does not interfere with the equipment's warranty period.

STANDARD ITEMS
Check all hand control commands
Check the armrests joints
Check and cleaning of all castors
Spindle drive motors and joints lubrication
Check the left and right bracket locks

OPTIONAL ITEMS
Check the auxiliary stainless-steel tray joints
Check of the armrests for blood donation

Allow only qualified Olsen technicians to perform installation and maintenance on your equipment and accessories.

Use only Olsen original parts and accessories. The use of non-original components can compromise the performance of the equipment, increasing its emissions or reducing its electromagnetic immunity.

Do not make adaptations, modifications or changes to the equipment or its components or accessories.

16 - Warranty Terms

The warranty period for this product is 12 months, counted from the equipment installation date, considering the 90 days legal warranty term, provided that the installation is performed within the period of 90 days from the product's Purchase Invoice issue and the fulfilled of other requirements of this certificate.

- 1** - The maximum term of storage is 3 months from the date of purchase of the product. If the storage period is exceeded, the guarantee is still in progress, even if the product is still stored.
- 2** - Upholstered parts are guaranteed for 6 months.
- 3** - Fuses, cables and transformers are not covered under warranty.
- 4** - The warranty is limited to repair or replacement of defective parts and does not cover defects originated by:
 - a) Non-compliance with the instructions for use and maintenance;
 - b) Falls, crashes and inadequate storage;
 - c) Action of nature agents;
 - d) Damage to upholstery, improper use of chemicals, exposure to unsuitable weather conditions, contact with tissues, leather, disposable gloves, paints, pigmented detergents, razors or any sharp instruments;
 - e) Damage to painted parts and plastic covers caused by improper use of chemicals or by contact with disposable gloves and sharp or piercing objects;
 - f) Connection to wrong voltage power supply.
- 5** - This warranty will cease when:
 - a) In the normal course of its period of warranty;
 - b) Make changes in the product not authorized by Olsen;
 - c) Adulterations in the purchase, installation or services document;
 - d) Installation or technical assistance made by a person not authorized by Olsen;
 - e) Failure to install the equipment for more than 90 days, counted from the purchasing date contained in the invoice;
 - f) By using non-genuine spare parts.
- 6** - The parts repairing or replacement during the warranty period will not extend the original expiry date.
- 7** - The expenses originated from the unit installation, scheduled preventive maintenance, travel, and hotel of the service staff involved in the calls for service for installation or units repair will run under the unit's owner responsibility and in accordance with the distributor's norms.
- 8** - The purchaser, after verifying the services performed in the installation and revision of the equipment, must date and sign the service order provided by the technician and keep along with his invoice of purchase of the equipment, failing which the warranty extension product when necessary.
- 9** - All service requests for warranty equipment must be made with the serial number of the equipment to be serviced and a copy of the purchase or installation document. If this information is not communicated, the service request will be made as not covered by the warranty.


17 - Message from the President

Olsen and its clients:
A successful relationship.

I have linked my name to the factory and to the Dental and Medical equipment which are currently produced and trade in more than 100 countries having in mind the responsibilities and long-term response to this initiative.


Our products are modern, innovative, durable and of low maintenance cost. These characteristics have been achieved thanks to our competent and dedicated team, which make me very proud for many reasons, as they are always giving the best of their creative capacity to our clients.


Our company will always be open to all those who prefer Olsen products, for any necessary information and technical assistance, but especially for comments regarding the relationship with customers. We expect this connection to always bring you satisfaction, resulting in more and more benefits to all of us.

Cesar Olsen

www.olsen.odo.br

+55 48 2106 6000

export3@olsen.odo.br

Olsen

Equipment made to last

Registration at Ministry of Health of Brazilian 10281300012

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