USER MANUAL AMBULATORY STRETCHER

R_x only

Attention: Carefully read this manual
before operating the Olsen's equipment!

Olsen
Equipment made to last

Index

1 - Introduction	5
2 - Package Contents	5
3 - Equipment Presentation	5
3.1 - Standards Items	5
3.2 - Optional Items	5
4 - Parts Identification	6
4.1 - Applied Parts	6
4.2 - Accessories and Detachable Parts	6
5 - Description and Operating of the Equipment	6
5.1 - Turning On the Equipment	7
5.2 - Emergency Battery (Optional)	7
5.2.1 - Battery Charge Warning Sounds	7
5.3 - Movement of the Equipment	7
5.3.1 - Manual Control	7
5.3.2 - Fixation of the Control	8
5.3.3 - Movement Handle	8
5.3.4 - Movements Interruption	8
5.3.5 - Positioning of the Patient, the Operator and Others	8
5.4 - Multi articulated Headrest and Auxiliary Pillow	8
5.5 - Protection Side Rail	9
5.6 - Castors with Locks	9
5.7 - Paper Roll Holder	9
6 - General Features	10
7 - Installation Requirements	10
7.1 - Pre-installation	10
7.2 - Electrical Installation	10
8 - Installation	11
8.1.1 - Equipment Positioning	11
8.2 - Olsen Certified Technical Assistance Network	11
9 - Cleaning and Disinfection	12
9.1 - Upholsteries and Plastic Covers	12
9.2 - Painted Parts	12
10 - Technical Features	13
10.1 - Electromagnetic Compatibility (EMC)	14
11 - Dimensional	16

12 - Symbology.....	17
13 - Important Notes.....	17
13.1 - General Cares - Compulsory Reading	18
13.2 - Disposal.....	19
13.3 - Transport and Storage.....	19
13.4 Contraindication	19
13.5 - Finishing.....	19
14 - Troubleshooting.....	20
15 - Preventive Review	20
16 - Warranty Terms	21
17 - Message from the President.....	22

1 - Introduction

Congratulations on the great choice!

You have acquired an equipment developed for the outpatient environment, designed to accommodate patients during consultations and general medical exams.

This manual is supplemented to the Ambulatory Stretcher Quick Operation Guide provided with your equipment and provides all the information you need to get the most out of your equipment, so read them carefully before using it.

2 - Package Contents

Standards Items:

- 1 Ambulatory Stretcher
- 1 Manual Control
- 2 Protection Side Rail
- 1 Auxiliary Pillow
- 1 Ambulatory Stretcher Quick Operation Guide

Optionals Items:

- 1 Emergency Battery

3 - Equipment Presentation


3.1 - Standards Items

- Automatic Table with 2 Commands
- Manual Control
- Oil-Free Bosch Motors
- Multi articulated Headrest
- Auxiliary Pillow
- 2 Swivel Castors with Individual Locking
- 2 Swivel Castors
- Protection Side Rail
- Paper Roll Holder
- Folding Step
- Movement Handle

3.2 - Optional Items

- Leather Upholstery
- Emergency Battery

4 - Parts Identification

- A - Folding Step
- B - Swivel Castors with Individual Locking
- C - Protection Side Rail
- D - Manual Control
- E - Movement Handle

- F - Multi articulated Headrest
- G - Upholstery of Stretcher
- H - Paper Roll Holder
- I - Electrical Panel
- J - Swivel Castors

4.1 - Applied Parts

The following items of the Ambulatory Stretcher are considered as applied parts to the patient:

- Upholsteries;
- Multi articulated Headrest;
- Auxiliary Pillow;
- Protection Side Rail.

4.2 - Accessories and Detachable Parts

The following items of the Ambulatory Stretcher are considered as accessories or detachable parts:

Detachable Parts:

- Multi articulated Headrest;
- Auxiliary Pillow.

Accessories:

- Emergency Battery.

5 - Description and Operating of the Equipment

The Ambulatory Stretcher is an automated equipment designed for patient accommodation in general medical procedures. With automatic stretcher height movement, it allows the professional to act in an ergonomic way, facilitating patient mobility, avoiding the use of stairs and avoiding sudden impacts.

It is equipped with 1 electrical motor to perform 2 individual movements, controlled by manual control.

The Ambulatory Stretcher features high-strength castors constructed of molded polyurethane that protect the floor and produce extremely low noise levels when being moved.

This equipment should not be used to patient transport.

Before starting use your equipment make sure that it is properly installed, connected to the mains, the circuit breaker is on and the mains is energized.

5.1 - Turning On the Equipment

The electrical panel is in the base frame of the Ambulatory Stretcher, below the footrest.

Being easily accessible, it provides protection fuses (3), On / Off switch (1) and On / Off indicator LED (2).

To turn On the equipment: press the On/Off switch (1) to the "I" position. The LED (2) next to the On/Off switch (1) will turn GREEN, indicating that the equipment is ready for use.

To turn Off the equipment: press the On/Off switch (1) to the "O" position.

The Ambulatory Stretcher has the optional emergency battery (section 5.2), which when available adds a battery charge indicator LED (4).

5.2 - Emergency Battery (Optional)

The emergency battery was developed exclusively for emergency situations where it is necessary to complete a procedure with no electricity on the mains to supply the equipment.

The battery is automatically activated when there is a lack of power, or when the equipment is disconnected from the mains. In this case the equipment emits the "Battery Mode" warning (section 5.2.1).

The battery provides autonomy approximately 3 hours of operation, a period that may vary according to the amount of commands executed by equipment.

The LED (1) indicates the status of the battery according to the colors shown:

- **GREEN:** battery ready for use (full charge).
- **BLUE:** battery in critical condition - Connect the equipment to the mains power supply as soon as possible.

While charging the battery, the blue LED is blinking.

Note: when the BLUE LED On, all commands are blocked. Operation only with the equipment connected to the mains. Upon reaching the maximum battery charge, the LED turns GREEN and unlocks the commands for operation in battery mode again.

5.2.1 - Battery Charge Warning Sounds

In addition to the LED battery charge indicator, the Ambulatory Stretcher uses warning sounds, as shown below:

- **"Battery Mode":** the equipment is disconnected from the mains or there is a lack of mains power.
- **"Low Battery":** the battery charge is less than 50% and must be recharged. While the equipment is not connected to the mains, the warning sound will be repeated every 5 minutes.
- **"Discharged Battery":** the battery charge is at a critical level. While the equipment is not connected to the mains, the warning sound will be repeated every 3 minutes. When this occurs, all the equipment commands are blocked.
- **"Full Battery":** the battery reaches the maximum charge level. The system emits this warning sound only once.


5.3 - Movement of the Equipment

5.3.1 - Manual Control

The Ambulatory Stretcher has 2 commands identified with icons that help the operation in an intuitive way. It is located on the side of the stretcher, on the same side of the step, near the headrest.

To raise the stretcher: set the control switch up (1).

To lower the stretcher: set the control switch down (2).

5.3.2 - Fixation of the Control

The manual control is attached to the Ambulatory Stretcher, but there is the possibility of installing the control on both sides of the stretcher. To do this, ask the accredited technician Olsen for the control to be installed on the side of your preference.

5.3.3 - Movement Handle

In order to adjust the positioning of the Ambulatory Stretcher, the movement handle is produced with steel SAE 1020 mechanical drawn tube, offering safety and stability.

To facilitate transport of the equipment, the movement handle is supplied folded (to down). It is therefore necessary to remove it and fix it in the correct position (to up).

5.3.4 - Movements Interruption

The up and down commands are continuous activation, i.e. they act as long as the control switch is in the selected movement or when it reaches its movement limit (minimum or maximum).

To interrupt these movements, simply position the control switch to the central position, which the equipment immediately stops moving.

5.3.5 - Positioning of the Patient, the Operator and Others

During the equipment movement, the operator and other persons close to him/her should position themselves on the equipment sides, respecting the minimum distance of 50 centimeters, remaining out of the equipment movement area and its components.

To ensure safe and proper patient positioning, the operator should advise him to keep his hands in sight during the entire equipment movement on the equipment upholstery. If the operator sees any situation that may generate any kind of risk to the patient, the stretcher must be stopped immediately.

 **The presence of the operator, persons or objects in the movement areas of the equipment and its components may cause damage to the equipment and / or impair its correct operation.**

 **All equipment's automatic movements must be supervised by the operator. If the operator sees any situation that may generate any kind of risk to the patient, the stretcher must be stopped immediately.**



5.4 - Multi articulated Headrest and Auxiliary Pillow

The headrest (2) can be adjusted manually by its multi articulation controlled by a knob (1).

To adjust the headrest position: turn the knob (1) counterclockwise to loosen it while adjusting. Tighten the knob (1) clockwise to lock in the set position.

To adjust the distance from the headrest to the stretcher: simply pull the headrest (2) away from the stretcher. Do not exceed the 17 cm limit between the headrest and the stretcher (3).

It accompanied by the auxiliary pillow (4) with 2 strips with Velcro® (5) for easy positioning and removal for cleaning.

5.5 - Protection Side Rail

Produced with steel SAE 1020 mechanical drawn tube, it offers safety, avoiding falls in case of lateral displacement of the patient.

Fastened to the equipment through locks also of steel, they have slip with self-lubricating system of high durability and simple operation.

To raise and lock the side rails: lift the railing up against the side of the equipment, and gently press the railing down to lock the sliding device at the end of the course.

To lower the side rails: pull it up to unlock the device, then pull the railing in the opposite direction to the stretcher and bring it down.



5.6 - Castors with Locks

With the Nylon structure and the sideband in Polyurethane, the castors of the equipment have high resistance and perfect symmetry in the displacement of the assembly. Developed to withstand heavier structures and provided with a locking system, they protect the floor and exhibit extremely low noise levels during travel.

The rear castors lock can be activated with the feet.

To lock the castor: press the lever (1) down.

To release the castor lock: push the lever (2) down.


5.7 - Paper Roll Holder

Positioned at the bottom of the stretcher, near the electrical panel, the paper roll holder expedites the preparation of the equipment for the next procedure. Rolls up to 100 mm in diameter can be used.

To place the paper roll on the holder: turn the holder's shaft (1) clockwise, by pulling it in the opposite direction to the stretcher. No need to remove the 2 shaft ends. Just remove one side only to put the disposable sheet roll.

Note: the disposable sheet roll does not accompany the stand.

6 - General Features

- **Upholsteries:** are assembled on a very resistant frame, with density foam 33 and coated with flexible laminated PVC or leather, which besides being resistant, providing easy asepsis for the whole.
- **Mechanical structure:** manufactured in SAE 1020 rolled steel profiles and welded by MIG process, to guarantee the strength and durability.
- **Electrostatic painting:** it is applied to all structural metal parts of the equipment. The special polyurethane paint provides the equipment with high durability coating and has antibacterial properties according to *JIS Z 2801: 2000*, where in 24 hours, the bacterial reduction is greater than 99.9%.
- **Plastic covers:** made of high-strength ABS with acrylic cover, the plastic covers do not require painting and allowing polishing for recovery in case of wear or scratches.
- **Electrical system:** the equipment can work in frequencies of 50 or 60 Hz and can be configured to be connected to the voltages of 118/127/220/230 Volts by an authorized technician. The maximum supply voltage of electronic boards, motors, and other controls is 24 V. The electrical system has On/Off switch and protection fuses.
- **Spindle drive motor:** all motors used in the equipment are produced by Robert Bosch of Brazil and have advantages such as low noise, the absence of oil reservoir, uniformity in displacement, reduction of energy consumption and low cost for maintenance. The motors also have protection systems that act on the mechanical and electrical parts of the set.

7 - Installation Requirements

7.1 - Pre-installation

The pre-installation must be guided by the Olsen authorized service to ensure that the environment is adequate the specified requirements for operation of the equipment (*chapter 8 - Installation*) and the positioning in which it is to be used. In this step the electrical connections must be prepared.

This equipment is not designed to be installed or operated in surgical center.

7.2 - Electrical Installation

The power grid must be single-phase and have specific grounding and 10A/30mA DR circuit breaker. The circuit breaker must exclusively supply the Ambulatory Stretcher and must be easy and quick to disconnect from the mains. If the power grid presents voltage variation, the installation of a surge protection device is required. The power switch must be turned Off at the end of the workday.

There is a table below for sizing the electrical installation.

Tension (V)	Wire Diameter (mm ²)	Distance (m)	Current (A)
118/127/220/230 V	2,5	Up to 20	10

This equipment should only be connected to a power source with protective grounding. There is risk of electric shock!

8 - Installation

The installation of the Ambulatory Stretcher must be done by an Olsen certified technician and consists of positioning the movement handle and auxiliary pillow, as well as checking the mains voltage which will supply the equipment and, if necessary, make the tension adjustment. In case of optional items have been purchased, these should also be installed at this time.

The technician will also verify if the equipment complies with the purchasing order, if it maintains its integrity, and will provide operation, cleaning and maintenance guidance.

When the installation is complete, observe the following items:

- When the On/Off switch is activated, the On/Off indicator LED lights green;
- The headrest adjustment is in perfect working order;
- The 2 commands of the stretcher work correctly;
- The protection side rail are smooth to slide and stay firm when locked;
- The movement handle is fixed perfectly;
- The upholsteries of stretcher, headrest and auxiliary pillow are intact;
- The castors and their locks working properly;

In the case of optional items have been purchased, check the items below according to the acquired optional item:

- The internal emergency battery cable has been installed and the commands operate only on battery power, without connecting the equipment to the mains.

For instructions on the correct operation of the equipment and its components, read *chapter 5 (Description and Operation of the Equipment)* of this manual.

The installation of the equipment must be carried out only by the technical service authorized by Olsen. The installation by unauthorized person will imply the loss of the guarantee.

8.1.1 - Equipment Positioning

When positioning the equipment, perform up and down movements checking for enough space to move the stretcher, perform procedures, and if access to the electrical panel (*section 5.1*) remains free.

8.2 - Olsen Certified Technical Assistance Network

To access the Olsen certified technical assistance network for installation and maintenance contact us by e-mail export3@olsen.odo.br or if you prefer +55 48 2106 6000.

9 - Cleaning and Disinfection

The whole sanitize process must be done with the use of gloves suitable for cleaning and protection, in addition to a mask and protective glasses, according to biosafety standards.

9.1 - Upholsteries and Plastic Covers

The plastic covers and upholsteries must be cleaned with a damp cloth containing neutral soap or detergent only. Olsen advises against using any chemical product to clean these parts, but in the case of disinfection products, it is important to check if it has suitable compatibility and specifications before use on these materials.

Never use hypochlorite or alcohol-based products.

9.2 - Painted Parts

It must be cleaned with a slightly damp cloth containing only soap or mild detergent.

Never use hypochlorite or alcohol-based products.

All equipment item mentioned in this chapter must be sanitized and sterilized (when appropriate) prior to use.

Olsen is not responsible for defects, deformities, stains or abnormalities caused by improper use of chemical products, contact with tissues, leather, disposable gloves, inks, pigmented detergents and other organic or synthetic products.

10 - Technical Features

Power supply: 118/127/220/230 VC~.

A certified technician when installing the equipment must select the voltage;

Note: all equipment is adjusted in the factory to 220 VC~.

Number of phases: single phase.

Frequency: 50/60 Hz.

Power 118/127 V: 260 VA.

Power 220/230 V: 290 VA.

Protection fuses:

- **For 220/230 VC~:** 2 fuses of 1,5 AH (5 x 20 mm).
- **For 118/127 VC~:** 2 fuses of 3 AH (5 x 20 mm).

Power cable specification (in compliance with requirements 6.1 and 6.2 of IEC 60601-1-2:2010):

- Flexible Cable PP Circular 500 V 3 x 1 mm 247-5 NM 53-C5.
- Tri-polar Male Plug 10 A - 250 V (NBR 14136).

Electrical shock protection type: Class I equipment, according to IEC 60601-1-1 and IEC 60601-1-2 standards.

Protection degree: Type B;

Operation mode: non-continuous.

- **Time On:** 30 s.
- **Time Off:** 5 min.

Operating environment conditions:

- **Temperature:** between 15° and 30° C.
- **Relative humidity:** 30% to 70% non-condensing.
- **Pressure:** 75 kPa to 106 kPa.

Harmful water penetration's protection:

- **Equipment:** IPX0.

Thermal protection of the transformer: aperture with 130° C ± 3%.

Lifting capacity (maximum patient weight + accessories): up to 200 kg.

Equipment weight:

- **Net:** 85 kg.
- **Gross:** 113 kg.

Seat height in relation to the floor:

- **Minimum:** 656 mm;
- **Maximum:** 899 mm.

10.1 - Electromagnetic Compatibility (EMC)

The Ambulatory Stretcher needs special attention regarding electromagnetic compatibility and must be installed and put into use in accordance with the electromagnetic compatibility information presented in this chapter.

Radio frequency (RF) communication equipment, portable and mobile, can affect the Ambulatory Stretcher.

Guidance and Manufacturer's Declaration - Electromagnetic Emissions		
The Ambulatory Stretcher is intended for use in the electromagnetic environment specified below. The customer or the user should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic Environment - Guidance
RF emissions CISPR 11	Grupo 1	The Ambulatory Stretcher uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference with nearby electronic equipment. The Ambulatory Stretcher is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings use for domestic purposes.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Non-applicable	
Emissions due to voltage fluctuations/scintillation emissions IEC 61000-3-3	Non-applicable	

Guidance and Manufacturer's Declaration - Electromagnetic Emissions - I			
The Ambulatory Stretcher is intended for use in the electromagnetic environment specified below. The customer or the user should assure that it is used in such an environment.			
IMMUNITY Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact	± 6 kV contact	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
	± 8 kV air	± 8 kV air	
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines	± 2 kV for power supply lines	Mains power quality should be that a typical commercial or hospital environment.
	±1 kV for input/output lines	±1 kV for input/output lines	
Surge IEC 61000-4-5	± 1 kV differential mode	± 1 kV differential mode	Mains power quality should be that a typical commercial or hospital environment.
	± 2 kV common mode	± 2 kV common mode	
Voltage dips, short interruptions and voltage variations on the power supply input lines IEC 61000-4-11	<5% U_T (dip > 95% in U_T) for 0.5 cycle	<5% U_T (dip> 95% in U_T) for 0.5 cycle	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Ambulatory Stretcher requires continued operation during power mains interruptions, it is recommended that the Ambulatory Stretcher be powered from an uninterruptible power supply or a battery.
	40% U_T (60% dip in U_T) for 5 cycles	40% U_T (60% dip in U_T) for 5 cycles	
	70% U_T (30% dip in U_T) for 25 cycles	70% U_T (30% dip in U_T) for 25 cycles	
	<5% U_T (dip > 95% in U_T) for 5 s	<5% U_T (dip> 95% in U_T) for 5 s	
Power frequency (50/60 Hz) Magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE U_T is the a.c. mains voltage prior to applications of the test level.			

Guidance and Manufacture's Declaration - Electromagnetic Immunity

The Ambulatory Stretcher is intended for use in the electromagnetic environment specified below. The customer or the user should assure that it is used in such an environment.

IMMUNITY Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Conducted RF IEC 61000-4-6	3 V _{rms} 150 kHz to 80 MHz	3V _{rms}	<p>Portable and mobile RF communications equipment should not be used closer to any part of the Ambulatory Stretcher, including cables, than the recommended separation distance calculated by the equation applicable to the transmitter frequency.</p> <p>Recommended separation distance</p> $d = [1,2]^{\frac{1}{2}}\sqrt{P}$ <p>80 MHz to 800 MHz</p> $d = [1,2]^{\frac{1}{2}}\sqrt{P}$ <p>800 MHz to 2,5 GHz</p> $d = [2,3]^{\frac{1}{2}}\sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	<p>Where <i>P</i> is the transmitter maximum output power rating in watts (W) according to the transmitter manufacturer and <i>d</i> is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a, should be less than the compliance level in each frequency range. ^b Interference may occur in the surroundings of equipment marked with the following symbol:</p> 

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

a	Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measure filed strength in the location in which the Ambulatory Stretcher is used exceeds the applicable RF COMPLIANCE LEVEL above, the Ambulatory Stretcher should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Ambulatory Stretcher.
b	Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended Separation Distances Between Portable and Mobile RF Communications Equipment and the Ambulatory Stretcher

The Ambulatory Stretcher is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Ambulatory Stretcher as recommended below, according to the maximum output power of the communications equipment.

Rated Maximum Output the Power of the Transmitter (W)	Separation Distance According to Frequency of Transmitter (m)		
	150 kHz to 80 MHz	80 kHz to 800 MHz	800 kHz to 2,5 GHz
	$d = [1,2]^{\frac{1}{2}}\sqrt{P}$	$d = [1,2]^{\frac{1}{2}}\sqrt{P}$	$d = [2,3]^{\frac{1}{2}}\sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,37	0,37	0,74
1	1,2	1,2	2,3
10	3,7	3,7	7,4
100	12	12	23


For transmitters rated at a maximum output power not listed above, the recommended separation distance *d* in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where *P* is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

11 - Dimensional

Measured in Millimeters

12 - Symbology

The symbology complies with standards IEC 60601-1 and IEC 60878.

	Ambulatory Stretcher Up		Ambulatory Stretcher Down		Type B Applied Part
	"OFF" (Power)		"ON" (Power)		Refer to the Instruction Manual
	Manufacturer		Protective Earth (Ground)		Earth (Ground)
	Operating Instructions		General Mandatory Action Sign		General Prohibition Sign
	Caution		General Warning Sign		Warning, Electricity
	Serial Number		Non-Sterile		Heaping Up
	Fragile, Handle with Care		Temperature Limitation		Humidity Limitation
	This Side Up		Keep Dry		Keep Away from Heat
	Stepping Prohibited		Alternating Current		
	Authorized Representative in the European Community				

13 - Important Notes

The reproduction and distribution of these instructions can only be made with prior permission from Olsen Indústria e Comércio S.A.

The technical features of the products described in this manual correspond to the time of its publication. Future technical improvements do not result in any right to update existing products.

The images presented in this manual are illustrative.

This equipment is designed to be free from interference from magnetic fields, external electrical influences, electrostatic discharges, pressure or pressure variation, provided that the equipment is transported, installed, operated and sanitized in accordance with the instructions for use contained in this manual.

13.1 - General Cares - Compulsory Reading

 Follow the instructions in *chapter 7* of this manual (*Installation Requirements*) to suit the electrical network where the equipment will be installed.

 Follow instructions for proper use of the equipment and its accessories as described in *chapter 5* (*Description and Operating of the Equipment*). Improper use could be harmful to equipment, which would not be covered by warranty.

 Follow the instructions in *chapter 9* (*Cleaning and Disinfection*) of this manual for daily cleaning of your equipment.

 Protect your equipment from direct exposure to sunlight. The equipment direct exposure to sunlight may cause premature aging of plastic covers and upholstery.

 Turn off the circuit breaker or unplug the equipment at the end of the working day.

 In case of manual control damage, discontinue use of the equipment, turn it off, and contact Olsen accredited assistance.

 The electric mains cable was developed for use only on the Ambulatory Stretcher. The use of this component in other equipment may compromise the emissions and electromagnetic immunity of these.

 Only the authorized technician can replace the equipment's mains cable and internal fuses.

 This equipment is not suitable for use in the presence of a flammable anesthetic mixture with air, O₂ or Nitrous Oxide. Also not suitable for use in environment rich in Oxygen.

 This equipment should only be operated by physicians and nurses for clinical procedures of ambulatory areas.

 The battery system has a charger integrated into the electrical system of the equipment, which guarantees the good operation and durability of the batteries, avoiding short circuits, oxidations and leaks. However, should any of these defects be found, discontinue use of the equipment, unplug it from the mains and contact the Olsen accredited technician.

 Use only the mains cable supplied with the equipment for mains connection. Use of cables other than those specified (*chapter 10 - Technical Features*) may result in increased emissions or reduced immunity of the Ambulatory Stretcher.

 Do not remove the plastic covers from the equipment. The risk of electric shock. Only the accredited technician is authorized to perform this procedure.

 In the case of damage to the electrical panel, adjacent plastic covers and motors plastic covers disconnect the equipment from the mains and contact the Olsen accredited service. The use of the equipment must be interrupted until maintenance is complete. The use of the equipment in these conditions offers a risk of electric shock!

 Do not perform maintenance or cleaning on the Ambulatory Stretcher while it is in use or turned On.

 Do not install or use any electrical equipment over or near the Ambulatory Stretcher. If necessary, the Ambulatory Stretcher must be checked to see if it is functioning normally in the configuration in which it will be used.

 Do not fuse replacement if it is possible to touch the patient, even if unintentionally.

 Do not perform the patient load or other procedures on the Ambulatory Stretcher if the castors are not properly locked.

13.2 - Disposal

Debris, residues and infectious materials resulting from the procedures performed on this equipment must be deposited in biological waste duly identified and in accordance with current legislation.

For proper disposal of this equipment and its components and accessories, we recommend that it be sent to specialized recycling companies to ensure the best destination of each component without harm to the environment.

13.3 - Transport and Storage

It is recommended that the equipment's transportation and storage to be made in its original package.

Transport carefully protecting the equipment from falls and impacts.

Protect from moisture, rain exposure and direct contact with liquids.

Keep it sheltered from the sun.

Do not heap more than 5 volumes.

Do not move or store the equipment on uneven surfaces.

Temperature range for transportation and storage: -10° to +45° C.

Moisture limits for transport and storage: 20% to 70%.

13.4 Contraindication

This equipment is contraindicated for any use other than that for which it is intended, or to be operated by unauthorized personnel.

13.5 - Finishing

At the end of the day, switch Off the equipment, unplug it from the mains and clean all upholsteries, plastic covers and metal parts, as described in *chapter 9* of this manual (*Cleaning and Disinfection*).

14 - Troubleshooting

For solving possible problems in a simple and practical way, just follow the instructions in the following table:

Item	Problem	Causes	Solutions
1	The equipment doesn't perform any command	1° - Equipment is not connected to the electric mains	1° - Connect the equipment to the main electricity
		2° - Electrical circuit breaker is switched Off	2° - Switch on the electrical circuit breaker
		3° - There is no power on the main electricity	3° - Call the Power Supply Company
		4° - Protection fuse is blown	4° - Contact the Olsen authorized service
2	The device only executes commands when using battery power	1° - Equipment is not connected to the electric mains	1° - Connect the equipment to the main electricity
		2° - Electrical circuit breaker is switched Off	2° - Switch on the electrical circuit breaker
		3° - There is no power on the main electricity	3° - Call the Power Supply Company
		4° - Protection fuse is blown	4° - Contact the Olsen authorized service
3	The machine does not execute commands when using battery power	1° - The indicator light is lit in BLUE color	1° - Recharge the battery until the LED lights GREEN
		2° - The indicator light is lit in GREEN color	2° - Contact the Olsen authorized service

If you are in doubt or find a problem with the equipment that is not mentioned in this chapter, stop using the equipment immediately and contact your authorized service center or contact us by e-mail export3@olsen.odo.br or by calling +55 48 2106 6000.

15 - Preventive Review

In order to extend your equipment's lifespan, Olsen has prepared a list of the main equipment items for which it recommends performing semi-annual preventive maintenance.

Performing a preventive maintenance by an accredited technician does not interfere with the equipment's warranty period.

STANDARDS ITEMS
Check all manual commands
Checking multi articulated headrest adjustments
Checking the sliding fit of the protection side rails
Checking the attachment of the movement handles
Check and cleaning of all castors
Spindle drive motors and joints lubrication

OPTIONAL ITEMS
Check the battery charge

Allow only qualified Olsen technicians to perform installation and maintenance on your equipment and accessories.

Use only Olsen original parts and accessories. The use of non-original components can compromise the performance of the equipment, increasing its emissions or reducing its electromagnetic immunity.

Do not make adaptations, modifications or changes to the equipment or its components or accessories.

16 - Warranty Terms

The warranty period for this product is 12 months, counted from the equipment installation date, considering the 90 days legal warranty term, provided that the installation is performed within the period of 90 days from the product's Purchase Invoice issue and the fulfilled of other requirements of this certificate.

- 1** - The maximum term of storage is 3 months from the date of purchase of the product. If the storage period is exceeded, the guarantee is still in progress, even if the product is still stored.
- 2** - Upholstered parts are guaranteed for 6 months.
- 3** - Fuses, cables and transformers are not covered under warranty.
- 4** - The warranty is limited to repair or replacement of defective parts and does not cover defects originated by:
 - a) Non-compliance with the instructions for use and maintenance;
 - b) Falls, crashes and inadequate storage;
 - c) Action of nature agents;
 - d) Damage to upholstery, improper use of chemicals, exposure to unsuitable weather conditions, contact with tissues, leather, disposable gloves, paints, pigmented detergents, razors or any sharp instruments;
 - e) Damage to painted parts and plastic covers caused by improper use of chemicals or by contact with disposable gloves and sharp or piercing objects;
 - f) Connection to wrong voltage power supply.
- 5** - This warranty will cease when:
 - a) In the normal course of its period of warranty;
 - b) Make changes in the product not authorized by Olsen;
 - c) Adulterations in the purchase, installation or services document;
 - d) Installation or technical assistance made by a person not authorized by Olsen;
 - e) Failure to install the equipment for more than 90 days, counted from the purchasing date contained in the invoice;
 - f) By using non-genuine spare parts.
- 6** - The parts repairing or replacement during the warranty period will not extend the original expiry date.
- 7** - The expenses originated from the unit installation, scheduled preventive maintenance, travel, and hotel of the service staff involved in the calls for service for installation or units repair will run under the unit's owner responsibility and in accordance with the distributor's norms.
- 8** - The purchaser, after verifying the services performed in the installation and revision of the equipment, must date and sign the service order provided by the technician and keep along with his invoice of purchase of the equipment, failing which the warranty extension product when necessary.
- 9** - All service requests for warranty equipment must be made with the serial number of the equipment to be serviced and a copy of the purchase or installation document. If this information is not communicated, the service request will be made as not covered by the warranty.


17 - Message from the President

Olsen and its clients:
A successful relationship.

I have linked my name to the factory and to the dental and medical equipment which are currently produced and trade in more than 100 countries having in mind the responsibilities and long-term response to this initiative.

Our products are modern, innovative, durable and of low maintenance cost. These characteristics have been achieved thanks to our competent and dedicated team, which make me very proud for many reasons, as they are always giving the best of their creative capacity to our clients.

Our company will always be open to all those who prefer Olsen products, for any necessary information and technical assistance, but especially for comments regarding the relationship with customers. We expect this connection to always bring you satisfaction, resulting in more and more benefits to all of us.

Cesar Olsen

www.olsen.odo.br

+55 48 2106 6000

export3@olsen.odo.br

Olsen

Equipment made to last

Registration at Ministry of Health of Brazilian 10281300011

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