# USER MANUAL PORTABLE DENTAL UNIT

Attention: Carefully read this manual  
before operating the Olsen's equipment!

**Olsen**  
Equipment made to last

# Index

1 - Introduction .....	7
2 - Package Content .....	7
3 - Equipment Presentation .....	7
3.1 - Chair .....	7
3.2 - Operating Light .....	7
3.3 - Water Unit .....	7
3.4 - Delivery System .....	7
3.5 - Optional Items .....	7
4 - Parts Identification .....	8
5 - Equipment Assembly .....	9
5.1 - Assembling Instructions .....	9
6 - General Features .....	15
6.1 - Valves .....	15
7 - Installation Requirements .....	15
7.1 - Pre-installation .....	15
7.2 - Electrical Installation .....	16
7.3 - Water for the Water Unit .....	16
7.4 - Water for the Tanks .....	16
7.5 - Sewage System .....	16
8 - Equipment Installation .....	17
9 - Equipment Description and Operating .....	19
9.1 - Before Turning On the Equipment .....	19
9.2 - Turning On the Equipment .....	19
9.3 - Electrical Panel .....	19
9.4 - Chair .....	19
9.4.1 - Multi-articulating Headrest .....	19
9.4.2 - Dental Chair Movements .....	20
9.4.3 - Removable Armrests .....	20
9.4.4 - Emergency Position .....	20
9.5 - Air Compressor .....	20
9.5.1 - Automatic Vapor Drain .....	21
9.5.2 - Air Tank Drainage .....	21
9.6 - Retractable Castors .....	21
9.7 - Tool Tray .....	21
9.8 - Stainless Steel Tray and Straight Arm .....	21

9.9 - Pedals.....	22
9.9.1 - Ejection Pedal .....	22
9.9.2 - Propulsion Pedal .....	22
9.10 - Evolution LED Operating Light .....	22
9.11 - Water Unit.....	22
9.11.1 - Water Tanks.....	22
9.12 - Instruments .....	23
9.12.1 - Couplings for Instruments .....	23
9.12.2 - Venturi Saliva Ejector .....	23
9.12.3 - 3-Way Syringe.....	23
9.12.4 - LED Ultrasonic Scaler (Optional) .....	24
9.12.5 - Ultrasonic scaler use precautions .....	24
9.12.6 - LED Curing Light (Optional) .....	25
9.12.7 - Precautions on the Curing Light Use .....	26
9.13 - Dental Stool.....	26
9.13.1 - Handpieces Kit .....	26
9.14 - Hydraulic Connection Set .....	27
10 - Cleansing and Disinfection .....	27
10.1 - Upholsteries and Plastic Covers .....	27
10.2 - Painted Parts .....	27
10.3 - Saliva Ejectors and Solids Collector .....	27
10.4 - Water Unit and Strainer.....	28
10.5 - Curing Light .....	28
10.6 - Ultrasonic Scaler .....	28
10.7 - Sterilization by Autoclave .....	28
11 - Disassembling.....	29
11.1 - Pre-disassembly Cleaning .....	29
11.2 - Starting Disassembling.....	29
12 - Technical Features .....	37
12.1 - Electromagnetic Compatibility (EMC) .....	38
13 - Dimensional .....	40
14 - Symbology.....	42
15 - Important Notes .....	43
15.1 - General Cares - Compulsory Reading.....	43
15.2 - Disposal .....	44
15.3 - Transport and Storage.....	44

15.4 - Contraindications .....	44
15.4.1 - Curing Light Contraindications .....	44
15.4.2 - Contraindications on the Scaler Use .....	45
15.5 - Daily Care .....	45
16 - Troubleshooting .....	45
17 - Preventive Review .....	46
18 - Warranty Term .....	48
19 - Message from the President .....	49

## 1 - Introduction

Congratulations on the great choice!

You have acquired equipment designed for dental use, built to be assembled and disassembled to facilitate transportation and installation.

The user manual presents the descriptive and technical information, installation instructions, maintenance, cleaning, precautions for use and regulatory guidelines. So before starting the operation of the equipment read carefully its instructions.

## 2 - Package Content

Check out the equipment package contents:

### Standard Items:

- 1 Portable Dental Unit Box
- 1 Portable Dental Unit Chair
- 1 Dental Stool
- 1 Straight Arm
- 1 Articulated Arm of the Operating Light
- 1 Operating Light
- 2 Armrests
- 1 Water Unit
- 1 Stainless Steel Tray
- 1 Castor's Adjustment Wrench
- 2 Water Tanks
- 1 3-Way Syringe's Tip Case
- 1 Kit of Handpieces
- 1 Hydraulic Connection Set
- 1 User Manual Portable Dental Unit

### Optional Items:

- 1 Coupling for Optical Fiber Handpiece
- 1 LED Curing Light
- 1 Ultrasonic Scaler with LED

## 3 - Equipment Presentation

### 3.1 - Chair

- Multi-articulating Headrest
- Backrest and Footrest with Mechanical Adjustment

### 3.2 - Operating Light

- Evolution Operating Light LED 30.000 Lux
- On/Off on the Operating Light Arm
- Side Handles
- Multifaceted Mirror with Polycarbonate Protection Shield

### 3.3 - Water Unit

- Translucent Spittoon Bowl
- Venturi Saliva Ejector of Ø6,5 mm

### 3.4 - Delivery System


- 3-Way Syringe
- Coupling for High-speed Handpiece
- Coupling for Low-speed Handpiece without Cooling

### 3.5 - Optional Items

- Coupling for Optical Fiber Handpiece
- LED Curing Light
- Ultrasonic Scaler with LED
- Porcelain Spittoon Bowl

## 4 - Parts Identification

### A - Water Unit

- 1 - Venturi Saliva Ejector
- 2 - Water Unit Valve
- 3 - Spittoon Bowl Water Spout
- 4 - Cup Filler Water Spout
- 5 - Cup Filler Valve

### B - Operating Light

- 6 - On/Off Operating Light

### C - Air Compressor

- 7 - Air Compressor Pressure Gauge

### D - Dental Stool

- 8 - Dental Stool's Backrest
- 9 - Backrest Adjustment Knob
- 10 - Dental Stool's Seat
- 11 - Dental Stool's Base

### E - Chair

- 12 - Backrest
- 13 - Multi-articulating Headrest
- 14 - Armrests
- 15 - Footrests
- 16 - Seat
- 17 - Progressive Pedal
- 18 - Saliva Ejector Pedal

### F - Delivery System

- 19 - Support of the Straight Arm
- 20 - Straight Arm
- 21 - Stainless Steel Tray
- 22 - Coupling for Low-speed Handpiece
- 23 - Coupling for High-speed Handpiece
- 24 - 3-Way Syringe

## 5 - Equipment Assembly

Before you start assembling the equipment, check if the voltage and frequency of the available electrical power are compatible with the equipment's voltage and frequency.

The place where the equipment will be assembled must have enough space for its assembly and operation (*chapter 13 - Dimensional*), and the surface must be regular and no inclination.


**This equipment has covers and straps for safety and security. When assembling the equipment, remove and retain all covers and straps as they will be required for disassembly and transportation of the equipment.**


**Do not connect this equipment to the mains with a voltage or frequency other than the specified. This will cause damage to equipment that will not be covered by warranty.**


### 5.1 - Assembling Instructions


**1** - Remove the 4 locking clips on the sides of the case and lower the latches (**1**).


**2** - Use the 4 side handles (**2**) of the metal case to lift it. Carefully remove the metal case.

**3** - Loosen the straps (**3**) that secure the dental stool frame and the tray frame to the light's arm, move the light arm away; Move the tray arm away (**4**); loosen the strap (**5**) that is attached to the equipment's frame and the dental stool's frame, then remove the dental stool frame (**6**). Loosen the strap from the packaging of the armrests (**7**) and detach it from the base of the dental stool.

**4** - Remove the tray's arm (**4**) from the operating light's column (**8**), rotate it so that the tray's frame (**9**) stays in the right position, then snap it back into the operating light's column (**8**).


**5** - Release the 2 Velcro® straps (**10**) that secure the base of the dental stool (**11**), loosen the knob (**12**) without removing it and remove the base of the dental stool.

**6** - Insert the dental stool's frame (**6**) into its base (**11**) as shown in the image on the side and tighten the attachment knob (**14**).


**The adjustment screw of the dental stool's backrest (13) is in the water tank's bag. Install the screw before using the dental stool.**

**7** - Release the strap that is attached to the backrest chair (15), remove the 2 side knobs (16); Lift the backrest and attach the right side knob, only to temporarily lock the backrest.


**8** - Install the armrest on the left side (17) by turning it clockwise. Install the left side knob (18). Repeat the process with the right armrest.

**9** - Loosen the knob (19) and raise the footrest (20); Remove the packaging from the pedals (21) and loosen the strap (22) that secures the hoses from the water unit.

**10** - Remove the packaging from the water tanks (23), locate the trim ring (24) and install it on top of the operating light's column (25).


**11** - Loosen the strap (26) and lift the protective case from the operating light's arm (27) to remove it from the holder.

Remove the operating light's arm package from the opposite side of the headlight.


**12** - Connect the operating light's power cable (28) fit the arm into the column (29) being careful not to damage the connector during the fitting.

**13** - Loosen the handle (30) and remove the water unit (31) passing it under the footrest, to install it in its holder on the left side.

**14** - Fit the water unit into the holder under the seat, fitting the jamb (32) and then tighten the handle (33) for attachment.

**15** - Remove the hoses from the packaging, pass them through the finish of the operating light's arm (34) and fasten them with the Velcro® strap. Position the handpieces in the holders (35).

**16** - Install the water spout of the cup filler (36) and spittoon bowl (37); Put water in the tanks (38) and install them by turning them clockwise for fasten.

**17** - Connect the air hose (blue) to the air outlet under the seat (39); Remove the stainless steel tray, instruments and syringe's tip from the package (40); Position the pedals (41) on the right side of the chair.

**18** - Collect all packaging, straps and other materials that come with the equipment and store them with the metal case.

**19** - If your equipment has a curing light (**43**), connect the curing light cable (**42**) to the electrical panel, and place the curing light on its holder in the water unit.

## 6 - General Features

- **Upholsteries:** are assembled on a very resistant frame, covered with soft foam and coated with flexible laminated PVC, it is easy to aseptic of the set;
- **Mechanical structure:** manufactured in SAE 1020 rolled steel profiles and welded by MIG process, giving the equipment resistance and low weight, facilitating transportation, installation, and maintenance.;
- **Electrostatic painting:** It is applied to all structural metal parts of the equipment. The special polyurethane paint provides the equipment with high durability coating and has antibacterial properties according to JIS Z 2801: 2000, in 24 hours, the bacterial reduction is greater than 99.9%.
- **Plastic covers:** made of high-strength ABS with acrylic cover, the fairings do not require painting, which facilitates surface polishing to repair small scratches or wear;
- **Electrical system:** The equipment can work in frequencies of 50 or 60 Hz and can be configured to be connected to the voltages of 127/220 Volts by an authorized technician. The maximum supply voltage of operating light and other electronic resources is 24 V. The electrical system has On/Off switch and protection fuses;

### 6.1 - Valves

The valves used in the equipment's hydro-pneumatic system have been developed to ensure its performance and durability, built-in nonferrous metal with chromium coating.

- **Command valves:** releases spray to the pneumatic handpiece, with a piston system to achieve stagnation and non-retraction of water;
- **Pressure-regulating valves:** they guarantee the propulsion pressure stability of the instruments, while also allowing adjustment to the instrument specifications. The factory setting is 2.2 bar air pressure for high-speed handpieces with a regulated flow rate of 35 liters per minute and 2.4 bar for low-speed handpieces with 55 liters of air per minute.

## 7 - Installation Requirements

### 7.1 - Pre-installation

The pre-installation is the suitability of the environment for the installation of the equipment and should be guided by authorized Olsen service. In order to use all the resources available on the Portable Dental Unit, it is necessary to provide electric power, water, and sewage access. If the place where the equipment will be installed does not have running water or sewage, it is still possible to use the equipment (without using the resources of the water unit) using only filtered drinking water in the tanks and power supply.

**This equipment is not designed to be installed or operated in a surgical center.**

## 7.2 - Electrical Installation

The power grid must be single-phase, specific grounding and 10 A/30 mA DR circuit breaker exclusive to the equipment. If the power grid presents voltage variation, the installation of a surge protection device is required. The power switch must be turned off at the end of the workday.

There is a table below for sizing the electrical installation.

Tension (V)	Wire Diameter	Distance (m)	Current (A)
127/220	2,5	Up to 20	10

**This equipment should only be connected to a power source with protective grounding. There is a risk of electric shock!**

## 7.3 - Water for the Water Unit

For supply to the water unit, the water network must be close to the equipment and easily accessible so that the operator can stop water flow when needed.

The water should have a pressure between 2.8 and 6.0 bar, with recommended pH (Hydrogen Potential) between 6.5 and 8.

It is recommended to use a filter before external equipment water supply to avoid clogging in the internal water system.

**If there is no water network available to supply the water unit, water, and other disinfectant materials must be provided for the sanitization of the spittoon bowl, in case it is to be used (chapter 10 - Cleaning and Disinfection).**

## 7.4 - Water for the Tanks

The water used in the tanks must be potable and filtered. The water tank supplies water to the 3-way-syringe and to the cooling system of the high-speed coupling. This system allows a constant air pressure balance between the handpieces because the water tank maintains the pressure balanced with the equipment's air pressure.


## 7.5 - Sewage System

The sewage system must present good hydraulic declivity (minimum  $-2^\circ$ ) and should preferably be installed underground. The nominal diameter of the tubing must be  $\varnothing$  40 mm.


**The equipment uses the sewage network to drain the Venturi saliva ejector, drain the spittoon bowl, and drain the cup filler. If it is not possible to connect the equipment to a sewage system, the sewer should be installed in a sealed container with a vent, to ensure safe drainage of material from the drain and water unit (tank, ejector Venturi and cup filler).**

## 8 - Equipment Installation


**1** - Locate the sewer connector (1) in the water tank packaging and the following parts:

- Venturi saliva ejector (2);
- Water connector (3);
- Hose clamp (4);
- Pressure clamp (5).

**2** - Identify the hoses coming from the water unit:

- Air of the saliva ejector (6);
- Saliva ejector drainage (7);
- Spittoon bowl drainage (8);
- Water inlet (9).


**3** - Connect the air hose to the Venturi saliva ejector.

**4** - At the sewer connector, identify the open drainage ducts for the Venturi saliva ejector.

**5** - Connect the Venturi valve drain hose to one of the open ducts of the sewer connector.

**6** - Connect the air hose (6) to the Venturi saliva ejector valve. Use the pressure clamp.

**7** - Connect the drainage hose from the spittoon bowl to the sewer connector.


**8** - Apply thread sealing tape to the water connector and install it at the water supply point.


**9** - Connect the water hose to the water connector. Use the hose clamp.

**10** - Install the sewer connector at the sewer access point (saliva ejector hose and drain hoses should already be connected).

**11** - After checking the compatibility of the mains with the equipment, connect the plug to the mains.

**12** - In the absence of a sewage system and/or drinking water supply, two reservoirs may be used:

- **Drinking water reservoir:** must be positioned at a minimum height of 2.5 m/8,2 ft from the equipment base. Its purpose is to supply water to the spittoon bowl and cup filler;
- **Sewage reservoir:** it must be positioned below the level of the equipment base, to ensure the drainage flow of the spittoon bowl and Venturi saliva ejector.

**Note:** Olsen does not provide the external hoses and reservoirs.

## 9 - Equipment Description and Operating

The Portable Dental Unit is equipment intended for dental use, designed to be used in adverse situations and can be easily assembled and disassembled.

It provides individual backrest and footrest movements through mechanical adjustment.

The footrest features clear plastic for upholstery protection.


Its seat is built on the compressor air reservoir that is fixed to the base of the equipment.

### 9.1 - Before Turning On the Equipment

- 1° - Verify that the equipment is correctly assembled according to the instructions in *chapter 5 (Equipment Assembly)* of this manual;
- 2° - Make sure that the equipment is mounted on a regular surface, without inclination and that it does not risk being easily displaced from its current position;
- 3° - Check if the water supply and the drainage installation are in accordance with the instructions of *chapter 8 (Equipment Installation)* of this manual;
- 4° - Check if the equipment is connected to the grid/electric power source compatible with its power supply (*chapter 12 - Technical Features*).

### 9.2 - Turning On the Equipment

- 1° - Press the On/Off switch (A) to the "I" position. The green LED (B) will light indicating that the equipment is switched On and the compressor will start up by filling the compressed air tank.
- 2° - In the water unit (*section 9.1.1*) locate the depressurizing valve (D) and close it;
- 3° - Remove water tanks (E) and make sure they are full. If necessary, fill them with filtered drinking water to the limit indicated on each tank and reconnect them to the platform. If desired place prophylactic products (low concentration);
- 4° - Open the depressurizing valve (D) and check that the water tanks (E) are tightly coupled and without leaks;
- 5° - Before starting to use the equipment, check the operation of all available instruments and controls and clean and sterilize the instruments (*chapter 10 - Cleaning and Disinfection*).

### 9.3 - Electrical Panel

Under the seat, on the left side, near the base of the backrest is the electrical panel, making available the following items:

- A - On/Off Switch
- B - Indicator On/Off LED
- C - Fuse Holders


### 9.4 - Chair

#### 9.4.1 - Multi-articulating Headrest

The headrest (B) is characterized by multiple adjustable positions that ensure more comfort for patients during procedures. The multi-articulating headrest (B) can also be adjusted to care for wheelchair patients (D).

**To adjust the headrest position:** turn the knob (A) counterclockwise to loosen it while adjusting. Turn knob (A) clockwise to lock in the set position.

**To adjust the distance from the headrest to the backrest (C):** pull the headrest (B) away from the backrest. Do not exceed the 13 cm limit between the headboard and the backrest.

## 9.4.2 - Dental Chair Movements

The Portable Dental Unit chair offers backrest and footrest adjustments.

The backrest adjustment mechanism is located under the seat, on the right side. To move it down, press the trigger (A) up while moving the backrest. Release the trigger (A) to lock the backrest in the desired position. To move the backrest up, it is not necessary to use the trigger (A).

The footrest adjustment mechanism (B) is on the left side, under the seat, under the operating light post. To move the footrest down, proceed in the same way as the backrest mechanism (A).

 **To facilitate the adjustment of the backrest, push the backrest upwards while pressing the trigger (A), then lower the backrest.**

 **The footrest adjustment mechanism has a locking knob, which should be used whenever the equipment is disassembled and/or transported. If the knob is tight the mechanism can not be used.**

## 9.4.3 - Removable Armrests

The armrests are designed to rest the patient's arms and can be easily removed, allowing a better approach of the user to perform the procedures.

## 9.4.4 - Emergency Position

Allows brain irrigation by gravity, with slow movements to the maximum negative level ( $-5^{\circ}$  related to the ground). The user only needs to adjust the backrest down to its maximum position.

## 9.5 - Air Compressor


The Portable Dental Unit has a compact 1 HP air compressor (A) incorporated in its structure and controlled by pressure switch (G), to supply the compressed air needs of pneumatic instruments and Venturi ejector.

The air compressor is fixed in the base of equipment and its compressed air reservoir (D) of the equipment is the seat column itself, which has a manometer (C), drain (B) and safety valve (H).

The compressed air outlet has an automatic drain (E) with valve (F) that allows regulating the airflow by turning the valve clockwise to reduce or close the flow and counterclockwise to open.

**Note:** the compressor air pressure is set at the factory according to the instrument's specifications.

 **The reducing airflow in the automatic drain valve may impair the performance of pneumatic and suction instruments.**

### 9.5.1 - Automatic Vapor Drain

Moisture removal filter from the compressed air line installed at the base of the equipment. Performs automatic drainage when inlet pressure is less than 30 PSI.

It has a valve that allows to reduce and close the passage of compressed air to the instruments.

### 9.5.2 - Air Tank Drainage

The air tank drain is a preventive procedure, which avoids the moisture accumulation, the automatic drain clogging and consequently a drop in the performance of pneumatic and ejection instruments.

**To drain the air tank:** with the compressor turned on, open the air tank drain valve by turning it counterclockwise. Leave the hose outlet free to let out the debris/water until the tank is clean.

**Note:** when the equipment remains assembled for a certain period, it is recommended to drain the tank daily and before disassembly.


### 9.6 - Retractable Castors

For easy transport, the base of the Portable Dental Unit has 3 retractable castors. They can be positioned with an L wrench provided in the tool tray (section 9.7).

**To retract the castor (A):** insert the L wrench (C) into the base (B) and turn it clockwise.

**To extend the castor (A):** turn the L wrench (B) counterclockwise.

**Note:** when the base is lifted the castors automatically retract.

### 9.7 - Tool Tray

The tool tray features compartments for storing the portable case's closing clips, the L wrench for moving the castors and additional fuses.

The tool tray is provided under the base of the equipment.

### 9.8 - Stainless Steel Tray and Straight Arm

The stainless-steel tray (C) is seated on a metal frame that has a handle (D) to facilitate its horizontal positioning through the straight arm (B).

The arm allows the positioning of the tray and the instruments, according to the need of the professional.

Next to the stainless-steel tray (C) are the 3-way syringe holder and holders for coupling for high and low-speed handpieces.

**To adjust arm (B) height:** turn knob (A) counterclockwise, position arm (B) as desired, then rotate clockwise to lock into position.


## 9.9 - Pedals

### 9.9.1 - Ejection Pedal

The ejection pedal (A) is responsible for activating the Venturi saliva ejector (B). This drive system provides less air consumption, making the device only work while the pedal (A) is depressed.

**To activate the Venturi saliva ejector (B):** press the pedal (A).

### 9.9.2 - Propulsion Pedal

The propulsion pedal has 2 drive buttons, one for the coupling of the pneumatic micromotor (C) and the other for the coupling for high-speed handpiece (D).

**To activate the micromotor coupling:** press button (C).

**To activate the high-speed coupling:** press button (D).

## 9.10 - Evolution LED Operating Light

The Evolution operating light features LED illumination, a multifaceted mirror for shadow elimination, with a switch on the articulated arm near the headlight (A).


It has articulated arm and side handles allowing easy adjustment of positioning.

The headlight has a polycarbonate protection shield and allows use by LED lighting with single light intensity.

**To turn the operating light On/Off:** press the button on the operating light arm (A).

Evolution LED technical characteristics:

- Light intensity: 30,000 Lux;
- Color temperature: 4500 K.

## 9.11 - Water Unit

The water unit has a spittoon bowl translucent (B) made of polyester resin, with the porcelain spittoon bowl as an optional.

The water unit has a push-button valve-activated cup filler (E), tank water valve (H), air pressurization valve (A) and the Venturi saliva ejector (I) with debris filter.

The spittoon bowl water spout (C) and the cup filler water spout (D) are detachable for an easy cleaning process.

In the water unit are also the water tanks (G), one for the 3-Way syringe and the other for the high-speed handpiece coupling.

To release the water from the spittoon bowl: open the valve (H).


**For water in the cup filler:** press the button (F). Water is released into the cup filler while the button is pressed.

### 9.11.1 - Water Tanks

The use of water tank (G) in the equipment allows the dentist to select the quality of the water used in his instruments.

The water in the tank is intended to supply the coupling for a high-speed handpiece and a 3-way syringe.

**Note:** for the best durability of these components, use good drinking water and replace it daily.

## 9.12 - Instruments

### 9.12.1 - Couplings for Instruments

The couplings (C) are produced with Borden (F) type coupling without cooling for pneumatic micromotor and with cooling (B) for high-speed rotation. These components have exclusive features:

- 1 - All construction in high-quality non-ferrous metal with chrome surface treatment;
- 2 - Disassembly for cleaning without the need for specific tools;
- 3 - The terminal with cooling has a spray flow control ring (E), does not have ledges or jambs, providing greater comfort to the user.

The optical fiber handpiece coupling (G) is optionally supplied with 3,1 V output for handpieces with optical fiber lighting.

**To operate the pneumatic instruments:** remove the instrument from the holder and press the button in the progressive pedal (section 9.9.2). The propulsion air (A) will be released to the turbine coupling, causing it to start up.

**To adjust the water flow of the cooling spray (B):** turn the adjusting ring (E) until the desired water flow is obtained. The green dots are water spray indicators (D).


**To adjust the maximum flow of water:** bring the points closer to each other.

**To decrease the flow of water:** rotate the ring so it distances one point from the other. The further away, the less water in the spray.


**Never use high-speed couplings at the micromotor coupling terminal.**

#### Borden without Cooling

#### Borden

### 9.12.2 - Venturi Saliva Ejector


Developed for saliva ejection, the Venturi ejector has cannula adapters (A and D) for fitting disposable cannulas or autoclavable metal cannulas. It also has a solid debris filter (E) that prevents them from being sent to the sewage system.

**To activate the suction:** remove the ejector (B) from the holder (C) and press the saliva ejector pedal (section 9.9.1).

**To deactivate the suction:** release the saliva ejector pedal (section 9.9.1) and place the ejector (B) in its holder (C).

The ejector has cannula adapters for Ø6.5 mm (A) or Ø9.5 mm (D) and is removable for cleaning.

For proper functioning of this device, it is necessary that the solid debris filter (E) is clean and the drain to the sewage system is properly installed and with proper inclination.

**Note:** the cannulas do not come with the Venturi saliva ejector.

### 9.12.3 - 3-Way Syringe

The 3-way syringe (E) has 3 functions: water jet, air jet, and spray jet.


**To emit air jet:** press the air button (A).

**To emit water jet:** press the water button (B).

**To emit spray jet:** press simultaneously the air (A) and water (B) buttons.

The syringe tip (C) is rotatable and detachable.

**To detach the tip:** press the ring (D) and pull the tip (C), to insert it press the ring (D) and snap it into place.

### 9.12.4 – LED Ultrasonic Scaler (Optional)

The scaler has an ultrasonic generator electronic board that, through the piezoelectric system, provides the tip vibration, giving a high-frequency condition. When the scaler is triggered, it starts the tip vibration, activates the LED and opens the water passage for cooling. Vibration power and water flow control for cooling can be adjusted through individual controls.

Dental applications of this instrument: scarification, root planning, root canal treatment, cavity, and periodontal preparation.

Before starting the scaler operation, install the tip on the transducer and connect the transducer (A) to the coupling (C) checking if are connected.

**To attach the transducer (A) to the coupling (C):** align the adjustment indicators (B) of the two parts and carefully insert them.

**To install the scaler tips:** insert the tip into the transducer (A) by threading it carefully, then insert the Torque wrench (D) into the tip and then rotate it clockwise until it is tight.

**To remove the tip:** turn it counterclockwise with the Torque wrench (D). After installing the tip, check for leakage between the transducer (A) and coupling (C).

**To install a file adapter (not supplied with scaler):** attach the adapter to the transducer (A) and screw it gently. Insert the file into the tip of the adapter, thread the tip and then use the Endo wrench (E) to tighten.




**To remove the adapter:** use the Endo wrench (E) to loosen it and then carefully unscrew.

**To activate the scaler:** push the propulsion/activation pedal (section 9.9.1).

**To adjust the vibration power:** turn the control knob (F) until the desired power is reached.

**To toggle functions:** push the control knob (F) in direction of the table for function General and Perio. Pull the knob (F) gently in the opposite direction of the table for Endo function.

**To adjust the water:** use the control knob (G).

#### Technical Features:

- Frequency: 28 kHz ± 3 kHz
- Power supply: 24 VCC
- Power Rating: 3 W~20 W
- Operation mode: Continuous
- Output force (half-displacement): <2 N
- Tip vibration displacement: ≤100 μm

### 9.12.5 – Ultrasonic scaler use precautions

- Check for vibration with the instrument outside the patient's oral cavity before use. If any abnormalities are found, discontinue use immediately and contact an Olsen accredited service center;
- Wear gloves to handle the Ultrasonic Scaler or its components;
- Use only Torque and Endo wrenches for attachment and removal of the tips. If the tip is not properly attached, it will lose vibration;
- This Ultrasonic instrument has been developed for professional dental use only and should not be used for any other purpose;
- The tip wears out with the use. Tip wear can cause power reduction. If this occurs, replace the tip;
- Do not sharpen or bend the tip. The tips may be damaged and not generate sufficient vibration during scarification;
- When working, the Ultrasonic scaler may affect computers and local area network (LAN) cables. During operation near a radio unit interference may be heard;
- Only use dental autoclave for sterilization of tips, switches, transducer and LED;
- Respect the maximum operating power of the tips used with the Ultrasonic Scaler. Use of power above recommended will damage tip and the scaler itself;

- Keep away from patients with cardiac pacemakers;
- During the Ultrasonic Scaler operation, the tip temperature may rise if spray is not used. Always use enough spray water to cool the tip and dental plane;
- Only use the Ultrasonic Scaler on the dental plane. Contact with skin, gums and mucous membranes may cause injury;
- Keep the Ultrasonic Scaler away from explosive substances and flammable materials;
- Do not subject the handpiece to a strong impact or drop it;
- Do not exceed the power-oriented for each tip. This may damage the dental plan and the tips;
- Do not use the Ultrasonic Scaler on metallic surfaces, ceramic, porcelain or resin prostheses;
- Do not touch or wet the back of the transducer, where electrical connections to the power cord are made. This may result in electric shock;
- Do not force the Ultrasonic Scaler cable when removing it. This may cause disconnection..


### 9.12.6 - LED Curing Light (Optional)

Device for polymerization of materials for restoration and tooth bleaching, curing light produces light radiation to solidify light-sensitive resins, with a short duration shot.

This instrument has a 20 seconds timer for continuous use and protection mode to prevent overheating.

#### Technical Features:

- **Light intensity with optical fiber tip:** 1000 - 1200 mW/cm<sup>2</sup>;
- **Dimensions:** 26 x 25 x 260 mm;
- **Net weight:** 135 g;
- **Wave length:** 420-480 nm;
- **Protection timer:** 200s.

Before you start using the curing light, install the fiber optic tip (A/B) and eyeshield (C).

**To install the tip:** first, attach the light hood (C) to the tip and then attach it to the curing light unit (A/B). The tip must be pushed all the way to the fitting end.

**To activate/deactivate the curing light:** press the On/Off button (E) twice. The green LED (D) indicates that the curing light is operating. The curing light unit will run for 20 seconds and will turn off automatically.

**To toggle the operating mode:** each tap on the selection button (I) switches the operation mode. Press the selection button (I) until the orange LED (F/G/H) indicates the desired mode.

#### The curing light operates in the following modes:

- **Continuous (6):** LED on the maximum power for 20 seconds;
- **Ramp (7):** LED activates at its minimum power, gradually increasing. In 5 seconds reaches its maximum power until the cycle of 20 seconds is finished;
- **Pulse (8):** LED keeps blinking for 20 seconds.

**Activation of protection mode:** If the curing lights triggered for several consecutive times without interval, the protection mode will be activated, locking the device for 20 seconds after each drive. To disable the protection mode, allow the device to stand for 4 minutes. The protection mode can be activated automatically from 9 successive starts.

#### The curing light uses a beep to indicate the following situations:

- 1 - Activation of selection button;
- 2 - Activation of the On/Off button;
- 3 - After operation for 10 seconds;
- 4 - After operation for 20 seconds.

### 9.12.7 - Precautions on the Curing Light Use

- Prohibited the use in patients who present biological reactions;
- Do not aim the light of the curing light directly to the eyes. The curing light produces optical radiation emitted by LED;
- Do not touch the tip directly on the curing material. This will prevent material adhering to the tip, impairing the device performance;
- Do not use the curing light without the light hood (C);
- The optical fiber tip should only be used on the teeth. Avoid touching the patient's gums, lips or skin;
- Use the curing light only on the teeth;
- After 40 seconds of continuous use, the tip end can reach 56° C;
- Consecutive activations of the curing light with an optical fiber tip may cause the tip end to heat, reaching a maximum temperature of 68° C.

### 9.13 - Dental Stool

The dental stool has upholstery with the same characteristics of the upholstery used in the chair.

The gas column (or piston) has high wear resistance and is self-lubricating, allowing it to be assembled and disassembled several times.

Castors and plastic covers are produced in high impact Nylon with high abrasion resistance without suffering abnormalities.

**To adjust the height of the backrest (B):** turn the knob (A) counterclockwise. Position the backrest to the desired height and tighten the knob to lock the backrest in position.

**To adjust the backrest to the seat (D):** turn the knob (C) counterclockwise to release the backrest. Adjust it to the desired position and turn the knob (C) to lock it.

**To adjust the seat height (D):** press the lever (E) up.

**To lower the seat:** sit on the seat (D) until the seat is at the desired height and release the lever (E).

**To raise the seat:** get up from the stool and release the lever (E) as soon as the seat (D) is at the desired height.

#### 9.13.1 - Handpieces Kit

The Portable Dental Unit receives as a series item the PRO PAX SU B2 kit, with the following items:

- 1 Pneumatic Micromotor EX 203C B2;
- 1 Contra Angle NAC EC;
- 1 Straight Piece EX 6B;
- 1 High-speed Pana-Max Standard Push-Button Pen, Single Spray.

Ergonomic instruments made of stainless steel, with ceramic bearings, offer high strength and durability with low noise levels (approximately 62 dB).

## 9.14 - Hydraulic Connection Set

The connection set is to facilitate the installation of the water hoses and the drainage of the water from the spittoon bowl, cup filler, and Venturi saliva ejector, with the drain connector, hose clamps, water connector and for Venturi saliva ejector.

## 10 - Cleansing and Disinfection

**The whole sanitize process must be done with the use of gloves suitable for cleaning and protection, in addition to a mask and protective glasses, according to biosafety standards.**

### 10.1 - Upholsteries and Plastic Covers

The plastic covers and upholsteries must be cleaned with a damp cloth containing neutral soap or detergent only. Olsen advises against using any chemical product to clean these parts, but in the case of disinfection products, it is important to check if it has suitable compatibility and specifications before use on these materials.

To clean the operating light polycarbonate protection shield, use only a soft cloth dampened with liquid glycerin.

**Never use hypochlorite or alcohol-based products.**

### 10.2 - Painted Parts

It must be cleaned with a slightly damp cloth containing only soap or mild detergent.

**Never use hypochlorite or alcohol-based products.**

### 10.3 - Saliva Ejectors and Solids Collector


Daily disinfect the saliva ejector's hoses with an appropriate disinfection product for PVC hoses. It is essential to read the instructions of the asepsis product to avoid misuse or overdosage, which can cause damage to the medium and long term to the hoses. Using a proper concentration of the cleaning product, suck up with the suction needed for the effectiveness of the process. With the product still in the hose, place the suction on the hose holder. After the necessary action period of the product, suck 1 L of water.

The solids collector must also be cleaned daily. The saliva ejectors' efficiency may be impaired if this filter is clogged. In the case of saliva ejector's performance reduction, clean its filters.

To clean the filters, follow the instructions below:

- 1° - Undo the solids collector cover (A);
- 2° - Remove the solids collector (B) for cleaning;
- 3° - After cleaning, re-assemble the solids collector.

To clean the water unit drain, pour aseptic liquid observing the drain in order to remove any debris from the inside of the hose.

## 10.4 - Water Unit and Strainer

The water spouts from both spittoon bowl (B) and cup filler (A) are removable to facilitate removal from spittoon bowl (E) for cleaning. The disinfection of the spouts should be done with 70% alcohol.

For spittoon bowl cleaning (E) remove spittoon bowl spout (B), spittoon bowl finish (C), strainer (D) and then remove spittoon bowl (E) for cleaning.

For removal of the strainer (D) use a tweezer or gloves to avoid direct contact with the waste.

With a soft sponge, running water and mild soap or detergent clean the strainer (D) and spittoon bowl (E). Do not use an abrasive sponge.

**All debris and contaminated materials should be disposed of in the biological waste.**

## 10.5 - Curing Light

To clean and disinfect the curing light tips, use gauze or disposable wipe moistened with soap and neutral detergent. Do not use alcohol, strong or abrasive alkaline detergents, detergents based on bleach, acetone, or other germicides.

The whitening tip cannot be autoclaved.

Curing light's body cleaning must be done with neutral detergent or 70% alcohol.

## 10.6 – Ultrasonic Scaler

The transducer (E), the LED (C), the tips and the scaler wrenches must be autoclaved for sterilization.


Prior to the transducer (E) sterilization, remove the tip of the scaler (if it is attached), LED removable cover (A), the light driver (B) and the finishing ring (D). These items should be sterilized with 70% alcohol.

For the proper functioning of the scaler, the LED (C) must be correctly inserted in the transducer (E) matching the LED's positive pole with the transducer's positive pole.

Only use a dental autoclave for sterilization of tips, wrenches, transducer (E) and LED (C).

**Following sterilization methods are forbidden:**

- Put in boiling water;
- Heat in stove, oven or microwave;
- Soak in disinfectants such as iodine, alcohol or glutaraldehyde.

## 10.7 - Sterilization by Autoclave

Prior to autoclaving, clean the items, removing any organic residue, both from the surface and internal ducts (if any). Then carefully dry each item, including the internal ducts, if possible, with compressed air.

Individually pack each item, with its own sterilized packaging for the autoclaving process. For steam autoclaving, the following values may be used (sterilization in accordance with ISO 17665):

- 130° C, 2 bars, 15 min;
- 120° C, 1 bar, 30 min;
- 134° C, 2,2 bars, 4 min.

**The following equipment's items are autoclavable:**

- Stainless steel trays;
- 3-way syringe's tip;
- Scaler Torque and Endo wrench;
- Scaler LED and transducer;
- Optical fiber tip of the curing light;

**Note:** the scaler transducer with LED cannot be autoclaved in contact with other materials.

To sterilize the scaler tips, set the autoclave to 132° C for 3 to 6 minutes and pressure of 30 PSI (2 bar) or as indicated on ISO 17665 standard.

**Regarding the resistance to the autoclaving process, the following items support:**

- Up to 1000 cycles: scaler tip, stainless steel tray, 3-way syringe tip and optical fiber curing light tip;
- Up to 600 cycles: transducer and LED of the scaler.

**All equipment items referenced in this chapter must be sanitized and sterilized (where applicable) prior to use.**

**Do not use any type of oil on the items to autoclave.**

**Olsen is not responsible for defects, deformities, spots on or abnormalities caused by improper use of chemical products, contact with tissues, leather, disposable gloves, inks, pigmented detergents and other organic or synthetic products.**

## **11 - Disassembling**

The process of disassembling the Portable Dental Unit equipment requires attention to detail so that the equipment can fit inside its metal case and not be damaged during transportation.

Before disassembling the equipment, it must be properly sanitized (*chapter 10 - Cleaning and Disinfection*), to avoid deterioration of the hoses and contamination of the equipment.

Empty and clean the equipment water tank separately to avoid cross-contamination.

Clean and drain the hoses to prevent leakage in the equipment from inside.


### **11.1 - Pre-disassembly Cleaning**

- 1° - Close the air valve of the water unit, remove the 2 water tanks and empty them;
- 2° - Reinstall the empty tanks and open the air valve;
- 3° - Press the water button on the 3-way-syringe until no more water comes out;
- 4° - Release the maximum flow of water in the high-speed coupling and actuate until no more water comes out;
- 5° - Remove the instruments and the syringe tip, sterilize them and store them in their packaging;
- 6° - Clean the saliva ejector, filter, spittoon bowl, spouts, and strainer according to the instructions in *chapter 10 (Cleansing and Disinfection)* of this manual;
- 7° - Drain the water unit drain completely and disconnect the drain or sewer connector. It is not necessary to disconnect the hoses from the sewer connector.
- 8° - Open the air compressor drain and drain all compressed air.

### **11.2 - Starting Disassembling**

After cleaning the equipment, disconnect the equipment from the power supply and then identify and separate all the packages and straps with Velcro® to follow the process of disassembly

- A - 1 Operating Light arm packing
- B - 1 Stainless steel tray and instrument packing
- C - 1 Packing of water tanks
- D - 2 Packaging 23 x 23 cm - pedals/hoses
- E - 1 Packing of Armrests

- F - 1 Strap 50 cm
- G - 2 Straps 60 cm
- H - 1 Strap 105 cm
- I - 4 Straps 40 cm

**1** - Loosen the Velcro® straps that hold the instruments hoses and wrap them close to the water unit, forming a circumference of approximately 20 cm; use the 23 x 23 package to store it in the water unit holder between the unit and the seat.

**2** - Disconnect the compressed air supply hose from the equipment.

**3** - Keep the 2 pedals in the other 23 x 23 package; place them in front of the automatic drain gauge.

**4** - Remove the water spouts from the spittoon bowl and the the cup filler and remove the tanks of the water unit, store them in the water tank package.

**5** - To remove the water unit, Release the water unit locking handle (1), store it with the hoses and instruments in the package; If the equipmente has a curing light and/or ultrasonic scaler, disconnect it and store it in the stainless steel tray and instrument package.

**6** - Pass the water unit under the seat and attach it to the holder (2) on the right side of the unit; lock the handle (3).

**7** - Attach the corrugated hose to the strap (4) attached to the base of chair.

**8** - Lower the footrest to the limit and lock the handle (8).

**9** - Remove the handle from the left side of the backrest (9), position the metal rod (10) in parallel with the armrest and rotate the armrest counterclockwise and remove it; repeat the operation with the right armrest.

**10** – Store the armrests in the armrests package (11); store the tray, instrument kit, and syringe tip on the packaging of the tray and instruments.

**11** - Lower the backrest, reposition the side handles (13) and attach the Velcro® to the packaging of the tray and instruments (12) behind the backrest as indicated above.


**12** - Lift the operating light's arm to remove it carefully; disconnect the internal cable from the headlight, store the operating light's arm in the package up to the headlight.

**13** - Attach the operating light's arm to the bracket base support by snapping it all the way up into the post light. (21); keep the finish (22) near the water tanks.

**14** - Remove the straight arm from its socket on the post light and reinsert it upside down; position the straight arm on the side of the chair with the tray structure (23) on the backrest.

**15** - Adjust the dental stool seat to the minimum height and its backrest to the maximum height, loosen the seat knob (24) and detach its base from the frame; fit the base into the holder under the seat of the chair and fasten the knob (25).

**16** - With the 40 cm strap secure the dental stool base on the seat adjustment lever (26) and make sure that the hoses and the electric cable will not be loose.

**17** - Fit the dental stool backrest under the knob (28), and on the straight arm (29); position the tray frame (30) on the dental stool seat; Remove the knob from the backrest of of the dental stool (31) and store together with the water tanks.


**18** - Position the operating light's arm with the headlight on the backrest chair as in the picture above, attach the dental stool seat with the 105 cm strap (32) and pass the strap through the opening in the seat chair frame (33). Attach the tray frame to the operating light's arm with the 40 cm strap (34).

**19** - Attach the armrests package in front of the base of the dental stool, with the Velcro® strap (35) behind the base.

**20** - Attach with Velcro® strap the packaging of water tanks on the operating light's arm.

**21** - To move the equipment, drive the castors in the base (37) according to section 9.6 of this manual.

**22** - Position the metal case on the equipment, save the user manual and lock the clamps. Use the side handles (38) to carry the equipment.

## 12.1 - Electromagnetic Compatibility (EMC)

The Portable Dental Unit needs special attention regarding electromagnetic compatibility and must be installed and put into use in accordance with the electromagnetic compatibility information presented in this chapter.

Radio Frequency (RF) communication equipment, portable and mobile, can affect the Portable Dental Unit.

Guidance and Manufacturer's Declaration - Electromagnetic Emissions		
The Portable Dental Unit is intended for use in the electromagnetic environment specified below. The customer or the user of the Portable Dental Unit should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group1	The Portable Dental Unit is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Not Applicable	
Emissions due to voltage fluctuations/scintillation emissions IEC 61000-3-3	Not Applicable	

Guidance and Manufacturer's Declaration - Electromagnetic Immunity- I			
The Portable Dental Unit is intended for use in the electromagnetic environment specified below. The customer or the user of the Portable Dental Unit should assure that it is used in such an environment.			
Immunity test	IEC 60601 Test Level	Compliance level	Electromagnetic environment guidance
Electrostatic Discharge (DES) IEC 61000-4-2	± 6 kV contact ± 8 kV air	Class A	Floors should be wood, concrete or ceramic. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	Class A	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV differential Mode ± 2 kV common mode	Class A	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	< 5% $U_T$ (dip > 95% in $U_T$ ) for 0.5 cycle 40% $U_T$ (60% dip in $U_T$ ) for 5 cycles 70% $U_T$ (30% dip in $U_T$ ) for 25 cycles < 5% $U_T$ (dip > 95% in $U_T$ ) for 5 s	Class A Class A Class A Class B	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Portable Dental Unit requires continued operation during power mains interruptions, it is recommended that the Portable Dental Unit be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) Magnetic field IEC 61000-4-8	3 A/m	Class A	Power frequency magnetic fields should be at levels characteristic of a typical commercial or hospital environment.
NOTE: $U_T$ is the supply voltage c.a. before the application of the test level.			

### Guidance and Manufacture's Declaration - Electromagnetic Immunity

The Portable Dental Unit is intended for use in the electromagnetic environment specified below. The customer or the user of the Portable Dental Unit should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 V <sub>rms</sub> 150 kHz to 80 MHz	3 V <sub>rms</sub>	Portable and mobile RF communications equipment should not be used closer to any part of the Portable Dental Unit, including cables, than the recommended separation distance calculated by the equation applicable to the transmitter frequency. <b>Recommended separation distance</b> $d = [1,2]^{\frac{2}{3}}\sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = [1,2]^{\frac{2}{3}}\sqrt{P} \quad 800 \text{ MHz to } 2,5 \text{ GHz}$ $d = [2,3]^{\frac{2}{3}}\sqrt{P}$ where <i>P</i> is the transmitter maximum output power rating in watts (W) according to the transmitter manufacturer and <i>d</i> is the recommended separation distance in meters (m).  Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey <sup>a</sup> , should be less than the compliance level in each frequency range <sup>b</sup> . Interference may occur in the surroundings of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

- a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM, and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measure filed strength in the location in which the Portable Dental Unit units are used exceeds the applicable RF COMPLIANCE LEVEL above, the Portable Dental Unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Portable Dental Unit.
- b

Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

### Recommended separation distances between portable and mobile RF communications equipment and the Portable Dental Unit

The Portable Dental Unit is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Portable Dental Unit can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Portable Dental Unit as recommended below, according to the maximum output power of the communications equipment




The rated maximum output power of the transmitter (W)	Separation distance according to the frequency of the transmitter (m)		
	150 kHz to 80 MHz $d = [1,2]^{\frac{2}{3}}\sqrt{P}$	80 MHz to 800 MHz $d = [1,2]^{\frac{2}{3}}\sqrt{P}$	800 MHz to 2,5 GHz $d = [2,3]^{\frac{2}{3}}\sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,37	0,37	0,74
1	1,2	1,2	2,3
10	3,7	3,7	7,4
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance *d* in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where *P* is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 - At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 - These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

Measured in millimeters

## 14 - Symbology

The symbology complies with IEC 60601-1 and IEC 60878.

	Backrest Up		Backrest Back		Footrest Down
	Footrest Up		Bowl Flush		Cup Filler
	Multifunction Syringe (Air-Water)		Dental Air Motor		Turbine
	Turbine with Illumination		Ultrasonic Scaler		Curing Light
	Spray-cooling		Saliva Ejector		Suction Handpiece
	Variability, for Rotating Movement		Dental Operating Light		Hand Control Valve
	Foot Switch		"OFF" (Power)		"ON" (Power)
	Non-Sterile		Caution		General Warning Sign
	Warning, Electricity		General Mandatory Action Sign		Operating Instructions
	Refer to the Instruction Manual		General Prohibition Sign		Stepping Prohibited
	Manufacturer		Level		Type B Applied Part
	Type BF Applied Part		Serial Number		Sterilizable Up to the Temperature Specified
	Protective Earth (Ground)		Earth (Ground)		Keep Away from Heat
	Humidity Limitation		Fragile, Handle with Care		Temperature Limitation
	Keep Dry		Heaping Up		This Side Up
	Alternating Current	 Authorized Representative in the European Community			

## 15 - Important Notes

The reproduction and distribution of these instructions can only be made with prior permission from Olsen Indústria e Comércio S.A.

The technical features of the products described in this manual correspond to the time of its publication. Future technical improvements do not result in any right to update existing products.

The images presented in this manual are illustrative, and the actual equipment may vary in color, size, and shape compared to the illustrations presented here.

This equipment is designed to be free from interference from magnetic fields, external electrical influences, electrostatic discharges, pressure or pressure variation if the equipment is transported, installed, operated and sanitized in accordance with the instructions for use contained in this manual.

### 15.1 - General Cares - Compulsory Reading

 Follow the instructions in *chapter 7* of this manual (*Installation Requirements*) to suit the electrical and hydraulic network where the equipment will be installed.

 Follow instructions for proper use of the equipment and its accessories as described in *chapter 5* of this manual (*Equipment Assembling*). Improper use could be harmful to equipment, which would not be covered by warranty.

 Clean the equipment according to the instructions in *chapter 10* (*Cleansing and Disinfection*) of this manual.

 Protect your equipment from direct exposure to sunlight. Direct exposure of the equipment to sunlight may cause premature aging of plastic covers and upholstery.

 Turn off the circuit breaker or disconnect equipment from the mains and close the office water supply at the end of work.

 This equipment does not have a battery pack.

 In case of damage to the pedals, discontinue the use of the equipment, switch it off, and contact an Olsen accredited service.

 The mains cable and the scaler transducer are designed for exclusive use on the Portable Dental Unit. The use of these components in other equipment might compromise their emissions and electromagnetic immunity.

 Use only the scaler cable and transducer supplied with the equipment. The use of a cable or transducer other than those specified (*chapter 12 - Technical Features*) may result in increased emissions or reduced electromagnetic immunity of the Portable Dental Unit.

 Only the authorized technician can replace this equipment's mains cable and internal fuses.

 This equipment is not suitable for use in the presence of a flammable anesthetic mixture with air, O<sub>2</sub> or Nitrous Oxide.

 This equipment should be operated by dental surgeons for oral and maxillo-facial exams and procedures only.

 Use only the mains cable supplied with the equipment for mains connection. The use of cables other than those specified (*chapter 12 - Technical Features*) may result in increased emissions or reduced immunity of the Portable Dental Unit.

 Do not remove the plastic covers from the equipment. There's a risk of electric shock! Only the accredited technician is authorized to perform this procedure.

In the case of damage to the electrical panel and adjacent plastic covers, disconnect the equipment from the mains and contact the Olsen accredited service. Interrupt the equipment use until maintenance is complete. The equipment used in these conditions offers a risk of electric shock!

Do not perform maintenance or cleaning on the Portable Dental Unit while it is in use or turned On.

Do not install or use any electrical equipment over or near the Portable Dental Unit. If necessary, the dental chair must be checked to see if it is functioning normally in the configuration in which it will be used.

Do not perform the following procedures if it is possible to touch the patient, even if unintentionally, during the procedure:

- Coupling or removal of low and high-speed handpiece from its coupling;
- Coupling or removal of the scaler from its coupling;
- Fuse replacement;
- Do not perform patient boarding or other procedures on the Portable Dental Unit if the castors are not retracted.

## 15.2 - Disposal

Debris, residues and infectious materials resulting from the procedures performed on this equipment must be deposited in biological waste duly identified and in accordance with current legislation.

For proper disposal of this equipment and its components and accessories, we recommend that it be sent to specialized recycling companies to ensure the best destination of each component without harm to the environment.

The draining of this equipment should not be disposed of in a common sewage system, complying with the requirements established by the health agencies.

## 15.3 - Transport and Storage

For transportation and storage, the equipment must be disassembled and protected inside its original metal case.

Transport carefully protecting equipment from falls and impacts.

Protect from moisture, rain exposure and direct contact with liquids.

Keep it sheltered from the sun.

Do not heap up more than 4 volumes.

Temperature range for transportation and storage: -10 °C to +45° C.

Moisture limits for transport and storage: 20% to 70%.

Do not move or store the equipment on uneven surfaces.

## 15.4 - Contraindications

This equipment is contraindicated for any use other than that for which it is intended or to be operated by unqualified personnel.

### 15.4.1 - Curing Light Contraindications

**The curing light use in cardiac, pregnant or child patients should be performed with precautions;**

### 15.4.2 - Contraindications on the Scaler Use

- It's forbidden the scaler use of hemophilic patients.
- Patients, dentists or dental assistants with pacemakers are prohibited from using or approaching the scaler during its use.
- The scaler uses in cardiac, pregnant or child patients should be performed with precautions.

## 15.5 - Daily Care

At the end of the workday, turn off the equipment and follow these instructions:

- Replace the water in the water unit tanks daily;
- Clean the equipment daily according to the instructions in *chapter 10 (Cleansing and Disinfection)* of this manual;
- Close the water supply register and turn off the equipment's power supply circuit breaker;
- Drain the compressor as instructed in *section 9.5 (Air Compressor)* of this manual.

## 16 - Troubleshooting

For solving possible problems in a simple and practical way, just follow the instructions in the following tables:

Item	Problem	Causes	Solutions
1	Operating light is not lighting	1°-The equipment is not connected to the mains	1°- Connect the equipment to the mains
		2°- Electrical circuit breaker is switched Off	2°- Switch On the electrical circuit breaker
		3°-There is no power on the mains	3°- Call the power supply company
		4°-LED is blown	4°-Contact Olsen certified assistance
		5°-The protection fuse is blown	5°-Contact Olsen certified assistance
2	Pneumatic handpiece (micro-motor/turbine) doesn't work or is weak	1°-The handpiece ducts are clogged	1°-The handpiece ducts are clogged
		2°-There's looseness on the handpiece coupling	2°-Attach the handpiece correctly
		3°-The air valve is not completely open	3°-Open completely the air valve
		4°-There's insufficient air pressure for the equipment	4°-Check the inlet air pressure
		5°-The air compressor is not working properly	5°-Call the air compressor assistance
		6°-There's blocking in the pneumatic system	6°-Contact Olsen certified assistance
3	The pneumatic handpiece has water leakage on the coupling	1°-There's looseness on the handpiece coupling	1°-Attach the handpiece correctly
		2°-The gasket is worn	2°-Replace the gasket
		3°-The gasket is not properly sealing the coupling	3°-Apply original handpiece seal
		4°-The coupling is worn	4°-Contact Olsen certified assistance
4	The pneumatic handpiece doesn't have water on the spray	1°-The coupling spray's adjustment is closed	1°-Align the coupling's green dots
		2°-The water tank is empty	2°-Fill the water tank
		3°- There's looseness on the handpiece coupling	3°-Attach the handpiece correctly
		4°-The air valve is not completely open	4°-Open the air valve
		5°-There's insufficient air pressure for the equipment	5°-Open the main air valve
		6°-The hydro-pneumatic system is locked	6°-Contact Olsen certified assistance

5	The saliva ejector is weak or loses suction during the procedure	1°-The saliva ejector's filter is clogged	1°-Clean the saliva ejector's filter
		2°-There's insufficient air pressure for the equipment	2°-Open the main air valve
		3°-The equipment's drain hose is obstructed	3°-Check if the sewer's drain pipes are crushed and release the flow
		4°-The sewer is clogged	4°-Check if there's water outflow on the spittoon bowl and clean the sewer if necessary
		5°-The hydro-pneumatic system is locked.	5°-Contact Olsen certified assistance
6	The curing light isn't working	1°-There's a problem with power supply	1°-Check Troubleshooting, item 1
		2° - Block of 10 consecutive activations	2°-Wait at least for 20 seconds to activate the curing light again
		3°-Possible overheating	3°-Contact Olsen certified assistance
7	The scaler vibrates little or nothing	1° - There's a problem with the power supply.	1° - Check Troubleshooting, item 1.
		2° - The tip is not properly attached to the transducer.	2° - Remove and attach the tip again using the proper wrench.
		3° - The tip is worn, or the scaler has a defect.	3° - Contact Olsen certified assistance.
8	The scaler is overheating	1° - The vibration control is not adjusted according to tip specification.	1° - Adjust the scaler power control according to the tip specification.
		2° - The tip is not properly attached to the transducer.	2° - Remove and attach the tip again using the proper wrench.
		3° - The tip is worn, or the scaler has a defect.	3° - Contact Olsen certified assistance.
9	The scaler has little or no water	1° - The water tank is empty.	1° - Fill the water tank.
		2° - The scaler water valve is closed.	2° - Open scaler water valve.
		3° - The tip has a hydro-pneumatic system.	3° - Contact Olsen certified assistance.
10	The turbine's optical fiber light doesn't work	1° - The equipment is not performing the movement.	1° - Check Troubleshooting, item 1.
		2° - The optical fiber's light is switched Off.	2° - Turn On the optical fiber's light switch.
		3° - The handpiece is not properly attached.	3° - Reattach the handpiece until the joint gasket is tighter.
		4° - There's a pneumatic or electronic problem.	4° - Contact Olsen certified assistance.

If you are in doubt or find a problem with the equipment that is not mentioned in this chapter, stop using the equipment immediately and contact your authorized service center or contact us by e-mail [export3@olsen.odo.br](mailto:export3@olsen.odo.br) or by calling +55 48 2106 6000.

## 17 - Preventive Review

In order to extend your equipment's lifespan, Olsen has created the extended warranty system by performing scheduled service maintenance.

During the service maintenance, the technician will evaluate the equipment's general maintenance condition, as well as the monitoring of component wear and if there is a need for lubrication.

The technician may suggest the replacement of the parts with wear and tear and will provide guidance on the daily care necessary for the proper equipment functioning.

The following tables list the items that must be checked by the technician:

DENTAL CHAIR/PEDALS
Check backrest adjustment device
Check footrest adjustment device
check headboard movement
Check suction pedal valve
Check the valves of the progressive pedal
Check automatic vapor drain

DELIVERY SYSTEM/ STAINLESS STEEL TRAY
Check the straight arm clamp
Check the pressure and coupling rings
Check the valves, vane, and holders
Check and lubricate syringe buttons
Check instrument hoses
Disassembly and lubrication of syringe buttons
Insertion wear testing (scaler)

WATER UNIT
Check the Venturi ejector suction
Lubricate of debris separator filter rings
Spittoon bowl ring disassembly, cleaning, and lubrication

OPERATING LIGHT
Check of joint movements
Check bulkhead and mirror
Check focus and intensity
Check the LED

DENTAL STOOL
Check the castors
Check piston and movement

CURING LIGHT
Check the intensity and operating programs
Check tip and support ring

BASE/CONNECTIONS
Check the electrical, water and air supply connections
Check the sewer connection (damper)
Check the castors

**Olsen recommends that the items in this chapter be reviewed every 180 days to prevent possible equipment failure or performance even after the equipment warranty expires.**

**The performance of preventive or corrective maintenance by an accredited technician does not interfere with the warranty period.**

**Allow only qualified Olsen technicians to perform installation and maintenance on your equipment and accessories.**

**Use only Olsen's original parts and accessories. The use of non-original components can compromise the performance of the equipment, increasing its emissions or reducing its electromagnetic immunity.**

**Do not make adaptations, modifications or changes to the equipment or its components or accessories.**

## **18 - Warranty Term**

The warranty period is of 12 months starting from the equipment's installation date since the installation be made during the 90 days after the purchasing date and under the terms of this certificate.

- 1** - The maximum storage period must be 3 months from the purchasing date. In the case of this period is exceeded, the warranty period will start counting even if the product is kept in storage.
- 2** - Upholstered parts are under warranty for 6 months.
- 3** - Operating light mirror, fuses, wires, and transformers are not covered by this warranty.
- 4** - The warranty is limited to the repair or replacement of defected parts and does not cover defects originated by:
  - a) Non-observance of the use and maintenance instructions;
  - b) Falls, crashes and inadequate storage;
  - c) Effects caused by the agents of nature;
  - d) Damage caused to the upholstered parts, due to improper use of chemical products, exposure to inadequate weather conditions, contact with fabrics, leather, disposable gloves, paints, pigmented detergents, cutting or piercing objects;
  - e) Damage to painted parts and plastic covers caused by improper use of chemical products or by contact with disposable gloves and cutting or piercing objects;
  - f) Connection with the wrong voltage power supply.
- 5** - This warranty will cease when:
  - a) The warranty period expires;
  - b) The user makes any changes to the product or uses inappropriate accessories;
  - c) The warranty certificate presents any scraping or changes;
  - d) Installation and technical service is done by an unauthorized technician;
  - e) The equipment is stored for more than 6 months after the invoice issuing date;
  - f) For the use of spare parts other than Olsen originals.
- 6** - The repairing or replacement of parts under the warranty period will not extend its original expiring date.
- 7** - The expenses resulted from the equipment's installation, scheduled maintenance, travel and accommodation costs of the service staff involved in the calls for service for the equipment's installation or repair will run under the owner's responsibility and in accordance with the distributor's norms.
- 8** - The customer, after checking the installation, must keep the purchasing invoice and the equipment's serial number until the warranty period ends, as in the case of service needed, the serial number, the invoice number, and purchasing or installation date will be required.
- 9** - On all calls for service during the warranty period, the equipment's serial number and a copy of the purchasing or installation invoice must be shown. In the case of this information is not presented, the service will not be done under warranty.

## 19 - Message from the President

Olsen and its clients:  
A successful relationship.

I have linked my name to the factory and to the dental and medical equipment that are currently produced and trade in more than 100 countries. Having in mind the responsibilities and long-term response to this initiative.

Our products are modern, innovative, durable and low maintenance costs. These characteristics have been achieved thanks to our competent and dedicated team, which makes me very proud for many reasons, as they are always giving the best of their creative capacity to our clients.

Our company will always be open to all those who prefer Olsen products, for any necessary information and technical assistance, but especially for comments regarding the relationship with customers. We expect this connection to always, bring you satisfaction, resulting in more and more benefits to all of us.

Cesar Olsen

[www.olsen.odo.br](http://www.olsen.odo.br)

***/// Olsen***

Equipment made to last

Registration at Ministry of Health 1028130009

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