









IMPORTANT

CLEANING AND DISINFECTION

FOLLOW OUR TIPS TO PRESERVE YOUR EQUIPMENT

WARNINGS, NOTES AND CAUTIONS

-  This summary highlights key information from the manual but does not replace reading the Cleaning and Disinfection section of the manual.
-  Protect the equipment from direct sunlight to avoid premature aging of its covers.
-  Turn off the equipment before cleaning and disinfection.
-  Cleaning and disinfection must be performed before each patient and at the end of the workday.
-  Wear protective gloves, a mask, and goggles during sanitation, following biosafety standards.
-  Whenever possible, use disposable barriers and replace them between patients to preserve the equipment's surfaces and finishes.

PRODUCTS

WHAT TO USE



- Microfiber cloth, cotton, soft fabric, paper towel, or wet wipe;
- Soft sponge (non-abrasive);
- Brush with nylon or soft bristles;
- Neutral detergent (e.g., PROLYSTICA® 2X CONCENTRATE by STERIS CORPORATION);
- Disinfect with 0.52% Hydrogen Peroxide (e.g., OXIVIR TB, EPA Reg. No. 70627-56);
- For leather upholstery, use 3M™ Leather Revitalizer;

WHAT NOT TO USE




- Abrasive sponges, hard bristle brushes, steel wool, or sharp objects that can scratch or damage surfaces;
- Chemical products like alcohol, acetone, thinner, corrosive acids, concentrated bleaches, or solvents, which may damage surfaces or compromise material integrity.

PURGING WATER AND AIR LINES (CLEANING THE HOSES)



- Place the pneumatic coupling in the water unit;
- Activate the coupling, purging the water and air supply for at least 30 seconds;
- Repeat steps 1 and 2 for each coupling used during treatment;
- For cleaning, disinfection, and sterilization of the **INSTRUMENTS, WATER SYSTEM, and SUCTION SYSTEM**, refer to the **CLEANING AND DISINFECTION** section of your equipment manual.

REMOVABLE COMPONENTS



- Moisten a sponge or cloth with water and neutral detergent, diluted according to the instructions in the user manual;
- Apply the solution to the entire surface, ensuring it remains in contact for at least 3 minutes;
- Using a soft sponge or brush, scrub the entire external and internal surface until all residue is completely removed;
- Use a cloth or sponge soaked in water to remove all detergent from the surface;
- Dry the component with a soft, clean cloth or paper towel;
- Disinfection can be performed with a cloth soaked in 0.52% hydrogen peroxide, such as OXIVIR TB EPA, following the manufacturer's instructions. Refer to the **CLEANING AND DISINFECTION** section of your equipment manual.

NON-REMOVABLE COMPONENTS



- Upholstery, plastic parts, painted parts, and metal surfaces;
- Dampen a cloth with water and neutral detergent, diluted according to the instructions in the equipment's user manual;
- Scrub the entire surface until all visible dirt is removed, then dry the surface with a dry cloth or paper towel;
- Dental light must be cleaned with a microfiber cloth, paper towel, or wet tissue after each procedure due to chemical attack and the appearance of stains resulting from phosphoric acid splashes on the reflective surface of the mirror;
- To ensure the cleanliness of PVC upholstery, it is recommended to clean only with water and neutral detergent. For leather upholstery, regular application of moisturizer is recommended, with 3M™ Leather Revitalizer being the preferred product.