Quality Quick Operation Guide - Description of Parts and Controls of Equipment

Before using the equipment, please read the *User Manual Dental Line* carefully.


The purpose of this guide is to identify commands, parts, and instruments that may be available to the equipment, including optional items. This guide does not define the equipment configuration that may vary according to the customer's request.

In case of loss of the *User Manual* and/or *Operation Guide*, please contact us by e-mail export3@olsen.odo.br informing the serial number of your equipment.

Delivery System

Before Turning On the Equipment:



- 1A - Water Tank
- 2A - Pressure Valves

- 1° - Check if the equipment is properly installed in accordance with the instructions in *chapter 6* of the *User Manual Dental Line*;
- 2° - Turn On the air compressor that supplies the equipment;
- 3° - Open the main water supply valve of the equipment;
- 4° - Turn On the equipment's power circuit breaker.

Turning On the Equipment:


- 1B - On/Off Switch with LED Indicator
- 2B - Fuse Holder

- 1° - Turn On the On/Off switch;
- 2° - Activate the "Seat Up" command on the chair base until the seat reaches its maximum height limit;
- 3° - Close the pressure valve;
- 4° - Remove the water tank and make sure it filled with water. If necessary, fill it with drinkable water up to the limit indicated in tank and reattach it again. If you wish, to place prophylactic products (low concentration);
- 5° - Open the pressure valve and check if the water tank is well attached and free of leaks;
- 6° - Before starting the use equipment, verify the operation of all instruments and commands available and perform the cleaning and sterilization of the instruments.

Chair Base

QUALITY

- 1 - Zero Position
- 2 - Work Position
- 3 - Backrest Up (↗)
- 4 - Seat Down (↓)
- 5 - Backrest Down (↙)
- 6 - Seat Up (↑)
- 7 - On/Off Operating Light
- 8 - Spittoon Bowl Flush
- 9 - Progressive/Activation Pedal

Commands

Operating Lights

Evolution Plus

LED Premium

Concept Plus

- 1 - Control Knob
- 2 - Control Rod
- 3 - Sensor

- 1° - To activate the Operating Light, press the Operating Light On/Off button at the base of the chair (section Commands - 7), this action also enables the sensor (3) and the control button/ rod (2);
- 2° - To change the lighting intensity of the Evolution Plus Operating Light, keep the command button - Commands section (7) - pressed for more than 2 seconds;
- 3° - To turn on or off the Concept Plus and LED Premium Operating Light, pass your hand in front of the sensor (3) at approximately 5 cm distance or hold the control rod (2), in the case of Concept Plus, or turn the control button (1), in the case of the Premium LED spotlight;
- 4° - To change the lighting intensity of the Concep Plus and Premium LED reflectors, keep your hand in front of the sensor (3) or hold the control rod (2), in the case of Concept Plus or turn the control knob (1), in the case of the Premium LED Operating Light.

Water Unit


- 1 - Cup Filler Control
- 2 - Spittoon Bowl Control
- 3 - Area for Instruments

The water unit may display water controls for the bowl flush and cup filler according to the equipment configuration. The cup filler water spout is removable.

To control the flow of water for cup filler (1) and spittoon bowl (2).

Assistant Module


- 1 - Instruments Holder
- 2 - Vacuum Pump Ejector
- 3 - Venturi Saliva Ejector
- 4 - Vortex Saliva Ejector

QUALITY


Quality Quick Operation Guide - Description of Parts and Controls of Equipment

Working Tables

Logic Working Table


Cross Flex Working Table

- 1 - Scaler with LED Light
- 2 - Coupling for Brushless Electrical Micromotor
- 3 - Coupling for Pneumatic Micromotor
- 4 - Instruments Control Panel
- 5 - Capacitive Handle
- 6 - Coupling for High-speed Handpiece
- 7 - Coupling for Optical Fiber Handpiece
- 8 - 3-Way Syringe
- 9 - Stainless Steel Tray
- 10 - Command Panel

The Logic working tables can present some configurations of command panels, varying according to available instruments.

Only move the working table through the capacitive handle

- 1 - Electrical Micromotor Speed Control
- 2 - Electrical Micromotor Direction Control
- 3 - Scaler Power Control and General Function Selector - Perio/Endo
- 4 - Scaler Water Control
- 5 - Brushless Micromotor Water Control
- 6 - Optical Fiber Water Control
- 7 - Thermo Comfort LED On/Off
- 8 - Fiber Optics On/Off


Flex Arm


The Flex Arm features 2 articulations (1) and pneumatic lock that are disabled by activating the pneumatic valve (2), allowing the working table movement.

Command Panels


Blue Touch Panel Models

Blue Touch Panel - X-Ray Viewer


Command Panel - Brushless Electrical Micromotor

- 1 - Seat Up/Position 1
- 2 - Backrest Up/Position 2
- 3 - Seat Down/Position 3
- 4 - Backrest Down
- 5 - Work Position
- 6 - Zero Position
- 7 - Spittoon Bowl Flush
- 8 - Cup Filler
- 9 - On/Off Operating Light
- 10 - On/Off Optical Fiber Light
- 11 - On/Off Air Jet System
- 12 - On/Off Thermo Comfort System
- 13 - X-Ray Viewer

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