Sprint Quick Operation Guide - Description of Parts and Controls of Equipment

Before using the equipment, please read the *User Manual Dental Line* carefully.


The purpose of this guide is to identify commands, parts, and instruments that may be available to the equipment, including optional items. This guide does not define the equipment configuration that may vary according to the customer's request.

In case of loss of the *User Manual* and/or *Operation Guide*, please contact us by e-mail export3@olsen.odo.br informing the serial number of your equipment.

Delivery System

Before Turning On the Equipment:

- 1A - Water Tank
- 2A - Pressure Valves



- 1° - Check if the equipment is properly installed in accordance with the instructions in *chapter 6* of the *User Manual Dental Line*;
- 2° - Turn On the air compressor that supplies the equipment;
- 3° - Open the main water supply valve of the equipment;
- 4° - Turn On the equipment's power circuit breaker.

Turning On the Equipment:


- 1B - On/Off Switch with LED Indicator
- 2B - Fuse Holder

- 1° - Turn On the On/Off switch;
- 2° - Activate the "Seat Up" command on the chair base until the seat reaches its maximum height limit;
- 3° - Close the pressure valve;
- 4° - Remove the water tank and make sure it filled with water. If necessary, fill it with drinkable water up to the limit indicated in tank and reattach it again. If you wish, to place prophylactic products (low concentration);
- 5° - Open the pressure valve and check if the water tank is well attached and free of leaks;
- 6° - Before starting the use equipment, verify the operation of all instruments and commands available and perform the cleaning and sterilization of the instruments.

Chair Base


Commands

- 1 - Zero Position
- 2 - Work Position
- 3 - Backrest Up (↗)
- 4 - Seat Down (↓)
- 5 - Backrest Down (↙)
- 6 - Seat Up (↑)
- 7 - On/Off Operating Light
- 8 - Spittoon Bowl Flush
- 9 - Progressive/Activation Pedal

Remote Foot Pedal

- 1 - Seat Up/Position 1
- 2 - Backrest Up/Position 2
- 3 - Seat Down/Position 3
- 4 - Backrest Down/Spit Position
- 5 - Work Position
- 6 - Zero Position
- 7 - Progressive/Activation Pedal
- 8 - Cup Filler
- 9 - On/Off Operating Light
- 10 - Spittoon Bowl Flush
- 11 - Metal Handle

To save a work position:

- 1° - Press the Zero Position button;
- 2° - Adjust the backrest to the desired position;
- 3° - Adjust the seat to the desired position;
- 4° - Press the Work Position button for 5 seconds. The equipment will emit 2 long beeps;
- 5° - Press command 1 (Seat Up) while the equipment beeps; The equipment will emit 2 short beeps to confirm the procedure;
- 6° - Press the Zero Position button.

Operating Lights

Evolution Plus

LED Premium

- 1 - Control Knob
- 2 - Control Rod
- 3 - Sensor

1° - To activate the Operating Light, press the Operating Light On/Off button at the base of the chair (section Commands - 7), this action also enables the sensor (3) and the control button/rod (2);

Concept Plus

2° - To change the lighting intensity of the Evolution Plus Operating Light, keep the command button - Commands section (7) - pressed for more than 2 seconds;


3° - To turn on or off the Concept Plus and LED Premium Operating Light, pass your hand in front of the sensor (3) at approximately 5 cm distance or hold the control rod (2), in the case of Concept Plus, or turn the control button (1), in the case of the Premium LED spotlight;

4° - To change the lighting intensity of the Concept Plus and Premium LED reflectors, keep your hand in front of the sensor (3) or hold the control rod (2), in the case of Concept Plus or turn the control knob (1), in the case of the Premium LED Operating Light.

Sprint Quick Operation Guide - Description of Parts and Controls of Equipment


Working Table

The Sprint working table can present some configurations of command panels, varying according to available instruments.

- 1 - Instruments Control Panel
- 2 - Instruments Holder
- 3 - Stainless Steel Tray
- 4 - Bicarbonate Reservoir
- 5 - X-ray Viewer


Instruments Panel

Command Panel


- 1 - Pneumatic Brake Valve*
- 2 - Scaler Power Control
- 3 - Scaler Water Control
- 4 - Optical Fiber Water Control
- 5 - Optical Fiber On/Off
- 6 - Electric Micromotor Speed Control
- 7 - Electric Micromotor Water Control
- 8 - Electric Micromotor Rotation Direction Inverter

Brushless Micromotor Command Panel


- 1 - Display
- 2 - Program
- 3 - Directional
- 4 - Gear Ratio
- 5 - Setting
- 6 - Manual Start
- 7 - Adjust Keys

Assistant Module

- 1 - Instruments Holder
- 2 - Vacuum Pump Ejector
- 3 - Venturi Saliva Ejector
- 4 - Vortex Saliva Ejector

Water Unit


- 1 - Area for Instruments
- 2 - Spittoon Bowl Manual Control
- 3 - Cup Holder Manual Control
- 4 - Spittoon Bowl Control with Timer
- 5 - Cup Holder Control with Timer

To activate the water control in the Spittoon Bowl, turn the Spittoon Bowl Manual Control (2). To activate the manual water control on the Cup Holder, keep the Cup Holder Manual Control (3) pressed.


When the Water Unit is equipped with timed controls, to release timed water in the spittoon bowl, press the Spittoon Bowl Control with Timer (4), to release timed water in the Cup Holder, press the Cup Holder Control with Timer (5). In addition to being available on the side of the Water Unit, the controls (4) and (5) may be available on the base control or on the remote pedal.

Working Table Arm

The **Straight Arm** allows 5 manual height adjustment levels.

To lower/raise the arm: holding the arm in the part closest to the equipment column, pull it firmly upward to disengage it. Adjust the lock (1) to the desired height and then push the arm down until it locks into place.


The **Articulated Arm** features 2 articulations between the delivery unit and the working table and mechanical lock system. It requires manual adjustment by a handle (2) to adjust the height from the floor.

The **Flex Arm** features 2 articulations (3) and pneumatic lock that are disabled by activating the pneumatic valve (4), allowing the working table movement.

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