USER MANUAL LAB CHAIR

Attention: Carefully read this manual
before operating the Olsen's equipment!

Olsen
Equipment made to last

Index

1 - Introduction	5
2 - Package Contents	5
3 - Equipment Presentation	5
3.1 - Standard Items	5
3.2 - Optional Items	5
4 - Parts Identification	6
4.1 - Applied Parts	6
4.2 - Accessories and Detachable Parts	6
5 - Description and Operating of the Equipment	6
5.1 - Turning On the Equipment	7
5.2 - Headrests	7
5.2.1 - Articulated Headrest	7
5.2.2 - Multiarticulated Headrest (Optional)	7
5.3 - Chair Base Commands	7
5.3.1 - Movements Interruption	8
5.3.2 - Positioning of the Patient, the Operator and Others	8
5.4 - Double Bracket Locks	8
5.5 - Upholstered Leg Rest	8
5.6 - Heels Support	9
5.7 - Folding Step	9
5.8 - Stainless Steel Tray	9
6 - General Features	10
7 - Installation Requirements	10
7.1 - Pre-installation	10
7.2 - Electrical Installation	10
8 - Installation	11
8.1 - Position	11
8.2 - Olsen Certified Technical Assistance Network	11
9 - Cleaning and Disinfection	11
9.1 - Upholsteries and Plastic Covers	11
9.2 - Painted Parts	11
9.3 - Sterilization by Autoclave	12
10 - Technical Features	13
10.1 - Electromagnetic Compatibility (EMC)	14
11 - Dimensional	16

12 - Symbology.....	17
13 - Important Notes.....	17
13.1 - General Cares - Compulsory Reading	18
13.2 - Disposal.....	18
13.3 - Transport and Storage.....	19
13.4 Contraindication.....	19
13.5 - Finalizando.....	19
14 - Troubleshooting.....	19
15 - Preventive Review	20
16 - Warranty Terms	21
17 - Message from the President.....	22

1 - Introduction

Congratulations on the great choice!

You have acquired an equipment developed for the outpatient and clinical environment, aimed primarily for laboratory testing and provide patient comfort during procedures.

This manual is supplemented to the Lab Chair Quick Operation Guide provided with your equipment and provides all the information you need to get the most out of your equipment, so read them carefully before using it.

2 - Package Contents

Check out the equipment package contents:

Standard Items:

- 1 Lab Chair
- 1 Upholstered Leg Rest Set
- 1 Heel Support Set
- 4 Bushings # 10 for Fixing the Chair
- 1 Lab Chair Quick Operation Guide

Optional Items:

- 1 Multiarticulated Headrest

3 - Equipment Presentation

3.1 - Standard Items

- Automatic Chair with 4 Commands
- Upholstered Leg Rest
- Joystick Control on Chair Base
- Articulated Headrest
- Stainless Steel Retractable Tray
- 2 Oil-Free Bosch Motors
- Folding Step
- Heel Support

3.2 - Optional Items

- Lether Upholstery
- Multiarticulated Headrest

4 - Parts Identification

- A - Articulated Headrest
- B - Backrest
- C - Upholstered Leg Rest
- D - Seat
- E - Heel Support
- F - Stainless Steel Tray
- G - Doble Bracket Locks
- H - Folding Step
- I - Joystick Control
- J - Electrical Panel

4.1 - Applied Parts

The following items of the Lab Chair are considered as applied parts to the patient:

Standard Items:

- Upholsteries;
- Articulated headrest;
- Heel support;
- Upholstered leg rest;

Optionals Items:

- Multiarticulated headrest.

4.2 - Accessories and Detachable Parts

The following items of the Lab Chair are considered as accessories or detachable parts or:

Detachable parts:

- Articulated headrest;
- Heel support;
- Upholstered leg rest.

Accessories and Optionals:

- Multiarticulated headrest.

5 - Description and Operating of the Equipment

Technical Name: Lab Chair

Trade Name: Lab Chair

The Lab Chair is an automated equipment for patient accommodation on clinical and laboratory examinations, being a multifunctional chair, ideal for gynecology, obstetrics, collection and blood tests.

It is equipped with two (2) electrical motors to perform four (4) independent movements. The motors are controlled by Joystick on the chair base.


Before starting to use your equipment make sure that it is properly installed, connected to the mains, with the circuit breaker On and the mains powered.

5.1 - Turning On the Equipment

Located on the chair base, under the footrest and easy access, the electrical panel provides access to the protection fuses (2) and the On/Off switch with LED (1).

To turn On the equipment: press the On/Off switch (1) to the "I" position. The LED on the On/Off switch (1) will turn green, indicating that the equipment is ready for use.

To turn Off the equipment: press the On/Off switch (1) to the "O" position.

5.2 - Headrests

5.2.1 - Articulated Headrest

The articulating headrest (1) is characterized by the easy position adjust, always keeping the same distance between the headrest (1) attaching mechanism and its bar (3).

To adjust the headrest position: turn the knob (2) counterclockwise to loosen it and then adjust to the desired position. Turn the knob (2) clockwise to lock in the set position.

To adjust the height of the headrest: move the bar (3) until you reach the position. The headrest bar allows adjustment with up to 17 cm from the backrest.


5.2.2 - Multiarticulated Headrest (Optional)

The multiarticulated headrest (4) is characterized by the multiple adjustable positions that guarantee more comfort for the patients during the clinical and surgical procedures. The multiarticulated headrest can also be adjusted for the care of wheelchair patients (7).

To adjust the position of the headrest (4): turn the knob (5) counterclockwise to loosen it and then adjust it to the desired position. Turn the knob (5) clockwise to lock in the set position.


To adjust the height of the headrest (4): move the bar (6) until you reach the position. The headrest bar allows adjustment with up to 17 cm from the backrest.

5.3 - Chair Base Commands

The joystick on the chair base allows the movement of the backrest and the seat. Its controls can be selected by foot, avoiding the cross contamination.

To perform chair movements: press the joystick at the points indicated on the side.

5.3.1 - Movements Interruption

All Lab Chair movements are of continuous activation, i.e., while the operator presses the control button, the equipment realizes the movement. The movement is interrupted immediately when the operator releases the control button.

5.3.2 - Positioning of the Patient, the Operator and Others

During the equipment movement, the operator and other persons close to him/her should position themselves on the equipment sides, respecting the minimum distance of 50 centimeters, remaining out of the equipment movement area and its components.

To ensure safe and proper patient positioning, the operator should advise him to keep his hands in sight during the entire equipment movement on the equipment upholstery or on the specific accessories for performing the procedure. If the operator sees any situation that may generate any kind of risk to the patient, the chair must be stopped immediately.

The presence of the operator, persons or objects in the equipment movement areas and its components may cause damage to the equipment and/or impair its correct operation.


All equipment's automatic movements must be supervised by the operator. If the operator sees any situation that may generate any kind of risk to the patient, the chair must be stopped immediately.

5.4 - Double Bracket Locks

The double bracket locks allow the coupling of 2 accessories, taking up less space on the side. The plastic handles (1) allow height and rotation on the shaft (distance the accessory to seat) adjustment,

Simply turn counterclockwise to loosen and clockwise to lock.




5.5 - Upholstered Leg Rest

Manufactured with the same upholstery of the Lab Chair, the upholstered leg rests (1) have a removable support (3). With the aid of the bracket locks, the upholstered leg rest (1) feature height and depth adjustment.

The upholstered leg rests (1) have 2 locking points (4).

The leg rest features upholsteries (1) fitted on the metal support (3), allowing a free spin on the fitting shaft (5).

To fit the legrest (1): loosen the knob (2) and fit the leg rest onto the metal support (3). Tighten the handle (2) to secure the leg rests (1).



5.6 - Heels Support

They offer support for the heels and are available with double bracket locks (3) that allow the simultaneous fixation of the heels support (1) and the leg rest (section 5.5).

The heel supports (1) can be adjusted for height, position of the heel supports (1) and rotation on their axis.


In addition to bracket locks adjustments, heel supports (2) can be adjusted for the heel pad distance according to the patient profile or procedure to be performed.

To fix the heel support (1): fit the support (1) into the holes (4) along the metal support (2) according to the desired distance.

5.7 - Folding Step

It assists in the patient's embarkation and allows the best approximation of the professional to the equipment since the step (1) can be repositioned folded.

5.8 - Stainless Steel Tray

Just below the seat upholstery is positioned the stainless-steel tray, to make it easier to perform procedures and collections.

The stainless-steel tray has an exclusive compartment and can be removed for cleaning and sterilization.

6 - General Features

- **Upholsteries:** are assembled on a very resistant frame, with soft foam and coated with flexible laminated PVC or leather, which besides being resistant, providing easy asepsis for the whole.
- **Mechanical structure:** manufactured in SAE 1020 rolled steel profiles and welded by MIG process, to guarantee the strength and durability.
- **Electrostatic painting:** it is applied to all structural metal parts of the equipment. The special polyurethane paint provides the equipment with high durability coating and has antibacterial properties according to JIS Z 2801: 2000, where in 24 hours, the bacterial reduction is greater than 99.9%.
- **Plastic covers:** made of high-strength ABS with acrylic cover, the plastic covers do not require painting and allowing polishing for recovery in case of wear or scratches.
- **Electrical system:** the equipment can work in frequencies of 50 or 60 Hz and can be configured to be connected to the voltages of 118/127/220/230 Volts by an authorized technician. The maximum supply voltage of electronic boards, motors, and other controls is 24 V. The electrical system has On/Off switch and protection fuses.
- **Spindle drive motor:** all motors used in the equipment are produced by *Robert Bosch of Brazil* and have advantages such as low noise, the absence of oil reservoir, uniformity in displacement, reduction of energy consumption and low cost for maintenance. The motors also have protection systems that act on the mechanical and electrical parts of the set.
- **Pressurized suspension:** the high pressure nitrogen cylinder relieves the load of the gear motor lifting system of the equipment in 100 kg. This means that a patient with this approximate weight is not acting on the electromechanical lifting components.

7 - Installation Requirements

7.1 - Pre-installation

The pre-installation must be guided by the Olsen authorized service to ensure that the environment is adequate the specified requirements for operation of the equipment (*chapter 8 - Installation*) and the positioning in which it is to be used. In this step, the electrical connections must be prepared to supply the equipment.

This equipment is not designed to be installed or operated in a surgical center.

7.2 - Electrical Installation

The power grid must be single-phase, has specific grounding and 10 A/3 mA DR circuit breaker. The circuit breaker must exclusively supply for the Lab Chair and must be easy and quick to disconnect from the mains. If the power grid presents voltage variation, the installation of a surge protection device is required. The power switch must be turned off at the end of the workday.

There is a table below for sizing the electrical installation:

Tension (V)	Wire Diameter (mm ²)	Distance (m)	Current (A)
118/127/220/230	2,5	Up to 20	10

This equipment should only be connected to a power source with protective grounding. There is a risk of electric shock!

8 - Installation

The Lab Chair installation must be done by Olsen certified technician consists in assembling legs rest, heel support, headrest and checking the mains voltage which will supply the equipment and, if necessary, make the tension adjustment. In the case of optional items have been purchased, these should also be installed at this time.

The technician will also verify if the equipment complies with the purchasing order, if it maintains its integrity, and will provide operation, cleaning and maintenance guidance.

When the installation is complete, observe the following items:

- Both right and left double bracket locks are installed with the upholstered leg rest and heels support;
- All joystick commands work perfectly;
- The On/Off switch, when pressed, turns on the device and lights up the green indicator LED;
- The stainless-steel tray slides perfectly;
- The upholsteries of backrest, seat, headrest and leg rest are intact.

In the case of optional items have been purchased, check the items below according to the acquired optional item:

- The multiarticulated headrest are installed correctly;
- Integrity of upholstery and finishes;

For instructions on the proper functioning of the equipment and its components, read Chapter 5 (Description and Operating of the Equipment) of this manual.

Installation should only be made by an authorized technician. Installation performed by an unauthorized person will result in loss of warranty.

8.1 - Position

To position the equipment, with the heels installed, put the backrest parallel to the floor and place the Lab Chair, and then take it to the place where it will be used. Perform the movements of raising and lowering the seat making sure that there is enough room for all movements.

8.2 - Olsen Certified Technical Assistance Network

To access the Olsen certified technical assistance network for installation and maintenance contact us by e-mail export3@olsen.odo.br or if you prefer +55 48 2106 6000.

9 - Cleaning and Disinfection

The whole sanitize process must be done with the use of gloves suitable for cleaning and protection, in addition to a mask and protective glasses, according to biosafety standards.

9.1 - Upholsteries and Plastic Covers

The plastic covers and upholsteries must be cleaned with a damp cloth containing neutral soap or detergent only. Olsen advises against using any chemical product to clean these parts, but in the case of disinfection products, it is important to check if it has suitable compatibility and specifications before use on these materials.

Never use hypochlorite or alcohol-based products.

9.2 - Painted Parts

It must be cleaned with a slightly damp cloth containing only soap or mild detergent.

Never use hypochlorite or alcohol-based products.

9.3 - Sterilization by Autoclave

The stainless steel tray can be sterilized by steam autoclaving at the following values (sterilization according to SS-374 resolution of 12/15/1995 and ISO 17665):

- a) 130° C, 2 bar, 15 minutes;
- b) 120° C, 1 bar, 30 minutes;
- c) 134° C, 2,2 bar, 4 minutes.

Note: these items withstand up to 1.000 autoclaving cycles.

All equipment items mentioned in this chapter must be sanitized and sterilized (when appropriate) prior to use.

Olsen is not responsible for defects, deformities, stains or abnormalities caused by improper use of chemical products, contact with tissues, leather, disposable gloves, inks, pigmented detergents and other organic or synthetic products.

10 - Technical Features

Power supply: 118/127/220/230 VC~.

A certified technician when installing the equipment must select the voltage;

Note: all equipment is adjusted in the factory to 220 VC~.

Number of phases: single phase.

Frequency: 50/60 Hz.

Power 118/127 V: 260 VA.

Power 220/230 V: 290 VA.

Protection fuses:

- For 220/230 VC~: 2 fuses of 1,5 AH (5 x 20 mm).
- For 118/127 VC~: 2 fuses of 3 AH (5 x 20 mm).

Power cable specification (in compliance with requirements 6.1 and 6.2 of IEC 60601-1-2:2010):

- Flexible Cable PP Circular 500 V 3 x 1 mm 247-5 NM 53-C5.
- Tri-polar Male Plug 10 A - 250 V (NBR 14136).

Electrical shock protection type: Class I equipment, according to IEC 60601-1-1 and IEC 60601-1-2 standards.

Protection degree: Type B;

Operation mode: non-continuous.

- Time On: 30 s.
- Time Off: 5 min.

Operating environment conditions:

- Temperature: between 15° and 30° C.
- Relative humidity: 30% to 70% non-condensing.
- Pressure: 75 kPa to 106 kPa.

Harmful water penetration's protection:

- Equipment: IPX0.

Thermal protection of the transformer: aperture with 130° C ± 3%.

Lifting capacity (maximum patient weight + accessories): up to 200 kg.

Equipment weight:

- Net: 90 kg (with accessories: 104 kg).
- Gross: 115 kg (with accessories: 129 kg).

The weight of accessories and optional items:

- Upholstered leg rest: 5,50 kg;
- Heels support: 8,6 kg;
- Multiarticulated headrest: 4,25 kg.

Seat height in relation to the floor:

- Minimum: 560 mm;
- Maximum: 810 mm.

10.1 - Electromagnetic Compatibility (EMC)

The Lab Chair needs special attention regarding electromagnetic compatibility and must be installed and put into use in accordance with the electromagnetic compatibility information presented in this chapter.

Radio frequency (RF) communication equipment, portable and mobile, can affect the Lab Chair.

Guidance and Manufacturer's Declaration - Electromagnetic Emissions		
The Lab Chair is intended for use in the electromagnetic environment specified below. The customer or the user should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic Environment - Guidance
RF emissions CISPR 11	Grupo 1	The Lab Chair uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference with nearby electronic equipment.
RF emissions CISPR 11	Class B	The Lab Chair is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings use for domestic purposes.
Harmonic emissions IEC 61000-3-2	Non-applicable	
Emissions due to voltage fluctuations/scintillation emissions IEC 61000-3-3	Non-applicable	

Guidance and Manufacturer's Declaration - Electromagnetic Emissions - I			
The Lab Chair is intended for use in the electromagnetic environment specified below. The customer or the user should assure that it is used in such an environment.			
IMMUNITY Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ±1 kV for input/output lines	± 2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on the power supply input lines IEC 61000-4-11	<5% U_T (dip > 95% in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (dip > 95% in U_T) for 5 s	<5% U_T (dip > 95% in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (dip > 95% in U_T) for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Lab Chair requires continued operation during power mains interruptions, it is recommended that the Lab Chair be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) Magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE U_T is the a.c. mains voltage prior to applications of the test level.			

Guidance and Manufacture's Declaration - Electromagnetic Immunity

The Lab Chair is intended for use in the electromagnetic environment specified below. The customer or the user should assure that it is used in such an environment.

IMMUNITY Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Conducted RF IEC 61000-4-6	3 V _{rms} 150 kHz to 80 MHz	3V _{rms}	<p>Portable and mobile RF communications equipment should not be used closer to any part of the Lab Chair, including cables, than the recommended separation distance calculated by the equation applicable to the transmitter frequency.</p> <p>Recommended separation distance</p> $d = [1,2]^{\frac{1}{2}}\sqrt{P}$ <p>80 MHz to 800 MHz</p> $d = [1,2]^{\frac{1}{2}}\sqrt{P}$ <p>800 MHz to 2,5 GHz</p> $d = [2,3]^{\frac{1}{2}}\sqrt{P}$ <p>Where <i>P</i> is the transmitter maximum output power rating in watts (W) according to the transmitter manufacturer and <i>d</i> is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a, should be less than the compliance level in each frequency range. ^b Interference may occur in the surroundings of equipment marked with the following symbol:</p> 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

a	Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measure field strength in the location in which the Lab Chair is used exceeds the applicable RF COMPLIANCE LEVEL above, the Lab Chair should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Lab Chair.
b	Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended Separation Distances Between Portable and Mobile RF Communications Equipment and the Lab Chair

The Lab Chair is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Lab Chair as recommended below, according to the maximum output power of the communications equipment.

Rated Maximum Output the Power of the Transmitter (W)	Separation Distance According to Frequency of Transmitter (m)		
	150 kHz to 80 MHz	80 kHz to 800 MHz	800 kHz to 2,5 GHz
	$d = [1,2]^{\frac{1}{2}}\sqrt{P}$	$d = [1,2]^{\frac{1}{2}}\sqrt{P}$	$d = [2,3]^{\frac{1}{2}}\sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,37	0,37	0,74
1	1,2	1,2	2,3
10	3,7	3,7	7,4
100	12	12	23



For transmitters rated at a maximum output power not listed above, the recommended separation distance *d* in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where *P* is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.



11 - Dimensional

Measured in Millimeters

12 - Symbology

The symbology complies with standards IEC 60601-1 and IEC 60878.

	Patient Chair Up		Patient Chair Down		Backrest Up
	Backrest Back		"ON" (Power)		"OFF" (Power)
	Non-Sterile		Caution		General Warning Sign
	Warning, Electricity		Operating Instructions		General Mandatory Action Sign
	Refer to the Instruction Manual		General Prohibition Sign		Stepping Prohibited
	Type B Applied Part		Manufacturer		Serial Number
	Sterilizable Up to the Temperature Specified		Protective Earth (Ground)		Earth (Ground)
	Keep Away from Heat		Humidity Limitation		Fragile, Handle with Care
	Temperature Limitation		Keep Dry		Heaping Up
	This Side Up		Alternating Current		
	Authorized Representative in the European Community				

13 - Important Notes

The reproduction and distribution of these instructions can only be made with prior permission from Olsen Indústria e Comércio S.A.

The technical features of the products described in this manual correspond to the time of its publication. Future technical improvements do not result in any right to update existing products.

The images presented in this manual are illustrative.

This equipment is designed to be free from interference from magnetic fields, external electrical influences, electrostatic discharges, pressure or pressure variation, provided that the equipment is transported, installed, operated and sanitized in accordance with the instructions for use contained in this manual.

13.1 - General Cares - Compulsory Reading

 Follow the instructions in *chapter 7* of this manual (*Installation Requirements*) to suit the electrical network where the equipment will be installed.

 Follow instructions for proper use of the equipment and its accessories as described in *chapter 5* (*Description and Operating of the Equipment*). Improper use could be harmful to equipment, which would not be covered by warranty.

 Follow the instructions in *chapter 9* (*Cleaning and Disinfection*) of this manual for daily cleaning of your equipment.

 Protect your equipment from direct exposure to sunlight. The equipment direct exposure to sunlight may cause premature aging of plastic covers and upholstery.

 Turn off the circuit breaker or unplug the equipment at the end of the working day.

 In case of joystick damage, discontinue use of the equipment, turn it off, and contact Olsen accredited assistance.

 The electric mains cable was developed for use only on the Lab Chair. The use of this component in other equipment may compromise the emissions and electromagnetic immunity of these.

 Only the authorized technician can replace the equipment's mains cable and internal fuses.

 This equipment is not suitable for use in the presence of a flammable anesthetic mixture with air, O₂ or Nitrous Oxide.

 This equipment should only be operated by physicians and nurses for clinical and outpatient procedures of the areas gynecology, obstetrics, collection and blood tests.

 Use only the mains cable supplied with the equipment for mains connection. Use of cables other than those specified (*chapter 10 - Technical Features*) may result in increased emissions or reduced immunity of the Lab Chair.

 Do not remove the plastic covers from the equipment. The risk of electric shock. Only the accredited technician is authorized to perform this procedure.

 In the case of damage to the electrical panel, adjacent plastic covers and motors plastic covers disconnect the equipment from the mains and contact the Olsen accredited service. The use of the equipment must be interrupted until maintenance is complete. The use of the equipment in these conditions offers a risk of electric shock!

 Do not perform maintenance or cleaning on the Lab Chair while it is in use or turned On.

 Do not install or use any electrical equipment over or near the Lab Chair. If necessary, the Lab Chair must be checked to see if it is functioning normally in the configuration in which it will be used.


 Do not fuse replacement if it is possible to touch the patient, even if unintentionally, during the procedure.

13.2 - Disposal

 Debris, residues and infectious materials resulting from the procedures performed on this equipment must be deposited in biological waste duly identified and in accordance with current legislation.

 For proper disposal of this equipment and its components and accessories, we recommend that it be sent to specialized recycling companies to ensure the best destination of each component without harm to the environment.

13.3 - Transport and Storage

-  It is recommended that the equipment's transportation and storage to be made in its original package.
-  Transport carefully protecting the equipment from falls and impacts.
-  Protect from moisture, rain exposure and direct contact with liquids.
-  Keep it sheltered from the sun.
-  Do not heap more than 5 volumes.
-  Do not move or store the equipment on uneven surfaces.
-  Temperature range for transportation and storage: -10° to +45° C.
-  Moisture limits for transport and storage: 20% to 70%.

13.4 Contraindication

-  This equipment is contraindicated for any use other than that for which it is intended, or to be operated by unauthorized personnel.

13.5 - Finalizando

At the end of the work day, observe the following instructions:

- Provide cleaning of equipment, cleaning of upholstery, upholstered leg rests, heels support and stainless steel tray, as indicated in *chapter 9* of this manual (*Cleaning and Disinfection*);
- Provide sterilization of the stainless steel tray, as indicated in *chapter 9* of this manual (*Cleaning and Disinfection*);
- Place the equipment in the embarkation position and turn the equipment On/Off switch Off;
- Disconnect the circuit breaker from the mains supplying the equipment.

14 - Troubleshooting

For solving possible problems in a simple and practical way, just follow the instructions in the following table:

Item	Problem	Causes	Solutions
1	The equipment doesn't perform any command.	1° - Equipment is not connected to the electric mains	1° - Connect the equipment to the main electricity
		2° - Electrical circuit breaker is switched Off	2° - Switch on the electrical circuit breaker
		3° - There is no power on the main electricity	3° - Call the Power Supply Company
		4° - Protection fuse is blown	5° - Contact the Olsen authorized service

If you are in doubt or find a problem with the equipment that is not mentioned in this chapter, stop using the equipment immediately and contact your authorized service center or contact us by e-mail export3@olsen.odo.br or by calling +55 48 2106 6000.

15 - Preventive Review

In order to extend your equipment's lifespan, Olsen has prepared a list of the main equipment items for which it recommends performing semi-annual preventive maintenance.

Performing a preventive maintenance by an accredited technician does not interfere with the equipment's warranty period.

STANDARD ITEMS
Check the articulated headrest
Check all joystick commands
Check the upholstered leg rests
Check the heels support
Check the stainless-steel tray
Spindle drive motors and joints lubrication
Check the left and right double bracket locks

OPTIONAL ITEMS
Check the multiarticulated headrest
Check of the armrests for blood donation

Allow only qualified Olsen technicians to perform installation and maintenance on your equipment and accessories.

Use only Olsen original parts and accessories. The use of non-original components can compromise the performance of the equipment, increasing its emissions or reducing its electromagnetic immunity.

Do not make adaptations, modifications or changes to the equipment or its components or accessories.

16 - Warranty Terms

The warranty period for this product is 12 months, counted from the equipment installation date, considering the 90 days legal warranty term, provided that the installation is performed within the period of 90 days from the product's Purchase Invoice issue and the fulfilled of other requirements of this certificate.

- 1** - The maximum term of storage is 3 months from the date of purchase of the product. If the storage period is exceeded, the guarantee is still in progress, even if the product is still stored.
- 2** - Upholstered parts are guaranteed for 6 months.
- 3** - Fuses, cables and transformers are not covered under warranty.
- 4** - The warranty is limited to repair or replacement of defective parts and does not cover defects originated by:
 - a) Non-compliance with the instructions for use and maintenance;
 - b) Falls, crashes and inadequate storage;
 - c) Action of nature agents;
 - d) Damage to upholstery, improper use of chemicals, exposure to unsuitable weather conditions, contact with tissues, leather, disposable gloves, paints, pigmented detergents, razors or any sharp instruments;
 - e) Damage to painted parts and plastic covers caused by improper use of chemicals or by contact with disposable gloves and sharp or piercing objects;
 - f) Connection to wrong voltage power supply.
- 5** - This warranty will cease when:
 - a) In the normal course of its period of warranty;
 - b) Make changes in the product not authorized by Olsen;
 - c) Adulterations in the purchase, installation or services document;
 - d) Installation or technical assistance made by a person not authorized by Olsen;
 - e) Failure to install the equipment for more than 90 days, counted from the purchasing date contained in the invoice;
 - f) By using non-genuine spare parts.
- 6** - The parts repairing or replacement during the warranty period will not extend the original expiry date.
- 7** - The expenses originated from the unit installation, scheduled preventive maintenance, travel, and hotel of the service staff involved in the calls for service for installation or units repair will run under the unit's owner responsibility and in accordance with the distributor's norms.
- 8** - The purchaser, after verifying the services performed in the installation and revision of the equipment, must date and sign the service order provided by the technician and keep along with his invoice of purchase of the equipment, failing which the warranty extension product when necessary.
- 9** - All service requests for warranty equipment must be made with the serial number of the equipment to be serviced and a copy of the purchase or installation document. If this information is not communicated, the service request will be made as not covered by the warranty.


17 - Message from the President

Olsen and its clients:
A successful relationship.

I have linked my name to the factory and to the Dental and Medical equipment which are currently produced and trade in more than 100 countries having in mind the responsibilities and long-term response to this initiative.

Our products are modern, innovative, durable and of low maintenance cost. These characteristics have been achieved thanks to our competent and dedicated team, which make me very proud for many reasons, as they are always giving the best of their creative capacity to our clients.

Our company will always be open to all those who prefer Olsen products, for any necessary information and technical assistance, but especially for comments regarding the relationship with customers. We expect this connection to always bring you satisfaction, resulting in more and more benefits to all of us.

Cesar Olsen

www.olsen.odo.br

+55 48 2106 6000
export3@olsen.odo.br

/// Olsen

Equipment made to last

Registration at Ministry of Health of Brazilian 10281300013

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