USER MANUAL

AUTOMATIC OPHTHALMOLOGY/ ENT.

Attention: Carefully read this manual
before operating the Olsen's equipment!

Olsen
Equipment made to last

Index

1 - Introduction	5
2 - Package Contents	5
3 - Equipment Presentation	5
3.1 - Standard Items	5
3.2 - Optional Items	5
4 - Parts Identification	6
4.1 - Applied Parts	6
4.2 - Accessories and Detachable Parts	6
5 - Equipment Description and Operating	7
5.1 - Turning On the Equipment	7
5.2 - Multiarticulate Headrest	7
5.3 - Chair Base Commands	7
5.3.1 - Movement Interruption	7
5.3.2 - Positioning of the Operator and the Patient	8
5.4 - Folding Step	8
5.5 - Operating Light (Optional)	8
5.5.1 - Concept LED (Optional)	8
5.5.2 - Surgical Operating Light (Optional)	8
5.6 - Water Unit (Optional)	9
5.7 - Venturi Saliva Ejector (Optional)	9
5.8 - Vacuum Pump Adapter (Optional)	9
6 - General Features	10
7 - Installation Requirements	10
7.1 - Pre-installation	10
7.2 - Electrical Installation	10
7.3 - Water for the Water Unit	10
7.4 - Sewage System	11
8 - Installation	11
8.1 - Equipment Positioning	11
8.2 - Olsen Certified Technical Assistance Network	11
9 - Cleaning and Disinfection	11
9.1 - Upholsteries and Plastic Covers	11
9.2 - Painted Parts	12
9.3 - Saliva Ejectors	12
9.4 - Water Unit and Strainer	12

9.5 - Sterilization by Autoclave.....	12
10 - Technical Features	13
10.1 - Electromagnetic Compatibility.....	14
11 - Dimensional.....	16
12 - Symbology.....	17
13 - Important Notes.....	18
13.1 - General Cares	18
13.2 - Disposal.....	19
13.3 - Transportation and Storage	19
13.4 - Contraindications.....	19
13.5 - Finishing.....	19
14 - Troubleshooting.....	20
15 - Preventive Review	20
16 - Warranty Terms	21
17 - Message from the President.....	22

1 - Introduction

Congratulations for the great choice!

You have acquired an equipment developed for the outpatient environment, settings primarily at the specialty for Ophthalmology, Otolaryngology and Radiologic specialties, in order to provide maximum comfort to the patient and facilitate the performance of procedures by the professional.

This manual is supplemented to the Automatic Ophthalmology/Ent Quick Operation Guide provided with your equipment and provides all the information you need to get the most out of your equipment, so read them carefully before using it.

2 - Package Contents

Check out the equipment package contents:

Standard Items:

- 1 Automatic Ophthalmology/Ent Chair
- 2 Armrests
- 4 Bushings (#10) for attaching the chair
- 1 Automatic Ophthalmology/Ent. Quick Operation Guide

Optional Items:

- Water Unit
- Vacuum Pump Adaptor Set
- Venturi Saliva Ejector

3 - Equipment Presentation

3.1 - Standard Items

- Automatic Chair with 4 Commands
- Swivel Armrests
- Joystick Control on Chair Base
- Multiarticulate Headrest
- 2 Oil-free Bosch Motors
- Folding Step

3.2 - Optional Items

- Water Unit
- Cup Filler in the Water Unit
- Vacuum Pump Adaptor Set
- Venturi Saliva Ejector
- Concept LED Operating Light
- Surgical LED Operating Light

4 - Parts Identification

- A - Multiarticulate Headrest
- B - Backrest
- C - Armrests
- D - Seat
- E - Folding Step
- F - Electric Panel
- G - Joystick Control

4.1 - Applied Parts

The following items of the Automatic Ophthalmology/Ent Chair are considered as applied parts to the patient:

Standard Items:

- Upholstery;
- Multiarticulate Headrest;
- Swivel Armrests.

Optionals Items:

- Cannula Venturi Saliva Ejector;
- Cannula Vacuum Pump;

4.2 - Accessories and Detachable Parts

The following items of the Automatic Ophthalmology/Ent Chair are considered as accessories or detachable parts:

Detachable parts:

- Multiarticulate headrest;
- Cannula vacuum pump;
- Cannula Venturi saliva ejector;

Accessories and Optional:

- Water unit;
- Vacuum pump set;
- Venturi saliva ejector;
- Surgical LED operating light;
- Concept LED operating light.

5 - Equipment Description and Operating

The Automatic Ophthalmology/Ent Chair is automated equipment for patient accommodation on clinical and outpatient exams especially for the Ophthalmology, Otolaryngology, and Radiology specialties.

It is equipped with 2 electric motors to perform 4 movements individually and independently, controlled by a joystick at the base of the chair.


Before starting to use your equipment make sure that it is properly installed and that the power supply meets the specifications of the equipment (*chapter 7 - Installation Requirements*).

5.1 - Turning On the Equipment

Located on the base of chair, with easy access, the electrical panel provides access to the protection fuse (2) and the LED On/Off switch (1).

To turn On the equipment: press the LED On/Off switch (1) to the "I" position. The switch will light green indicating that the equipment is turned On and ready for use.

To turn Off the equipment: press the LED On/Off switch (1) to the "O" position.

5.2 - Multiarticulate Headrest

The multiarticulate headrest (1) allows multiple adjustable positions that guarantee greater comfort for patients during the clinical and surgical procedures. Adjust the multiarticulate headrest (1) to treat wheelchair patients (4).

To adjust the headrest (1): turn the knob (2) counter-clockwise, adjust the headrest (1) to the desired position and then turn the knob (2) clockwise to lock it in position.


To adjust the distance from the headrest (1) to the backrest: Pull the headrest (1) away from the backrest. Do not exceed the 10 cm limit between the headboard and the backrest (3).

5.3 - Chair Base Commands

The joystick at the base of the equipment enables backrest and seat movement. Your commands can be selected with your feet, avoiding cross contamination.

To perform chair movements: press the joystick at the points shown beside.

5.3.1 - Movement Interruption

All Automatic Ophthalmology/Ent Chair movements are continuously activated, i.e.; they only act when the operator keeps the command pressed. Movement is stopped immediately when the operator releases the control button.

5.3.2 - Positioning of the Operator and the Patient

During the equipment movement, the operator and other people close to him/her should position themselves on the sides of the equipment, respecting the minimum distance of 50 cm, remaining outside the movement area of both the equipment and its components.

The patient's safe and proper positioning should be guided by the operator so that the patient remains with his or her hands and feet within the operator's field of vision, on the equipment padding, or on the specific accessories used for performing the procedure. The operator must guide the patient to remain in this position during the whole equipment movement. If the operator identifies any situation that may create a risk, the equipment movement must be stopped immediately.


The presence of the operator, persons or objects in the equipment's and its components movement areas may cause damage to the equipment and/or impair its correct operation.

All movement of equipment must be supervised by the operator. If the operator finds any situation that could pose a risk to the patient, he should immediately stop the movement of the equipment.

5.4 - Folding Step

Assists in patient boarding and allows the best approach of the professional to the equipment as the step (1) can be repositioned folded.

5.5 - Operating Light (Optional)

The Automatic Ophthalmology/Ent Chair features lighting options through operating lights with articulated arms that can be attached to frame, to the floor or ceiling.

5.5.1 - Concept LED (Optional)


The Concept LED operating light has polycarbonate protection shield, LED illumination and cyclic and gradual intensity control.

The head light has articulated arm and side handles for position adjustment.

The lighting is generated by white color LED, with luminous intensity of 8.000 up to 30.000 lux and a focus of approximately 10 x 5 cm.

To turn the operating light On/Off: press the operating light On/Off switch (1) on the chair base.

To change the light intensity: keep the operating light On/Off switch (1) pressed for more than 2 seconds.

5.5.2 - Surgical Operating Light (Optional)


It features 6 LED spotlights positioned and deployed in order to provide a wide area of illumination for outpatient procedures.

To turn the surgical operating light On/Off: use the switch (4).

To move the surgical operating light, use the metal handle (2) or the side handles (3). The handle (2) can be removed for autoclaving by turning it counterclockwise.

It is recommended that the side handles (3) be packed with a disposable material (ex.: plastic film) to avoid cross-contamination and damage to the equipment's plastic covers.

This operating light has a luminous intensity of 30.000 Lux.

The operating light Concept LED and the Surgical operating light do not have a lamp.

5.6 - Water Unit (Optional)


The water unit provides an easy connection to the water and sewage network (9) and a spittoon bowl (2), which helps on the procedures. The water unit is 90° swivel allowing your approach or out the patient. It has connection box (7) for water and sewage hoses (9), and electrical connection, providing socket (8) with a 10 A fuse protected outlet.

Optional instrument holders (6) can be provided with water units such as Venturi saliva ejector and vacuum pump adapter.

To obtain water for the bowl flush: use the adjustment valve (1).

To obtain water for the cup filler (3): press the button (2).

The water spouts (3 and 4) are detachable for easy cleaning, and can be adjusted to direct the water flow on the spittoon bowl (2).

5.7 - Venturi Saliva Ejector (Optional)

The Venturi saliva ejector (2) uses the compressed air to generate ejection. It has a solids filter (5) that prevents these solids from being thrown into the sewer system. It allows the attachment of 6.5 mm (1) and 9.5 mm (4) cannulas and may have, below the water unit a valve to adjust the ejection flow.


To activate the suction: remove the ejector (2) from the holder (3).

To deactivate the suction: place the ejector (2) in its holder (3).

For the proper functioning of ejector, it is necessary that the solids filter (5) is clean and the drain to the sewage system is properly installed and properly tilted.

Note: for the filter (5) cleaning information refer to the *chapter 9 (Cleaning and Disinfection)* of this manual.

The cannulas and the odontology compressor do not accompany the ejector.

5.8 - Vacuum Pump Adapter (Optional)

The kit can be installed in the water unit by providing the suction system of the vacuum pump as one more functionality integrated with the equipment. It features an body (3) with suction flow control and an adapter (1) for use with 6.3 mm cannula.



The vacuum pump connection hose is built-in the corrugated hose of the water unit and the holder (2) is provided with a device for the vacuum pump automatic switch-off.

The vacuum pump and the cannulas are not part of this device and can be purchased separately.

To activate the vacuum pump: remove the adapter body (3) from the holder (2). The vacuum pump (3) will start automatically. The suction control is made by rotating the adapter body (3), where the blue dots (4) alignment indicates the total opening of the ejection and, the totally opposite points indicate the closing of ejection.

To deactivate the vacuum pump: place the adapter (3) from the holder (2).

The adapter tip is removable (1), allowing the use of cannulas of 11 mm.

6 - General Features

- **Upholstery:** the upholstery is mounted on a very resistant frame, with soft foam and coated with flexible laminated PVC seamless, providing easy for the aseptic assembly.
- **Metal frame:** Olsen equipment is manufactured in SAE 1020 rolled steel profiles and welded by MIG process, giving strength and low weight to the set, facilitating transportation, installation, and maintenance of the equipment.
- **Electrostatic painting:** it is applied to all equipment's structural metal parts. The paint is produced from polyurethane and gives the equipment a high strength coating. The paint has an antimicrobial property, with performance in accordance with *JIS Z 2801: 2000*, where, in a 24-hour period, the bacteria reduction is greater than 99.9%.
- **Plastic covers:** made of high resistance ABS with acrylic coating does not require painting, so small scratches and wear can be repaired by polishing.
- **Electrical system:** the equipment can operate at frequencies of 50 or 60 Hz. The electrical system can be adjusted to connect to 118/127/220/230 V. If necessary, this setting is performed by an authorized technician during the equipment assembly. The maximum internal voltage for the power supply of electronic boards, motors, and other controls is 24 V. The electrical system has On/Off switch and protection fuses.
- **Spindle drive motors:** one of the most important features of the Olsen units lies on the seat and backrest movement system, which counts on Bosch electric motors with the spindle drive system. Besides being very compact, they are oil-free, work on a very low level of noise with smooth movements, low energy consumption, and low-cost maintenance.
- **Pressurized suspension:** the high pressure nitrogen cylinder relieves the load of the gear motor lifting system of the equipment in 100 kg. This means that a patient with this approximate weight is not acting on the electromechanical lifting components.

7 - Installation Requirements

7.1 - Pre-installation

The pre-installation should be guided by Olsen authorized assistance to ensure that the environment is adequate to receive the equipment, in accordance with the temperature, relative humidity and pressure required for the equipment correct operation (*chapter 10 - Technical Features*), as well as positioning in which it should be used. In this step, the electrical connections must be prepared to supply the equipment.

This equipment is not designed to be installed or operated in a surgical center.

7.2 - Electrical Installation

The power grid must be single-phase and have specific grounding and DR circuit breaker. The circuit breaker must exclusively supply the Automatic Ophthalmology/Ent Chair and must be easy and quick to disconnect from the mains. If the power grid presents voltage variation, the installation of a surge protection device is required. The power switch must be turned off at the end of the workday.

There is a table below for sizing the electrical installation:

Tension (V)	Wire Diameter (mm ²)	Distance (m)	Current (A)
118/127/220/230	2,5	20	10

This equipment should only be connected to a power source with protective grounding. There is a risk of electric shock!

7.3 - Water for the Water Unit

If your equipment has a water unit, it is essential that the water network has an easy access main valve to stop water flow, because the water main valve should be closed daily at the end of the workday.

The water pressure should be between 2.8 to 4.0 bar, with ideal pH (Hydrogenionic Potential) between 6.5 to 8.

It is recommended to use a filter before the external water supply of the equipment to prevent clogging in the internal water system.

7.4 - Sewage System

If your equipment has a water unit, the sewage system must have good hydraulic declivity (minimum of -2°) and should be preferably installed underground. The nominal diameter of the tubing must be of Ø40mm.

8 - Installation

The Automatic Ophthalmology/Ent Chair installation must be done by Olsen certified technician consists of the equipment assembly, as well as the compatibility check between the mains voltage and the equipment supply voltage and if necessary make the appropriate adjustment, before connecting the equipment to the network.

The technician will also verify if the equipment complies with the order if it maintains its integrity, and will provide operation, cleaning, and maintenance guidance.

When the installation is complete, observe the following items:

- The multiarticulate headrest is properly installed;
- The 4 joystick commands are working perfectly;
- The On/Off switch, when activated, turns on the equipment and turns on the green indicator LED;
- The upholsteries of backrest, seat, headrest and armrest are intact.

In the case of optional items have been purchased; check the items below according to the acquired optional item:

- If it has the water unit, it is correctly installed, has good water flow, adequate drainage of the spittoon bowl and no leaks.
- The accessories provided with the water unit are well fitted on the supports, have no leaks and have their perfect drive;
- The operating light turns On and Off normally and presents a smooth arm movement.

For instructions on proper operation of the equipment and its components, read Chapter 5 (Equipment Description and Operation) of this manual.

The installation of the equipment should only be performed by factory authorized personnel. Installation by unauthorized persons will void the warranty!

8.1 - Equipment Positioning

To position the equipment, recline the backrest and place the Automatic Ophthalmology/Ent Chair where it will be used. Perform the up and down movements of the seat and the backrest, making sure there is enough room to perform the movements.

8.2 - Olsen Certified Technical Assistance Network

To access Olsen certified technical assistance network for installation and maintenance, with the equipment serial number, contact us via e-mail at export3@olsen.odo.br or if you prefer by phone, call +55 48 2106 6000.

9 - Cleaning and Disinfection

It is indispensable to use gloves and mask, according to biosafety standards for the cleaning procedures of the equipment.

9.1 - Upholsteries and Plastic Covers

The plastic covers and upholsteries must be cleaned with a damp cloth containing neutral soap or detergent only. Olsen advises against using any chemical product to clean these parts, but in the case of disinfection products, it is important to check if it has suitable compatibility and specifications before use on these materials.

Never use hypochlorite or alcohol-based products.

9.2 - Painted Parts

It must be cleaned with a slightly damp cloth containing only soap or mild detergent and disinfection with alcohol 70%.

Never use hypochlorite or alcohol-based products.

9.3 - Saliva Ejectors

Daily disinfect the saliva ejector's hoses with an appropriate disinfection product for PVC hoses. It is highly recommended to read carefully the manufacturer's instructions to avoid damage to the devices.


Using a proper concentration of the cleaning product, make slow ejection with the saliva ejector until reaching the adequate quantity to clean the hoses. After exposing the hoses to the necessary time, make ejection of 1 liter of water to decrease the reacting chemical effects on the material.

The solids collector filter must also be cleaned weekly. The saliva ejectors efficiency may be impaired if this filter is clogged. In the case of saliva ejector's performance reduction, clean its filter. To clean the filters, follow the instructions below:

1° - Detach the covers (1);

2° - Remove the strainer (2) for cleaning;

3° - After cleaning, assemble the set again (1 and 2).

To handle the filters, it is essential the use of gloves and mask according to biosecurity standards.

9.4 - Water Unit and Strainer

The water spout (2) is detachable to facilitate the spittoon bowl (1) removal for cleaning.

To clean the spittoon bowl (1) remove the water spout (2), then the spittoon bowl's finishing (3) and the strainer (4) and then the spittoon bowl (1).

The cleaning process of the strainer (4) and the spittoon bowl (1) can be made with running water and neutral soap or detergent. Use non-abrasive sponge.

All contaminated debris and materials should be disposed of in biological waste.

All items of the equipment referenced in this chapter must be sanitized and sterilized (if applicable) prior to use.

9.5 - Sterilization by Autoclave

The Surgical operating light metal handle can be sterilized by steam autoclaving. To do this, wipe the part with a slightly damp cloth and use a specific envelope for steam sterilization.

Steam autoclave up to 135° C for 20 minutes at a pressure of 2 bars or 121° C for 40 minutes at a pressure of 1 bar. This item resists up to 1,000 autoclaving cycles.

Olsen is not responsible for defects, deformities, stains or abnormalities caused by improper use of chemical products, contact with tissues, leather, disposable gloves, inks, pigmented detergents and other organic or synthetic products.

10 - Technical Features

Power supply: 118/127/220/230 V.

A certified technician when installing the equipment must select the voltage.

Note: all equipments are adjusted in the factory to 220 V~.

Number of phases: single phase.

Frequency: 50/60 Hz.

Power 118/127 V: 260 VA

Power 220/230 V: 290 VA

Protection fuses:

- For 220/230 V~: 1,5 A H (5 x 20 mm);
- For 118/127 V~: 3 A H (5 x 20 mm).

Power cable specification (in compliance with requirements 6.1 and 6.2 of IEC 60601-1-2:2010):

- Flexible Cable PP Circular 500 V 3 x 1 mm 247-5 NM 53-C5;
- Tri-polar Male Plug 10 A - 250 V.

Electrical shock protection type (IEC 60601-1-1 and IEC 60601-1-2 standards): I Class.

Protection degree: B Type.

Operation mode: non-continuous

- Time On: 30 s.
- Time Off: 5 min.

Operating environment conditions:

- Temperature: between 15° C and 28° C;
- Pressure: 75 kPa to 106 kPa;
- Relative humidity: 30% to 70% non-condensing.

Harmful water penetration's protection: IPX0.

Thermal protection of the transformer: aperture with 130° C ± 3%.

Lifting capacity (maximum patient weight): up to 200 kg.

Equipment weight:

- Net: 85,0 kg;
- Maximum with accessories: 110,0 kg.

Weight of the optional items:

- Operating light LED: 7,55 kg;
- Surgical operating light: 8,20 kg;
- Water unit: 2,33 kg;
- Venturi saliva ejector: 1,10 kg;
- Vacuum pump adaptor: 0,78 kg.

Seat height from the floor:

- Minimum: 530 mm;
- Maximum: 780 mm.

10.1 - Electromagnetic Compatibility

The Automatic Ophthalmology/Ent Chair needs special attention regarding electromagnetic compatibility and must be installed and put into use in accordance with the electromagnetic compatibility information presented in this chapter.

Radio frequency (RF) communication equipment, portable and mobile, can affect the Automatic Ophthalmology/Ent Chair.

Guidance and manufacturer's declaration - electromagnetic emissions

The Automatic Ophthalmology/Ent Chair is intended for use in the electromagnetic environment specified below. The customer or the user of the Automatic Ophthalmology/Ent Chair should assure that it is used in such an environment.

Emission test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group1	The Automatic Ophthalmology/Ent Chair uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference with nearby electronic equipment. The Automatic Ophthalmology/Ent Chair is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings use for domestic purposes.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Not applicable	
Emissions due to voltage fluctuations/scintillation emissions IEC 61000-3-3	Not applicable	

Guidance and manufacturer's declaration - electromagnetic Immunity- I

The Automatic Ophthalmology/Ent Chair is intended for use in the electromagnetic environment specified below. The customer or the user of the Automatic Ophthalmology/Ent Chair should assure that it is used in such an environment.

IMMUNITY test	IEC 60601 Test Level	Compliance level	Electromagnetic environment guidance
Electrostatic Discharge (DES) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	±1 kV differential mode ±2 kV common mode	Mains power quality should be that a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	< 5% UT (dip > 95% in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles < 5% UT (dip > 95% in UT) for 5 s	< 5% UT (dip > 95% in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles < 5% UT (dip > 95% in UT) for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Automatic Ophthalmology/Ent Chair requires continued operation during power mains interruptions, it is recommended that the Automatic Ophthalmology/Ent Chair be powered from an uninterruptible power supply or a battery
Power frequency (50/60 Hz) Magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE UT is the a.c. mains voltage prior to applications of the test level.

Guidance and manufacture's declaration - electromagnetic immunity

The Automatic Ophthalmology/Ent Chair is intended for use in the electromagnetic environment specified below. The customer or the user of the Automatic Ophthalmology/Ent Chair should assure that it is used in such an environment.

IMMUNITY test	IEC 60601 Test Level	Compliance level	Electromagnetic environment guidance
Conducted RF IEC 61000-4-6	3 V _{rms} 150 kHz to 80 MHz	3V _{rms}	Portable and mobile RF communications equipment should not be used closer to any part of the Automatic Ophthalmology/Ent Chair, including cables, than the recommended separation distance calculated by the equation applicable to the transmitter frequency. Recommended separation distance $d = [1,2]^2 \sqrt{P}$ $d = [1,2]^2 \sqrt{P}$ $d = [2,3]^2 \sqrt{P}$ where <i>P</i> is the transmitter maximum output power rating in watts (W) according to the transmitter manufacturer and <i>d</i> is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range. ^b Interference may occur in the surroundings of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	

NOTE 1 At 80MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

a	Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measure field strength in the location in which the Automatic Ophthalmology/Ent Chair is used exceeds the applicable RF compliance level above, the Automatic Ophthalmology/Ent Chair should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Automatic Ophthalmology/Ent Chair.
b	Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended Separation Distances Between Portable and Mobile RF Communications Equipment and the Automatic Ophthalmology/Ent Chair

The Automatic Ophthalmology/Ent Chair is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Automatic Ophthalmology/Ent Chair can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Automatic Ophthalmology/Ent Chair as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of the transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d = [1,2]^2 \sqrt{P}$	150 kHz to 80 MHz $d = [1,2]^2 \sqrt{P}$	150 kHz to 80 MHz $d = [2,3]^2 \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,37	0,37	0,74
1	1,2	1,2	2,3
10	3,7	3,7	7,4
100	12	12	23
0,01	0,12	0,12	0,23
0,1	0,37	0,37	0,74



For transmitters rated at a maximum output power not listed above, the recommended separation distance *d* in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where *P* is the maximum output power rating of the transmitter in Watts (W) according to the transmitter manufacturer.

NOTE 1 - At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 - These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

11 - Dimensional

Measured in Millimeters

12 - Symbology

According to IEC 60601-1 and IEC 60878 standards.

	Patient Chair Up (Dental)		Patient Chair Down (Dental)		Backrest Back
	Backrest Up		Bowl Flush		Cup Filler
	Saliva Ejector		Suction Handpiece		Suction Handpiece with Hand Control Valve
	Non-Sterile		Hand Control Valve		Manufacturer
	"OFF" (Power)		"ON" (Power)		Caution
	General Warning Sign		Warning, Electricity		Operating Instructions
	General Prohibition Sign		Refer to the Instruction Manual		General Mandatory Action Sign
	Level		Sterilizable Up to the Temperature Specified		Serial Number
	Stepping Prohibited		Type B Applied Part		Keep Away from Heat
	Protective Earth (Ground)		Earth (Ground)		Temperature Limitation
	Humidity Limitation		Fragile, Handle with Care		This Side Up
	Keep Dry		Heaping Up		Alternating Current

EC REP	Authorized Representative in the European Community
---------------	--

13 - Important Notes

The reproduction and distribution of these instructions can only be made with prior permission from Olsen Indústria e Comércio S.A.

The technical features of the products described in this manual correspond to the time of its publication. Future technical improvements do not result in any right to update existing products.

The images presented in this manual are illustrative.

This equipment is designed to be free from interference from magnetic fields, external electrical influences, electrostatic discharges, pressure or pressure variation, provided that the equipment is transported, installed, operated and sanitized in accordance with the instructions for use contained in this manual.

13.1 - General Cares

 Follow the instructions in *chapter 7* of this manual (*Installation Requirements*) to suit the electrical network where the equipment will be installed.

 Follow the instructions for proper use of the equipment and its accessories as described in *chapter 5* (*Equipment Description and Operation*). Improper use may result in damage to equipment and would not covered under warranty.

 Follow the instructions in *chapter 9* (*Cleaning and Disinfection*) of this manual for your equipment daily cleaning.

 Before starting office activities, check compressor condition, if any.

 Protect your equipment from direct exposure to sunlight. Direct exposure of the equipment to sunlight may cause premature aging of plastic covers and upholstery.

 Turn off the circuit breaker or disconnect equipment from the mains and close the office water supply (when equipped with a water unit) at the end of the day work.

 In the event of damage to the joystick on the chair base, stop the equipment use, turn it off, and contact Olsen accredited assistance.

 The mains cable was designed for exclusive use on the Automatic Ophthalmology/Ent Chair. Use of this component in other equipment may compromise their emissions and electromagnetic immunity.

 Only the authorized technician can replace the equipment's mains cable and internal fuses.

 This equipment is not suitable for use in the presence of a flammable anesthetic mixture with air, O₂ or Nitrous Oxide. As well it is not suitable for use in Oxygen rich environments.

 This equipment should be operated only by physicians, and nurses who work in the areas of ophthalmology, otolaryngology and radiology.

 Use only the mains cable supplied with the equipment for mains connection. Use of cables other than those specified (*chapter 10 - Technical Features*) may result in increased emissions or reduced immunity of the Automatic Ophthalmology/Ent Chair.

 Do not remove the plastic covers from the equipment. There's risk of electric shock! Only the accredited technician is authorized to perform this procedure.

 In the case of damage to the electrical panel, adjacent plastic covers, and motor's plastic covers disconnect the equipment from the mains and contact the Olsen accredited service. The equipment use must be interrupted until the maintenance is complete. The equipment use in these conditions offers a risk of electric shock.!

 Do not perform cleaning or maintenance of the equipment while it is in use, with the patient or switched on.

 Do not install or use any electrical equipment over or near the Automatic Ophthalmology/Ent Chair. If necessary, the Automatic Ophthalmology/Ent Chair must be checked to see if it is functioning normally in the configuration in which it will be used.

 Do not perform the fuse replacement if it is possible to touch the patient, even if unintentionally, during the procedure.

13.2 - Disposal

 Debris, residues and infectious materials resulting from the procedures performed on this equipment must be deposited in biological waste duly identified and in accordance with current legislation.

 For proper disposal of this equipment and its components and accessories, we recommend that it be sent to specialized recycling companies to ensure the best destination of each component without harm to the environment.

13.3 - Transportation and Storage

 It is recommended that the equipment transportation and storage be made in its original packaging.

 Transport carefully protecting equipment from falls and impacts.

 Protect from moisture, rain exposure and direct contact with liquids.

 Keep it sheltered from the sun.

 Do not heap up more than 5 volumes.

 Do not move or store the equipment on uneven surfaces.

 Temperature range for transportation and storage: -10° C to +45° C.

 Moisture Limits for Transport and Storage: 20% to 70%.

13.4 - Contraindications

 This equipment is not intended for any use other than that for which it is intended, or to be operated by non-authorized personnel.

13.5 - Finishing

At the end of the workday, observe the following instructions:

- Provide equipment hygiene, cleaning the upholstery and its options, if available, as indicated in chapter 9 of this manual (*Cleaning and Disinfection*);
- Put the seat in the lowest position and turn off the equipment On/Off switch;
- Disconnect the electrical circuit breaker from the mains supply to the equipment.

14 - Troubleshooting

For solving possible problems in a simple and practical way, just follow the instructions in the following tables:

Item	Problem	Causes	Solutions
1	The chair doesn't perform any command	1° - The equipment is not connected to the main electricity	1° - Connect the equipment to the main electricity
		2° - Electrical circuit breaker is switched Off	2° - Switch On the electrical circuit breaker
		3° - There is no power on the main electricity	3° - Call the power supply company
		4° - The protection fuse is broken	4° - Contact Olsen certified assistance
2	Operating light is not lighting	1° - Equipment is not connected to the mains electricity	1° - Check troubleshooting <i>item 1</i>
		2° - The LED is burned	2° - Contact Olsen certified assistance
3	Ejector is weak or loses suction during procedure	1° - The ejector filter is obstructed	1° - Clean the ejector filter
		2° - Insufficient air pressure for the equipment	2° - Open the air register of the equipment
		3° - Drain hose obstruction	3° - Release bent/wrinkled hose
		4° - The sewage presents obstruction	4° - Provide sewage clearance
		5° - Blockage in the hydropneumatic system	5° - Contact Olsen certified assistance

If you are in doubt or find a problem with the equipment that is not mentioned in this chapter, stop using the equipment immediately and contact your authorized service center or contact us by e-mail export3@olsen.odo.br or by calling +55 48 2106 6000.

15 - Preventive Review

In order to extend your equipment lifespan, Olsen has prepared a list of the main equipment items for which it recommends performing semi-annual preventive maintenance.

Performing a preventive maintenance by an accredited technician does not interfere with the equipment's warranty period.

STANDARD ITEMS
Checking the multi-articulate headrest
Checking all joystick commands
Checking the swivel armrests
Checking the spindle drive motors and joints
OPTINOS ITEMS
Checking suction of Venturi ejector and/or Vacuum Pump
Lubrication of rings of the debris separator filter
Bowl ring disassembly, cleaning and lubrication
Checking Concept LED or Surgical operating light joint movements
Check the Surgical operating light and intensity of the Concept LED operating light
Check LED's
Checking the electrical, water, air and vacuum connections of the junction box
Checking the sewer connection (damper) of the junction box

Olsen recommends that the items in this chapter be reviewed every 180 days to prevent possible equipment failure or performance loss even after the equipment warranty expires.

Performing preventive or corrective review by an accredited technician does not interfere with the warranty period of the equipment.

Allow only qualified Olsen technicians to perform installation and maintenance on your equipment and accessories.

Use only genuine Olsen parts and accessories. The use of non-genuine parts, accessories or other components may compromise equipment performance by increasing its emissions or reducing its electromagnetic immunity.

Do not make any adaptations or changes to the equipment or its components or accessories.

16 - Warranty Terms

The warranty period is of 12 months starting from the equipment installation date only if the installation be made during 90 days after the purchasing date and under the terms of this certificate.

- 1** - The maximum storage time must be 3 months from the purchasing date. In case this time is exceeded, the warranty period will start even if the product is kept in storage.
- 2** - Upholstered parts are under warranty for 6 months.
- 3** - Multifaceted mirrors, fuses, wires and transformers are not covered by this warranty.
- 4** - The warranty is limited to the repair or replacement of defected parts and does not cover defects originated by:
 - a) Non-observance of the use and maintenance instructions;
 - b) Effects caused by the agents of nature;
 - c) Falls, crashes and an inadequate storage;
 - d) Damages caused by inadequate use of chemical products, mainly if not indicated on this user´s manual;
 - e) Contact with tissues, leather, disposable gloves, painting, pigmented detergents, razors, any sharp handpieces, etc. that can change the equipment original features;
 - f) A connection to a wrong voltage power supply.
- 5** - This warranty will cease when:
 - a) The warranty period expires;
 - b) The user makes any changes to the product or use of inappropriate accessories;
 - c) The warranty certificate presents any scraping or modifications;
 - d) The installation and technical service is done by an unauthorized technician;
 - e) Interrupting or failing to perform the maintenance previously scheduled;
 - f) The equipment is storage for more than 6 months after the invoice issue date.
- 6** - The repairing or parts replacement during the warranty period will not extend its original expiring date.
- 7** - The expenses originated from the equipment installation, scheduled maintenance, travel and hotel of the service staff involved in the calls for service for installation or repair of the equipment will run under the owner's responsibility and in accordance with the distributor's norms.
- 8** - The owner, after checking the installation when it is done, must keep the purchasing invoice and the serial numbers of each equipment until the end of the warranty period, as in the case of need of service the serial numbers, the invoice number and installation or purchasing date will be required.
- 9** - On all calls for service during the warranty period, the unit serial numbers and a copy of the purchasing invoice or installation form must be presented. In the case of this information is not presented the service will be done out of the warranty.


17 - Message from the President

Olsen and its clients:
A successful relationship.

I have linked my name to the factory and to the Dental and Medical equipment that are currently produced and trade in more than 100 countries having in mind the responsibilities and long-term response to this initiative.

Our products are modern, innovative, durable and of low maintenance cost. These characteristics have been achieved thanks to our competent and dedicated team, which make me very proud for many reasons, as they are always giving the best of their creative capacity to our clients.

Our company will always be open to all those who prefer Olsen products, for any necessary information and technical assistance, but especially for comments regarding the relationship with customers. We expect this connection always brings you satisfaction, resulting in more and more benefits to all of us.

Cesar Olsen

www.olsen.odo.br

+55 48 2106 6000
export3@olsen.odo.br

/// Olsen

Equipment made to last

Registration at Ministry of Health of Brazilian 10281300013

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