USER MANUAL

PERFORMANCE GYNECO TABLE

Rx only

Attention: Carefully read this manual
before operating the Olsen's equipment!

Olsen
Equipment made to last

Index

1 - Introduction	5
2 - Package Contents	5
3 - Equipment Presentation	5
3.1 - Standard Items	5
3.2 - Optional Items	5
4 - Parts Identification	6
4.1 - Applied parts	6
4.2 Accessories and Detachable Parts	6
5 - Description and Operating of the Equipment	7
5.1 - Turning On the Equipment	7
5.2 - Headrest and Auxiliary Pillow (Optional)	7
5.3 - Armrests (Optional)	7
5.4 - Remote Foot Pedal	8
5.4.1 - Movements Interruption	8
5.4.2 - Positioning of the Patient, the Operator and Others	8
5.5 - Emergency position	8
5.6 - Paper Roll Holder (Optional)	9
5.7 - Sliding Extender	9
5.8 - Upholstered Leg Rest	9
5.9 - Left and Right Bracket Locks	10
5.9.1 - Accessories on Bracket Locks	10
5.10 - Interaction Between Accessories	10
5.11 - Double Bracket Locks	11
5.12 - Colposcope Holder (Optional)	11
5.13 - Heels Support (Optional)	11
5.14 - Knees Support (Optional)	12
5.15 - Rear Castors Lock	12
6 - General Features	12
7 - Installation Requirements	13
7.1 - Pre-installation	13
7.2 - Electrical Installation	13
8 - Installation	13
8.1 Position	13
8.2 - Olsen Certified Technical Assistance Network	14
9 - Cleaning and Disinfection	14

9.1 - Upholsteries and Plastic Covers	14
9.2 - Painted Parts.....	14
9.3 - Sterilization by Autoclave.....	14
10 - Technical Features	15
10.1 - Electromagnetic Compatibility (EMC).....	16
11 - Dimensional.....	18
12 - Symbology.....	19
13 - Important Notes.....	19
13.1 - General Cares - Compulsory Reading	20
13.2 - Disposal.....	21
13.3 - Transport and Storage.....	21
13.4 Contraindication.....	21
14 - Troubleshooting.....	21
15 - Preventive Review	22
16 - Warranty Terms	23
17 - Message from the President.....	24

1 - Introduction

Congratulations on the great choice!

You have acquired an equipment developed for the outpatient environment, aimed primarily at the specialty of gynecology, aiming to provide maximum comfort to the patient and facilitate the performance of procedures by the professional.

This manual is supplemented to the Performance Gyneco Table Quick Operation Guide provided with your equipment and provides all the information you need to get the most out of your equipment, so read them carefully before using it.

2 - Package Contents

Check out the equipment package contents:

Standard Items:

- 1 Performance Gyneco Table
- 1 Left Bracket Lock
- 1 Right Bracket Lock
- 2 Upholstered Leg Rest
- 1 Performance Gyneco Table Quick Operation Guide

Optional Items:

- 1 Auxiliary Pillow
- 2 Upholstered Armrest
- 1 Colposcope Holder
- 2 Heel Support
- 2 Knee Support
- 1 Left Additional Bracket Lock
- 1 Right Additional Bracket Lock

3 - Equipment Presentation

3.1 - Standard Items

- Automatic Table with 6 Commands
- Headrest
- Interchangeable Seat Extension with Built-In Stainless Tray
- 2 Swivel Castors with Individual Locking
- 2 Fixed Castors
- 3 Oil-free Bosch Motors
- Upholstered and Interchangeable Leg Rest
- Remote Foot Pedal Commands

3.2 - Optional Items

- Auxiliary Pillow
- Upholstered and Interchangeable Armrests
- Paper Roll Holder
- Colposcope Holder
- Heels Support
- Knee Support
- Left Additional Bracket Lock
- Right Additional Bracket Lock

4 - Parts Identification

- A - Remote Foot Pedal
- B - Rear Castor with Lock
- C - Electrical Panel
- D - Left Bracket Lock
- E - Backrest Side Rail
- F - Paper Roll Holder
- G - Auxiliary Pillow
- H - Headrest
- I - Backrest

- J - Seat
- K - Leg Rest
- L - Armrests
- M - Extender Upholstery
- N - Seat Side Rail
- O - Right Bracket Lock
- P - Built-In Stainless-Steel Tray
- Q - Folding Step
- R - Free Front Castor

4.1 - Applied parts

The following items of the Performance Gyneco Table are considered as applied parts to the patient:

Standard Items:

- Upholstery;
- Upholstered ler rest;
- Seat extension;
- Extender upholstery;

Optional Items:

- Auxiliary pillow;
- Upholstery armrests;
- Knee support;
- Heel support;

4.2 - Accessories and Detachable Parts

The following items of the Performance Gyneco Table are considered as accessories or detachable parts or:

Detachable parts:

- Headrest;
- Upholstered leg rest.

Accessories:

- Auxiliary pillow;
- Additional left and right bracket locks;
- Upholstered armrest;
- Colposcope Holder;
- Heel support
- Knee support

5 - Description and Operating of the Equipment

Technical Name: Examination Table

Trade Name: Performance Gyneco Table

The Performance Gyneco Table is an automated equipment for patient accommodation on clinical and outpatient exams and gynecological procedures.

It is equipped with 3 electrical motors to perform 6 individual movements, controlled by remote foot pedal.

It performs clinical positioning for an emergency condition, allowing brain's irrigation by gravity, with slow movements to the maximum negative scale (-5° from horizontal) and can be moved easily because it has high resistance casters built in molded polyurethane that protect the floor and have very low noise level in the displacement. It has a foldaway step to facilitate patient access to the equipment.

This equipment should not be used to transport patients.


Before starting to use your equipment make sure that it is properly installed, connected to the mains, with the circuit breaker On and the mains powered.

5.1 - Turning On the Equipment

On the electrical panel, set the On/Off switch (1) to the "I" position. The LED (2) next to the switch will light green while the equipment beeps with 1 long beep, indicating that the equipment is ready for use.

The fuse holder (3) of the equipment is also located on the electrical panel.

To turn Off the equipment: press the On/Off switch (1) to the "O" position.


5.2 - Headrest and Auxiliary Pillow (Optional)

The headrest (1) was designed to be easily removed from the backrest according to the procedure need.

It has 2 fixing rods (2) for the fitting.

It can be accompanied by the auxiliary pillow (3) has 2 strips with Velcro® (4) for easy positioning and removal for cleaning.



5.3 - Armrests (Optional)

The armrests (1) have upholstery in the same color and characteristics as the Performance Gyneco Table.

It has adjustable rods (2) that allow the movements on side rails (3) for depth adjustment. It can also be used to assist in patient's lateral containment when he is lying on the Performance Gyneco Table.

To adjust the armrest (1): pull the metal lever (5) counterclockwise to loosen it.

Position the armrest (1) and turn the metal lever (5) clockwise to fasten it on the rail.


Note: the metal lever (4) moves only in the area indicated in the above figure. To turn it the metal lever (5) clockwise, push it upwards causing the lever to pass through the metallic knob's axis (4), and then rotate it again.

5.4 - Remote Foot Pedal

The remote foot pedal has the 6 commands of the Performance Gyneco Table identified with symbols that allow an intuitive operation. Conceived within the latest concepts of ergonomics and biosafety, it features foot-operated command and metal handle (7), that allows its drive and movement to be done with the feet, minimizing the risks of cross contamination.

- 1 - Seat Up
- 2 - Seat Down
- 3 - Backrest Down
- 4 - Backrest Up
- 5 - Recline Seat
- 6 - Incline Seat

5.4.1 - Movements Interruption

All Performance Gyneco Table movements are of continuous activation, i.e., while the operator presses the control button, the equipment realizes the movement. The movement is interrupted immediately when the operator releases the control button.

5.4.2 - Positioning of the Patient, the Operator and Others

During the equipment movement, the operator and other persons close to him/her should position themselves on the equipment sides, respecting the minimum distance of 50 centimeters, remaining out of the equipment movement area and its components.

To ensure safe and proper patient positioning, the operator should advise him to keep his hands in sight during the entire equipment movement on the equipment upholstery or on the specific accessories for performing the procedure. If the operator sees any situation that may generate any kind of risk to the patient, the chair must be stopped immediately.

All equipment automatic movements must be supervised by the operator.

The presence of the operator, persons or objects in the equipment movement areas and its components may cause damage to the equipment and/or impair its correct operation.

All equipment's automatic movements must be supervised by the operator. If the operator sees any situation that may generate any kind of risk to the patient, the chair must be stopped immediately.

5.5 - Emergency Position

Press the Backrest Down command (section 5.4) in the remote foot pedal, until the end of the course. The equipment will adjust the backrest position in -5° related to the floor.


The seat can be raised for lifting lower limbs.

5.6 - Paper Roll Holder (Optional)

Positioned on the back side of the backrest, this holder can be used with rolls up to 100 mm in diameter.


To put the sheet roll on the holder: turn the holder's shaft (1) clockwise, by pulling it in the opposite direction to the backrest. No need to remove the 2 shaft ends. Just remove 1 side only to put the disposable sheet roll.

Note: the roll of disposable sheet is not included.

5.7 - Sliding Extender

Located under the seat, it offers alternatives to facilitate the procedures through the extender upholstery (1) and stainless-steel tray (4). Both can be easily removed for cleaning and the tray can be sterilized by steam autoclave.

To use these items: pull the sliding extender (2) until its structure engages on the latches (3).

To use the stainless tray (4): remove the extender upholstery (1).

To use the upholstery (1): position it on the latches (3).

To retract the sliding extender (1): position the extender's tray (4) and upholstery (1) in its compartments and push the sliding extender (2) up to the end of the stroke.


5.8 - Upholstered Leg Rest

Manufactured with the same upholstery of the Performance Gyneco Table, they feature soft foam to provide greater comfort to the patient. With the aid of the bracket locks, the upholstered leg rest (1) feature height and depth adjustment.

The leg rest features upholsteries (1) fitted on the metal support (2), allowing a free spin on the fitting shaft (3).

The fitting shaft has joint (4), which gives the piece the inclination adjustment.

To adjust the inclination: position the upholstery on the requested position.

5.9 - Left and Right Bracket Locks

Used for attaching accessories to the side rails. They are identified according to the side of the seat to be installed. The bracket lock for the left side (1) receives the letter "L" and the right one (2) receives the letter "R".

To adjust the height (6) or turning (7) of the accessory, turn the plastic handle (3) counterclockwise. Adjust the accessory to the desired position and tighten the plastic handle (3) by turning it clockwise.

To adjust the position of the rails (8), pull the metal lever (4) counterclockwise by loosening the bracket on the side rail. Adjust the accessory to the desired position and turn the metal lever (4) clockwise.

Note: the metal lever (4) moves only in the area indicated in the figure to the side. To turn it the metal lever (4) clockwise, push it upwards causing the lever to pass through the metallic knob's axis (5), and then rotate it again.

5.9.1 - Accessories on Bracket Locks

Accessories which can be installed on the Performance Gyneco Table by using bracket locks:

Through the left and right bracket locks, these accessories can be adjusted in height, depth and rotation on the axis.

5.10 - Interaction Between Accessories

The Performance Gyneco Table allows the attendance of patients both lying down and sitting. In each of these conditions, the installation of the accessories to the side rails must be carefully observed to avoid damages to procedures, equipment and patient.

When moving the equipment, ensure that the accessories do not collide with one another or with other elements of the equipment, like figure be side.

The positioning of the accessories on the side rails must be made according to the needs to carry out each procedure, but the interaction between the accessories, during the equipment's movement it's operator's responsibility.


5.11 - Double Bracket Locks

The double bracket locks allow the coupling of 2 accessories, taking up less space on the side rails. The plastic handles (1) allow height and rotation on the shaft (distance the accessory to seat) adjustment, while the metallic lever (3) allows the bracket lock adjustment on the side rails.

Simply turn counterclockwise to loosen and clockwise to lock.

To adjust the bracket locks position on the side rail, pull the metal lever (3) counterclockwise. Adjust the attachment to the desired position and press the metal lever (3) clockwise.

Note: the metal lever (3) moves only in the area indicated in the figure above. To turn it the metal lever (3) clockwise, push it upwards causing the lever to pass through the metallic knob's axis (2), and then rotate it again.


5.12 - Colposcope Holder (Optional)

The Colposcope holder is attached to the metal support (4) of the upholstered leg rest or knee support, optimizing the professional's work space. This accessory can receive Colposcopes with pin fitting up to 22.5 mm diameter (7/8").

It has arm (3) with free joint, allowing better positioning of the Colposcope (5).

To attach the Colposcope: fit it and turn the knob (1) clockwise by locking it in the desired position.

To adjust the height (6): turn the knob (2) to loosen the regulatory bar, place the Colposcope at the requested height and turn the knob (2) clockwise to lock in position.

5.13 - Heels Support (Optional)



They offer support for the heels and are available with double bracket locks (3) that allow the simultaneous fixation of the heels support (1) and the leg (section 5.8) rest or knees support (section 5.14).

The heels support (2) can be adjusted for height, heel pads position and movement on the side rails.

In addition to bracket locks adjustments, heel supports (2) can be adjusted for the heel pad distance according to the patient profile or procedure to be performed.

To fix the heel support (2): fit the support (2) into

the holes (4) along the metal support (1) according to the desired distance.

5.14 - Knees Support (Optional)

They allow patient's accommodation with bent legs, providing more space to perform procedures. With the bracket locks is possible to adjust the knee support height and the depth.

5.15 - Rear Castors Lock

The rear castors lock can be activated with the feet.

To lock the castor: press the lever (1) down.

To release the castor lock: push the lever (2) down.

The lock must remain active during the equipment use. Only disable the lock to move the equipment.

6 - General Features

- **Upholsteries:** are assembled on a very resistant frame, with soft foam and coated with flexible laminated PVC or leather, which besides being resistant, providing easy asepsis for the whole.
- **Mechanical structure:** manufactured in SAE 1020 rolled steel profiles and welded by MIG process, to guarantee the strength and durability.
- **Electrostatic painting:** it is applied to all structural metal parts of the equipment. The special polyurethane paint provides the equipment with high durability coating and has antibacterial properties according to JIS Z 2801: 2000, where in 24 hours, the bacterial reduction is greater than 99.9%.
- **Plastic covers:** made of high-strength ABS with acrylic cover, the plastic covers do not require painting and allowing polishing for recovery in case of wear or scratches.
- **Electrical system:** the equipment can work in frequencies of 50 or 60 Hz and can be configured to be connected to the voltages of 118/127/220/230 Volts by an authorized technician. The maximum supply voltage of electronic boards, motors, and other controls is 24 V. The electrical system has On/Off switch and protection fuses.
- **Spindle drive motor:** all motors used in the equipment are produced by *Robert Bosch of Brazil* and have advantages such as low noise, the absence of oil reservoir, uniformity in displacement, reduction of energy consumption and low cost for maintenance. The motors also have protection systems that act on the mechanical and electrical parts of the set.
- **Pressurized suspension:** the high pressure nitrogen cylinder relieves the load of the gear motor lifting system of the equipment in 100 kg. This means that a patient with this approximate weight is not acting on the electromechanical lifting components.

7 - Installation Requirements

7.1 - Pre-installation

The pre-installation should be guided by Olsen authorized assistance to ensure that the environment is adequate to receive the equipment, in accordance with the temperature, relative humidity and pressure required for the equipment correct operation (*chapter 11 - Technical Features*), as well as positioning in which it should be used. In this step, the electrical connections must be prepared to supply the equipment.

This equipment is not designed to be installed or operated in a surgical center.

7.2 - Electrical Installation

The power grid must be single-phase, has specific grounding and 10 A/3 mA DR circuit breaker. The circuit breaker must exclusively supply for the Performance Gyneco Table and must be easy and quick to disconnect from the mains. If the power grid presents voltage variation, the installation of a surge protection device is required. The power switch must be turned off at the end of the workday.

There is a table below for sizing the electrical installation:

Tension (V)	Wire Diameter (mm ²)	Distance (m)	Current (A)
118/127/220/230	2,5	20	5,0

This equipment should only be connected to a power source with protective grounding. There is a risk of electric shock!

8 - Installation

The Performance Gyneco Table installation must be done by Olsen certified technician consists in assembling legs rest and checking the mains voltage which will supply the equipment and, if necessary, make the tension adjustment. In the case of optional items have been purchased, these should also be installed at this time.

The technician will also verify if the equipment complies with the purchasing order, if it maintains its integrity, and will provide operation, cleaning and maintenance guidance.

When the installation is complete, observe the following items:

- Both right and left bracket locks are installed with the upholstered leg rest;
- All remote foot pedal commands work perfectly;
- The sliding extender has the stainless-steel tray and the upholstery and slides perfectly;
- The upholsteries of backrest, seat, headrest and leg rest are intact.

In the case of optional items have been purchased, check the items below according to the acquired optional item:

- The armrests are installed on the side rails;
- The Colposcope holder was installed and performs all the joints perfectly;
- The knee supports were installed on their respective bracket locks;
- The heel supports were installed on their respective bracket locks;
- Integrity of upholstery and finishes;
- Joints without rigidity when moving;
- The fitting of accessory is secure.

For instructions on the proper functioning of the equipment and its components, read *Chapter 5 (Description and Operating of the Equipment)* of this manual.

Installation should only be made by an authorized technician. Installation performed by an unauthorized person will result in loss of warranty.

8.1 - Position

To position the equipment put the backrest parallel to the floor and raise the seat as far as possible, and then take it to the place where it will be used. Perform the movements of raising and lowering the seat, making sure that there is enough room for all movements.

8.2 - Olsen Certified Technical Assistance Network

To access the Olsen certified technical assistance network for installation and maintenance contact us by e-mail export3@olsen.odo.br or if you prefer +55 48 2106 6000.

9 - Cleaning and Disinfection

The whole sanitize process must be done with the use of gloves suitable for cleaning and protection, in addition to a mask and protective glasses, according to biosafety standards.

9.1 - Upholsteries and Plastic Covers

The plastic covers and upholsteries must be cleaned with a damp cloth containing neutral soap or detergent only. Olsen advises against using any chemical product to clean these parts, but in the case of disinfection products, it is important to check if it has suitable compatibility and specifications before use on these materials.

Never use hypochlorite or alcohol-based products.

9.2 - Painted Parts

It must be cleaned with a slightly damp cloth containing only soap or mild detergent.

Never use hypochlorite or alcohol-based products.

9.3 - Sterilization by Autoclave

The stainless steel tray can be sterilized by steam autoclaving at the following values (sterilization according to SS-374 Resolution of 12/15/1995 and ISO 17665):

- a) 130° C, 2 bar, 15 minutes;
- b) 120° C, 1 bar, 30 minutes;
- c) 134° C, 2,2 bar, 4 minutes.

Note: these items withstand up to 1.000 autoclaving cycles.

All equipment items mentioned in this chapter must be sanitized and sterilized (when appropriate) prior to use.

Olsen is not responsible for defects, deformities, stains or abnormalities caused by improper use of chemical products, contact with tissues, leather, disposable gloves, inks, pigmented detergents and other organic or synthetic products.

10 - Technical Features

Power supply: 118/127/220/230 VC~.

A certified technician when installing the equipment must select the voltage;

Note: all equipment is adjusted in the factory to 220 VC~.

Number of phases: single phase.

Frequency: 50/60 Hz.

Power 118/127 V: 200 VA.

Power 220/230 V: 250 VA.

Protection fuses:

- For 220/230 VC~: 1,5 A H (5 x 20 mm).
- For 118/127 VC~: >3 A H (5 x 20 mm).

Power cable specification (in compliance with requirements 6.1 and 6.2 of IEC 60601-1-2:2010):

- Flexible Cable PP Circular 500 V 3 x 1 mm 247-5 NM 53-C5.
- Tri-polar Male Plug 10 A - 250 V (NBR 14136).
- Tri-polar Female Plug 10 A - 250 V (IEC 60083/75).

Electrical shock protection type: Class I equipment, according to IEC 60601-1-1 and IEC 60601-1-2 standards.

Protection degree: B Type;

Operation mode: non-continuous.

Time On: 30 s.

Time Off: 5 min.

Operating environment conditions:

- **Temperature:** between 15° and 30° C.
- **Pressure:** 75 kPa to 106 kPa.
- **Relative humidity:** 30% to 70% non-condensing.

Harmful water penetration's protection:

- **Equipment:** IPX0.
- **Pedal:** IPX1.

Thermal protection of the transformer: aperture with 130° C ± 3%.

Lifting capacity (maximum patient weight + accessories): up to 200 kg.

Equipment weight:

- **Net:** 155,7 kg (with accessories: 177,4 kg).
- **Gross:** 182 kg (with accessories: 207 kg).

The weight of accessories:

- Leg rest: 4,20 kg.
- Colposcope holder: 2,90 kg.
- Heel support: 5,10 k.
- Knee support: 4,20 kg.

10.1 - Electromagnetic Compatibility (EMC)

The Performance Gyneco Table needs special attention regarding electromagnetic compatibility and must be installed and put into use in accordance with the electromagnetic compatibility information presented in this chapter.

Radio frequency (RF) communication equipment, portable and mobile, can affect the Performance Gyneco Table.

Guidance and Manufacturer's Declaration - Electromagnetic Emissions		
The Performance Gyneco Table is intended for use in the electromagnetic environment specified below. The customer or the user should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic Environment - Guidance
RF emissions CISPR 11	Grupo 1	The Performance Gyneco Table is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings use for domestic purposes.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Non-applicable	
Emissions due to voltage fluctuations/scintillation emissions IEC 61000-3-3	Non-applicable	

Guidance and Manufacturer's Declaration - Electromagnetic Emissions - I			
The Performance Gyneco Table is intended for use in the electromagnetic environment specified below. The customer or the user should assure that it is used in such an environment.			
IMMUNITY Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact	± 6 kV contact	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
	± 8 kV air	± 8 kV air	
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines	± 2 kV for power supply lines	Mains power quality should be that a typical commercial or hospital environment.
	±1 kV for input/output lines	±1 kV for input/output lines	
Surge IEC 61000-4-5	± 1 kV differential mode	± 1 kV differential mode	Mains power quality should be that a typical commercial or hospital environment.
	± 2 kV common mode	± 2 kV common mode	
Voltage dips, short interruptions and voltage variations on the power supply input lines IEC 61000-4-11	<5% U_T (dip > 95% in U_T) for 0.5 cycle	<5% U_T (dip > 95% in U_T) for 0.5 cycle	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Performance Gyneco Table requires continued operation during power mains interruptions, it is recommended that the Performance Gyneco Table be powered from an uninterruptible power supply or a battery.
	40% U_T (60% dip in U_T) for 5 cycles	40% U_T (60% dip in U_T) for 5 cycles	
	70% U_T (30% dip in U_T) for 25 cycles	70% U_T (30% dip in U_T) for 25 cycles	
	<5% U_T (dip > 95% in U_T) for 5 s	<5% U_T (dip > 95% in U_T) for 5 s	
Power frequency (50/60 Hz) Magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE U_T is the a.c. mains voltage prior to applications of the test level.			

Guidance and Manufacture's Declaration - Electromagnetic Immunity

The Performance Gyneco Table is intended for use in the electromagnetic environment specified below. The customer or the user should assure that it is used in such an environment.

IMMUNITY Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Conducted RF IEC 61000-4-3	3 V _{rms} 150 kHz to 80 MHz	3V _{rms}	Portable and mobile RF communications equipment should not be used closer to any part of the Performance Gyneco Table, including cables, than the recommended separation distance calculated by the equation applicable to the transmitter frequency. Recommended separation distance $d = [1,2]^{\frac{2}{3}}\sqrt{P}$ 80 MHz to 800 MHz $d = [1,2]^{\frac{1}{3}}\sqrt{P}$ 800 MHz to 2,5 GHz $d = [2,3]^{\frac{2}{3}}\sqrt{P}$ Where P is the transmitter maximum output power rating in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range. ^b Interference may occur in the surroundings of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

a	Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measure filed strength in the location in which the Performance Gyneco Table is used exceeds the applicable RF COMPLIANCE LEVEL above, the Performance Gyneco Table should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Performance Gyneco Table.
b	Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended Separation Distances Between Portable and Mobile RF Communications Equipment and the Performance Gyneco Table

The Performance Gyneco Table is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Performance Gyneco Table as recommended below, according to the maximum output power of the communications equipment.

Rated Maximum Output the Power of the Transmitter (W)	Separation Distance According to Frequency of Transmitter (m)		
	150 kHz to 80 MHz	80 kHz to 800 MHz	800 kHz to 2,5 GHz
	$d = [1,2]^{\frac{2}{3}}\sqrt{P}$	$d = [1,2]^{\frac{1}{3}}\sqrt{P}$	$d = [2,3]^{\frac{2}{3}}\sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,37	0,37	0,74
1	1,2	1,2	2,3
10	3,7	3,7	7,4
100	12	12	23


For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

11 - Dimensional

Measured in Millimeters

12 - Symbology

The symbology complies with standards IEC 60601-1 and IEC 60878.

	Seat Up		Seat Down		Recline Seat
	Backrest Up		Backrest Down		Incline Seat
	Non-Sterile		«Off» (Power)		“On” (Power)
	Caution		General Warning Sign		Warning, Electricity
	Operating Instructions		Refer to the Instruction Manual		General Mandatory Action Sign
	General Prohibition Sign		Stepping Prohibited		Manufacturer
	Type B Applied Part		Serial Number		Alternating Current
	Earth (Ground)		Protective Earth (Ground)		Sterilizable Up to the Temperature Specified
	Humidity Limitation		Temperature Limitation		Fragile, Handle with Care
	Keep Away from Heat		Keep Dry		This Side Up
	Heaping Up	 Authorized Representative in the European Community			

13 - Important Notes

The reproduction and distribution of these instructions can only be made with prior permission from Olsen Indústria e Comércio S.A.

The technical features of the products described in this manual correspond to the time of its publication. Future technical improvements do not result in any right to update existing products.

The images presented in this manual are illustrative.

This equipment is designed to be free from interference from magnetic fields, external electrical influences, electrostatic discharges, pressure or pressure variation, provided that the equipment is transported, installed, operated and sanitized in accordance with the instructions for use contained in this manual.

13.1 - General Cares - Compulsory Reading

 Follow the instructions in *chapter 7* of this manual (*Installation Requirements*) to suit the electrical network where the equipment will be installed.

 Follow instructions for proper use of the equipment and its accessories as described in *chapter 5* (*Description and Operating of the Equipment*). Improper use could be harmful to equipment, which would not be covered by warranty.

 Follow the instructions in *chapter 9* (*Cleaning and Disinfection*) of this manual for daily cleaning of your equipment.

 Protect your equipment from direct exposure to sunlight. The equipment direct exposure to sunlight may cause premature aging of plastic covers and upholstery.

 Turn off the circuit breaker or unplug the equipment at the end of the working day.

 In case of remote foot pedal damage, discontinue use of the equipment, turn it off, and contact Olsen accredited assistance.

 The electric mains cable was developed for use only on the Performance Gyneco Table. The use of this component in other equipment may compromise the emissions and electromagnetic immunity of these.

 Only the authorized technician can replace the equipment's mains cable and internal fuses.

 This equipment is not suitable for use in the presence of a flammable anesthetic mixture with air, O₂ or Nitrous Oxide.

 This equipment should only be operated by physicians and nurses for clinical, outpatient, and gynecological procedures.

 This equipment is not intended for any use other than that for which it is intended or to be operated by non-authorized personnel.

 Use only the mains cable supplied with the equipment for mains connection. Use of cables other than those specified (*chapter 10 - Technical Features*) may result in increased emissions or reduced immunity of the Performance Gyneco Table.

 Do not remove the plastic covers from the equipment. The risk of electric shock. Only the accredited technician is authorized to perform this procedure.

 In the case of damage to the electrical panel, adjacent plastic covers and motors plastic covers disconnect the equipment from the mains and contact the Olsen accredited service. The use of the equipment must be interrupted until maintenance is complete. The use of the equipment in these conditions offers a risk of electric shock!

 Do not perform maintenance or cleaning on the Performance Gyneco Table while it is in use or turned On.

 Do not install or use any electrical equipment over or near the Performance Gyneco Table. If necessary, the Performance Gyneco Table must be checked to see if it is functioning normally in the configuration in which it will be used.

 Do not perform the following procedures if it is possible to touch the patient, even if unintentionally, during the procedure:

- Fuse replacement;
- Do not perform the patient load or other procedures on the Performance table if the castors are not properly locked.

13.2 - Disposal

Debris, residues and infectious materials resulting from the procedures performed on this equipment must be deposited in biological waste duly identified and in accordance with current legislation.

For proper disposal of this equipment and its components and accessories, we recommend that it be sent to specialized recycling companies to ensure the best destination of each component without harm to the environment.

13.3 - Transport and Storage

It is recommended that the equipment's transportation and storage to be made in its original package.

Transport carefully protecting the equipment from falls and impacts.

Protect from moisture, rain exposure and direct contact with liquids.

Keep it sheltered from the sun.

Do not heap more than 4 volumes.

Do not move or store the equipment on uneven surfaces.

Temperature range for transportation and storage: -10° to +45° C.

Moisture limits for transport and storage: 20% to 70%.

13.4 Contraindication

This equipment is contraindicated for any use other than that for which it is intended, or to be operated by unauthorized personnel.

14 - Troubleshooting

For solving possible problems in a simple and practical way, just follow the instructions in the following table:

Item	Problem	Causes	Solutions
1	The equipment doesn't perform any command.	1° - Equipment is not connected to the electric mains	1° - Connect the equipment to the main electricity
		2° - Electrical circuit breaker is switched Off	2° - Switch on the electrical circuit breaker
		3° - There is no power on the main electricity	3° - Call the Power Supply Company
		4° - Protection fuse is blown	5° - Contact the Olsen authorized service

If you are in doubt or find a problem with the equipment that is not mentioned in this chapter, stop using the equipment immediately and contact your authorized service center or contact us by e-mail export3@olsen.odo.br or by calling +55 48 2106 6000.

15 - Preventive Review

In order to extend your equipment's lifespan, Olsen has prepared a list of the main equipment items for which it recommends performing semi-annual preventive maintenance.

Performing a preventive maintenance by an accredited technician does not interfere with the equipment's warranty period.

STANDARD ITEMS
Check all remote foot pedal commands
Check the sliding extender
Check the leg rest
Check and cleaning of all castors
Spindle drive motors and joints lubrication
Check the left and right bracket locks

OPTIONAL ITEMS
Check Colposcope holder joints
Check the double bracket locks
Check the armrests
Check the knees support bracket lock
Check heels support bracket lock
Check the shaft of paper roll holder lubrication

Allow only qualified Olsen technicians to perform installation and maintenance on your equipment and accessories.

Use only Olsen original parts and accessories. The use of non-original components can compromise the performance of the equipment, increasing its emissions or reducing its electromagnetic immunity.

Do not make adaptations, modifications or changes to the equipment or its components or accessories.

16 - Warranty Terms

The warranty period for this product is 12 months, counted from the equipment installation date, considering the 90 days legal warranty term, provided that the installation is performed within the period of 90 days from the product's Purchase Invoice issue and the fulfilled of other requirements of this certificate.

- 1** - The maximum term of storage is 3 months from the date of purchase of the product. If the storage period is exceeded, the guarantee is still in progress, even if the product is still stored.
- 2** - Upholstered parts are guaranteed for 6 months.
- 3** - Fuses, cables and transformers are not covered under warranty.
- 4** - The warranty is limited to repair or replacement of defective parts and does not cover defects originated by:
 - a) Non-compliance with the instructions for use and maintenance;
 - b) Falls, crashes and inadequate storage;
 - c) Action of nature agents;
 - d) Damage to upholstery, improper use of chemicals, exposure to unsuitable weather conditions, contact with tissues, leather, disposable gloves, paints, pigmented detergents, razors or any sharp instruments;
 - e) Damage to painted parts and plastic covers caused by improper use of chemicals or by contact with disposable gloves and sharp or piercing objects;
 - f) Connection to wrong voltage power supply.
- 5** - This warranty will cease when:
 - a) In the normal course of its period of warranty;
 - b) Make changes in the product not authorized by Olsen;
 - c) Adulterations in the purchase, installation or services document;
 - d) Installation or technical assistance made by a person not authorized by Olsen;
 - e) Failure to install the equipment for more than 90 days, counted from the purchasing date contained in the invoice;
 - f) By using non-genuine spare parts.
- 6** - The parts repairing or replacement during the warranty period will not extend the original expiry date.
- 7** - The expenses originated from the unit installation, scheduled preventive maintenance, travel, and hotel of the service staff involved in the calls for service for installation or units repair will run under the unit's owner responsibility and in accordance with the distributor's norms.
- 8** - The purchaser, after verifying the services performed in the installation and revision of the equipment, must date and sign the service order provided by the technician and keep along with his invoice of purchase of the equipment, failing which the warranty extension product when necessary.
- 9** - All service requests for warranty equipment must be made with the serial number of the equipment to be serviced and a copy of the purchase or installation document. If this information is not communicated, the service request will be made as not covered by the warranty.


17 - Message from the President

Olsen and its clients:
A successful relationship.

I have linked my name to the factory and to the Dental and Medical equipment which are currently produced and trade in more than 100 countries having in mind the responsibilities and long-term response to this initiative.

Our products are modern, innovative, durable and of low maintenance cost. These characteristics have been achieved thanks to our competent and dedicated team, which make me very proud for many reasons, as they are always giving the best of their creative capacity to our clients.

Our company will always be open to all those who prefer Olsen products, for any necessary information and technical assistance, but especially for comments regarding the relationship with customers. We expect this connection to always bring you satisfaction, resulting in more and more benefits to all of us.

Cesar Olsen

www.olsen.odo.br

+55 48 2106 6000
export3@olsen.odo.br

/// Olsen

Equipment made to last

Registration at Ministry of Health of Brazilian 10281300011

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